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These Instructions for Use are not  
intended for distribution in the USA.

Instruments and implants  
approved by the AO Foundation.

**DePuy Synthes**

COMPANIES OF **JOHNSON & JOHNSON**

### Table of contents

1. Introduction
2. Intended use
3. Usage
4. Cooling with irrigation liquid
5. Implant removal and teeth cutting with cutting tools
6. Locating of the cutting tools or fragments of cutting tools
7. Reprocessing instructions
8. Visual inspection and disposal of blunt or damaged reusable cutting tools
9. Warranty
10. Explanation of symbols

## 1. Introduction

Please read these Instructions for Use carefully prior to usage. It is essential that the surgeon and operating theater staff are fully conversant with the appropriate surgical technique for the cutting tools and associated implant.

The cutting tools are designed for use by physicians and trained medical personnel.

Synthes cutting tools are made of either stainless steel, tungsten carbide or are coated with nickel diamond. Each of them are designed for a specific surgical use.

These Instructions for Use cover the following Synthes cutting tools:

- Saw Blades, Rasps and Broaches
- Burrs
- Drill Bits (Radiolucent Drive Drill Bits and Drill Bits with Stop)

These instructions do not apply to Anspach cutting tools and Piezoelectric System cutting tips.

Detailed ordering information on the cutting tools can be found in the corresponding brochures of the cutting tools.

To ensure proper operation of the Synthes Power Tool, only use original Synthes cutting tools. Non Synthes cutting tools can considerably reduce the operating life of the system.

For further information on the Power Tool, please refer to the Instructions for Use of the corresponding system.

## 2. Intended use

The cutting tools for Power Tools are designed for use in surgery of the skeleton, i.e. for cutting, drilling, reaming and burring bone.

## Synthes Cutting Tools

Overview of the applications of the different types of cutting tools:

		To be used on		
		Bone	Teeth	Metal
Type of cutting tool	Saw Blades	×		
	Rasps	×		
	Broaches	×		
	Drill Bits (Radiolucent Drive Drill Bits and Drill Bits with Stop)	×		
	Burrs (stainless steel)	×		
	Burrs (tungsten carbide)*	×	×	×
	Burrs (stainless steel with a nickel diamond coating)*	×	×	×

\*For further information, please refer to chapter 7.

### 3. Usage

- The surgeon has to evaluate if the Power Tool and cutting tools are suitable for an application, based on power limitation of the machine, attachment and cutting tool regarding bone strength/anatomical situation as well as handling of the machine, attachment and cutting tool regarding bone size. In addition, the contraindications of the implant have to be respected. Please refer to the corresponding "Technique Guides" of the implant system used.
- Synthes recommends using a new cutting tool for each operation and discarding it after use. Cutting tools intended for single-use (⊗) must not be reused. Further information can be found in chapter 9.
- All cutting tools may only be used in combination with the appropriate Power Tool. Ensure prior to use that the coupling geometry of the cutting tool is compatible with the coupling geometry of the attachment or Power Tool.
- Cutting tools should only be used in combination with a well maintained and reprocessed Power Tool. Check proper operation of the Power Tool before using it on the patient. Refer to the Instructions for Use of the corresponding Power Tool.

- Cutting tool article numbers and dimensions are stated on the packaging label.
- A letter code on the burrs indicates the maximum speed recommended for the respective burr. Details about these codes can be found in the Instructions for Use of the corresponding Power Tool.
- Cutting tools should be carefully inspected before use to ensure that they are fully functional. Scratches or dents can result in breakage. Dullness of cutting edges can result in poor functionality. Damaged cutting tools should be replaced to prevent potential patient injury such as loss of metal fragments into the surgical site.
- Before use, make sure that the cutting tools are fixed correctly. Cutting tools, which are not inserted to their full chucking length, may cause high centrifugal forces and vibrations.
- For any additional information regarding the use of the cutting tool combined with a specific Power Tool, please refer to the Instructions for Use of the corresponding system.
- In order to ensure a smooth cut, the Power Tool has to be running before the cutting tool comes into contact with bone. Excessive contact pressure should be avoided as this might lead to:
  - damage to the cutting tool surfaces
  - splintering of the teeth of the cutting tool
  - thermal necrosis caused by excessive heat build-up
  - reduced lifespan of cutting tools
- Do not use cutting tools for working on metal (e.g. implants) and teeth, except for burrs made of tungsten carbide or with a nickel diamond coating. For additional details please refer to chapter 7.
- When cutting blocks are used, only saw blades, which do not exceed the thickness of the cutting block slots, should be used. Prior to use on the patient, free mobility of the saw blade has to be assured. Please refer to the corresponding "Technique Guide" of the implant system used.
- Jamming, leverage and bending of the saw blade must be avoided in the cutting block slot. The saw blade will otherwise overheat, causing it to jam. Other risks include fracture and thermal necrosis.

## Synthes Cutting Tools

- Contact with the cutting block, guiding device or other metallic objects (e.g. instruments) has to be avoided. Otherwise, damage to the cutting tools can occur and lead to injury of bone and surrounding tissue.

### Precautions

- Do not use the cutting tool if the package is open and/or damaged.
- Do not use the cutting tool after its expiration date.
- Do not use the cutting tool if damage is apparent or if the markings are not readable.
- Should the cutting tool be corroded, do not use it anymore.
- Ensure that the cutting edge of the cutting tool does not damage the sterile barrier (e.g. sterile cover) or gloves.
- It is recommended that a second cutting tool is available during application in case the cutting tool is dropped.
- The use of irrigation liquid is recommended to cool the cutting tools and prevent heat necrosis. Further information can be found in chapter 6.
- All operating room personnel must wear eye protection.
- If a product drops on the floor, fragments may split off. This represents a danger for the patient and user as these fragments may be sharp and/or unsterile fragments may enter the sterile field or hit the patient.

### Warning

- Improper use of the cutting tool may lead to injury to the patient, operator or other persons.
- Improper use of the cutting tool may lead to premature wear and destruction of the instrument.
- Cutting tools with a nickel diamond coating or containing stainless steel may elicit allergic reactions in patients with hypersensitivity to nickel. Ensure prior to use that the patient does not suffer from this allergy.

#### **4. Cooling with irrigation liquid**

Always irrigate adequately during cutting to prevent overheating of the saw blade or bone. Insufficient cooling may lead to increased heat generation (e.g. the saw blade teeth may become clogged with debris). At worst, this may cause irreversible damage to the bone due to thermal necrosis.

The Electric Pen Drive and Air Pen Drive are specifically equipped for irrigation:

- The console of the Electric Pen Drive System is equipped with an irrigation pump.
- For the Air Pen Drive System, an irrigation control unit is available.
- Specific irrigation nozzles and irrigation tube sets are available for the Electric and Air Pen Drive attachments.

For other Power Tools, a syringe may be used for irrigation.

For further information, please refer to the Instructions for Use of the corresponding Power Tool.

#### **5. Implant removal and teeth cutting with cutting tools**

##### **Nickel diamond coated or tungsten carbide burs**

The small torx cutting tools with a nickel diamond coating or made of tungsten carbide are indicated for the cutting, shaping and smoothing of bones, teeth and metal.

##### **Precautions to consider when removing implants or cutting teeth**

- Implants should only be removed with cutting tools if there is no other solution for implant removal.
- Only use nickel diamond coated or tungsten carbide cutting tools.
- Remove all particles by continuous flushing and vacuuming.
- Soft tissue must be well covered.
- Observe the material composition of the implant.

#### **6. Locating of the cutting tools or fragments of cutting tools**

Synthes cutting tools are designed and manufactured to perform within the scope of their intended use. However, if a metallic instrument (e.g. stainless steel, tungsten carbide, etc.) breaks during use, a medical imaging device (e.g. CT, Radiation Device, etc.) can aid on locating fragments and/or components of the cutting tool.

## Synthes Cutting Tools

### 7. Reprocessing instructions

Synthes recommends using a new cutting tool for each operation and discarding them after use. Performing cuts with a new and sharp cutting tool is faster, more precise, prevents infections caused by residues and generates less heat. This leads to a shorter surgery time, a reduced risk of bone necrosis and a better result.

#### Precautions

- Before first use, non-sterile cutting tools have to run through the complete reprocessing procedure.
- Special caution must be taken with cutting tools due to the risk of injury and infection of the personnel.

#### Single-use cutting tools

Cutting tools intended for single-use (⊗) must not be reused. Reuse or reprocessing (e.g. cleaning and resterilization) may compromise the structural integrity of the device and/or lead to device failure, which may result in patient injury, illness or death.

All single-use cutting tools (⊗) have to be discarded after the operation.

#### Precaution

- Cutting tools with a nickel diamond coating or made out tungsten carbide cannot be reprocessed since these burrs cannot be cleaned properly and the coating can chip off if used several times. They must therefore be discarded after each use.

#### Warning

Reuse or reprocessing of single-use device may create a risk of contamination, for example due to the transmission of infectious material from one patient to another. This could result in injury or death of the patient or user.

#### Reusable cutting tools

Reusable cutting tools may be reprocessed several times. However, cutting tools are frequently exposed to high mechanical loads and shocks during use and should not be expected to last indefinitely.

Consult the clinical processing instructions (Clinical Processing of Cutting Tools, 036.000.499). The reprocessing instructions are not applicable to cutting tools labeled as single-use items (⊗).

## **8. Visual inspection and disposal of blunt or damaged reusable cutting tools**

Inspect reusable cutting tools for blunt or damaged teeth using a magnifying glass (minimum enlargement 1:10) after each use and replace if necessary.

Check especially for:

- splintering and blunt (rounded) teeth, flutes, etc.
- damage to the cutting tools
- bent and corroded cutting tools

### **Precautions**

- Blunt or damaged cutting tools must no longer be used.
- Cutting tools must not be reground.
- The used cutting tools are to be disposed of according to internal hospital policies and procedures and/or current national recommendations.

### **Warning**

Failure to follow these precautions may lead to damage to the cutting tool and increased risk for the patient, the operator or third parties.

## **9. Warranty**

The manufacturer does not accept liability for damage resulting from improper use and failure to follow the Instructions for Use.

Synthes has not tested the compatibility of the Synthes cutting tools with Power Tools provided by other manufacturers and assumes no liability in such cases.

For further information please contact your local DePuy Synthes sales representative.

## **10. Explanation of symbols**

The symbols which are listed at the end of these Instructions for Use are shown on either the cutting tools or the labeling.

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60121204 SE\_542547 AB 03/2016 © Synthes, Inc. or its affiliates
Subject to modifications Synthes is a trademark of Synthes, Inc. or its
Synthes GmbH