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DePuy Ireland UC
Loughbeg
Ringaskiddy
Co. Cork Ireland

DePuy Synthes

PART OF THE JOHNSON-&JOHNSON FAMILY OF COMPANIES

Instruments and implants
approved by the AO Foundation

Instructions for Use

Please read these instructions for use, any enclosed information and relevant surgical techniques carefully before use. Ensure that you are familiar with the appropriate surgical technique.

The “Vertecem V+ Syringe Kit” is a sterile-packed, single-use syringe kit used for the injection of DePuy Synthes PMMA-based cement into cancellous bone.

Since the “Vertecem V+ Syringe Kit” is used in technically demanding procedures, the surgeon should be familiar with the use of the syringe kit as well as with the rest of the instrumentation and the procedure.

Product Name and Materials

“Vertecem V+ Syringe Kit”. One kit contains the following parts:

| Product name | Number of units | Materials |
|----------------------------------|------------------------|--|
| Stop-cock | 1 unit | PA12 (Polyamide) POM (Polyoxymethylen) |
| 2 mL syringes, white plungers | 5 units | PA12 (Polyamide) ABS (Acrylonitrile Butadiene Styrene) Silicone |
| 1 mL syringes, blue plungers | 8 units | PA12 (Polyamide) ABS (Acrylonitrile Butadiene Styrene) Silicone |

Intended use

The “Vertecem V+ Syringe Kit” is solely intended for the application of PMMA-based bone cement with the intent to augment cancellous bone. Refer to the corresponding directions regarding indications, contraindications, compatibility, use, precautions, warnings and side effects of the PMMA-based bone cement and access solution used in conjunction with the “Vertecem V+ Syringe Kit”. For information on compatibility with other devices or systems, consultation with a DePuy Synthes representative is recommended.

Side effects

There are no known side effects or interactions associated with the use of the “Vertecem V+ Syringe Kit” together with PMMA-based bone cement.

Sterile device

STERILE EO Sterilized using ethylene oxide

Store instruments in their original protective packaging and do not remove them from the packaging until immediately before use. Prior to use, check the product expiration date and verify the integrity of the sterile packaging. If all packaging appear to be intact, remove product(s) from the package using aseptic technique. Do not use if the package is damaged or the expiration date is passed.

Do not re-sterilize

Single-use device

Do not re-use

The “Vertecem V+ Syringe Kit” is intended for single-use only and must not be re-used. Re-use or reprocessing (e.g. cleaning and/or resterilization) may compromise the structural integrity of the device and/or lead to device failure which may result in patient injury, illness or death.

Furthermore, reuse or reprocessing of single-use devices may create a risk of contamination e.g. due to the transmission of infectious material from one patient or person to another. This could result in injury or death of the patient or user.

Device-specific warnings

- The purpose of the milliliter (mL) scale on the side of the “Vertecem V+ Syringe Kit” syringes is to give physicians an indication of how much PMMA-based bone cement has been injected. This scale is for information purposes only.

Precautions

- The “Vertecem V+ Syringe Kit” is used in technically demanding procedures. Therefore only physicians familiar with the proper technique and instrumentation used for the delivery and use of PMMA-based bone cement should use the “Vertecem V+ Syringe Kit”.
- The “Vertecem V+ Syringe Kit” should only be used to deliver PMMA-based bone cement in conjunction with the use of fluoroscopy imaging equipment capable of delivering high-quality images.
- Ensure a good fit between “Vertecem V+ Syringe Kit” syringe and “Vertecem V+ Syringe Kit” stop-cock/used access solution, but make sure to be on axis and avoid using excessive force when coupling them. They are both made of plastic and could otherwise break.

- After filling a “Vertecem V+ Syringe Kit” syringe with cement, grip it by its wings. Gripping it anywhere else could transfer heat from your hand to the cement, shortening the working time of the PMMA-based bone cement.
- Always fill all “Vertecem V+ Syringe Kit” syringes with PMMA-based bone cement directly after mixing.
- Always make sure the PMMA-based bone cement has reached the preferred viscosity prior to injection.
- To halt the flow of cement out of the “Vertecem V+ Syringe Kit” syringe at any given time, simply stop pushing the plunger.
- Make sure to always inject cement in a slow, controlled fashion.
- If cement flow is hindered at any time during injection, halt, investigate and correct the reason for the flow hindrance.
- Do not use tools or excessive force to operate the “Vertecem V+ Syringe Kit”.
- Caution: Federal law restricts this device to sale in the US by or on the order of a physician.

Combination of medical devices

The “Vertecem V+ Syringe Kit” is tested and validated with the DePuy Synthes PMMA-based bone cement and DePuy Synthes augmentation hardware only. DePuy Synthes has not tested compatibility with devices provided by other manufacturers and assumes no liability in such instances. For information on compatibility with other devices or systems, consultation with a DePuy Synthes representative is recommended. Refer to corresponding device’s or system’s directions for use.

Operating instructions

1. Prepare PMMA-based bone cement according to the corresponding directions of use and connect the "Vertecem V+ Syringe Kit" stop-cock using the side with no funnel. Ensure a good fit between "Vertecem V+ Syringe Kit" stop-cock and cement mixing device but make sure to be on axis and avoid using excessive force when attaching them – they are both made of plastic and could otherwise break.
2. Next air has to be removed from the system. Hold the cement mixing device vertically and gently turn its handle clockwise. Note: To remove air, rotate the handle clockwise in a controlled manner – do not push it.
3. You will see the mixer piston advancing inside the translucent cartridge and eventually a steady flow of cement will move into the "Vertecem V+ Syringe Kit" stop-cock. As soon as the cement starts to exit the "Vertecem V+ Syringe Kit" stop-cock at the funnel end, close it by turning the lever ("off") 90°, toward the cement mixing device. This should stop the cement flow.
4. Attach a "Vertecem V+ Syringe Kit" syringe to the stop-cock (funnel side). Ensure a good fit between "Vertecem V+ Syringe Kit" stop-cock and "Vertecem V+ Syringe Kit" syringe but make sure to be on axis and avoid using excessive force when attaching them – they are both made of plastic and could otherwise break.
5. Open the "Vertecem V+ Syringe Kit" stop-cock by turning the lever 90° to the left, back to its original position. Turn the cement mixer handle clockwise in a slow and controlled fashion to fill the "Vertecem V+ Syringe Kit" syringe. PMMA-

based bone cements can be highly viscous and need time to flow through the “Vertecem V+ Syringe Kit” stop-cock and into the “Vertecem V+ Syringe Kit” syringe. Make sure to fill the “Vertecem V+ Syringe Kit” syringe completely. Note: To transfer cement, rotate the handle clockwise in a controlled manner – do not push it.

6. As soon as the “Vertecem V+ Syringe Kit” syringe is filled, stop the cement flow by again closing the “Vertecem V+ Syringe Kit” stop-cock.
7. Disconnect the filled “Vertecem V+ Syringe Kit” syringe and attach the next one. Make sure to disconnect the “Vertecem V+ Syringe Kit” syringe on axis, otherwise the connection could break. Repeat steps 4 through 7 until all “Vertecem V+ Syringe Kit” syringes are filled. Always fill all “Vertecem V+ Syringe Kit” syringes immediately after mixing.
8. Test the cement viscosity by ejecting cement from one of the “Vertecem V+ Syringe Kit” syringes and evaluating it. Cement viscosity is an important factor in the safe delivery of cement to the intended operative site, as the risk for leakage decreases with increasing viscosity.
9. Once the viscosity is deemed appropriate by the responsible physician, the “Vertecem V+ Syringe Kit” syringe is attached to an access system and cement application can commence. Ensure a good fit between “Vertecem V+ Syringe Kit” syringe and access solution but make sure to be on axis and avoid using excessive force when attaching them – the “Vertecem V+ Syringe Kit” syringes are made of plastic and could otherwise break.

10. Exercise control when injecting PMMA-based bone cement. Always push the “Vertecem V+ Syringe Kit” syringe piston on axis, otherwise the connection between “Vertecem V+ Syringe Kit” syringe and access solution could break.

The force necessary to inject cement increases with time, as the cement sets. The required force also increases with increased syringe size. It is therefore advised to start using the 2 mL “Vertecem V+ Syringe Kit” syringes first and the 1 mL “Vertecem V+ Syringe Kit” syringes later.

Troubleshooting

If the force needed to inject cement becomes too high, this indicates that the maximum working time for the combination of DePuy Synthes PMMA-based bone cement and access solution used has been reached.

Device-related storage and handling information

Keep dry

Keep away from sunlight

Additional device-specific information

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© Synthes GmbH 2017. All rights reserved. **60051048 SE_235619 AF 01/17**
Synthes GmbH