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**SISTEMA DE CIMENTO VERTECEM V+**

## Instructions for Use

Please read these instructions for use, any enclosed information and relevant surgical techniques carefully before use. Ensure that you are familiar with the appropriate surgical technique.

The Vertecem V+ Cement Kit consists of two sterile packed components; a mixing device pre-filled with polymer powder and a glass ampoule containing monomer fluid. Also included in the kit is a lid for use during the mixing and transfer of bone cement into a suitable application or delivery system. For information on compatibility with other devices or systems, consultation with a DePuy Synthes representative is recommended.

### Article number

07.702.016S

### Article name

Vertecem V+ Cement Kit

### Material(s)

Containers for storage, transport and use:

Container	No. of units	Material(s)
Mixing device	1 unit	Polypropylene high density polyethylene terephthalate (contains powder component of bone cement)
Monomer ampoule	1 unit	Glass (contains liquid compo- nent of bone cement)
Mixing and transfer lid	1 unit	Polypropylene UHMW Polyethylene

Bone cement:

Component	Quantity	Constituent(s)
Polymer powder	26 g	<i>Polymethyl methacrylate/acrylate</i> <i>Zirconium dioxide</i> <i>Hydroxyapatite</i> <i>Benzoyl peroxide</i>
	44.6%	
	40.0%	
	15.0%	
	0.4%	
Monomer liquid	10 ml	<i>Methyl methacrylate, stabilized with</i> <i>60 ppm Hydroquinone</i> <i>N,N-dimethyl-p-toluidine</i>
	99.35%	
	0.65%	

### Intended use

The Vertecem V+ Cement Kit bone cement component is a radiopaque, injectable polymethyl methacrylate cement, indicated for use in spine pathologies being treated using vertebroplasty or kyphoplasty procedures. It

should only be used with devices or systems with which it has been tested and validated. For information on compatibility, consultation with a DePuy Synthes representative is recommended.

### Indications

- Vertebral compression fractures (VCF)
- Progressive compression fractures of one or multiple vertebrae with subsequent development of increasing kyphosis
- Patients with persisting instability after a vertebral fracture
- Combined procedures with internal fixation in osteoporosis
- Osteolysis
- Metastatic lesions
- Aggressive type of haemangioma

### Contraindications

- Infections
- Patients with clotting disorders
- Patients with severe cardiac and / or pulmonary insufficiency
- Patients with known hypersensitivity or allergy to any of the components of Vertecem V+ bone cement (see “Material(s)” above)
- Vertebra plana or circumstances where safe percutaneous access to the vertebra cannot be guaranteed
- Unstable vertebral fractures with involvement of the posterior wall in standalone vertebral augmentation procedures (e.g. vertebroplasty)
- Previous damage to the pedicle wall (transpedicular access)
- Lesions that feature narrowing of the spinal canal (more than 20%), including fracture or neoplasm, with or without myelopathy
- Retropulsing vertebral fragments with myelopathy
- A satisfactory response to conservative treatment
- Asymptomatic stable vertebral fractures
- Vertecem V+ bone cement is contraindicated for use in arthroplasty procedures

### Side effects

- Serious adverse events, some with fatal outcome, associated with the use of polymethyl methacrylate application include (but are not necessarily limited to) myocardial infarction, cardiac arrest, cerebrovascular accident, pulmonary embolism and anaphylaxis.
- The most frequent adverse reactions reported in conjunction with the application of polymethyl methacrylate are transitory drop in blood pressure, thrombophlebitis, hemorrhage and hematoma, superficial or deep wound infection, bursitis, short-term cardiac irregularities and heterotopic bone formation.
- Other potential adverse events reported with the application of polymethyl methacrylate include (but are not necessarily limited to) pyrexia, haematuria, dysuria, bladder fistula, transitory worsening of pain due to heat released during the polymerization, nerve entrapment and dysphagia due

- to extrusion of the bone cement beyond its intended application, adhesions and stricture of the ileum due to heat released during polymerization.
- Potential adverse events associated with vertebroplasty and kyphoplasty procedures include (but are not necessarily limited to) pneumonia, intercostal neuralgia, collapse of a vertebra adjacent to the one injected due to an osteoporotic disease, pneumothorax, cement extrusion into soft tissue, fracture of a pedicle, rib fracture in patients with diffuse osteopenia, especially during thoracic vertebroplasty procedure due to the significant downward force exerted during needle insertion, compression of the spinal cord with paralysis or loss of feeling and cement leakage into the vertebral disc. Deep wound infection is a serious complication and may require removal of the implant. It may lay dormant and may not manifest itself for several days or weeks postoperatively.
  - Adverse patient reactions affecting the cardiovascular system have been associated with the use of bone cements for arthroplasty. The absorption of methyl methacrylate into the vascular system may be associated with acute hypotensive effects. Hypotensive reactions have occurred shortly after implantation. They have lasted for a few minutes. Some patients have progressed to cardiac arrest. For this reason, patients should be monitored for any change in blood pressure during and immediately following the application of Vertecem V+ Cement Kit bone cement.
  - Polymethyl methacrylate may cause hypersensitivity among high-risk patients, which can result in an anaphylactic reaction.
  - Avoid uncontrolled or excessive bone cement injection, as this may cause cement leakage with severe consequences such as tissue damage, paraplegia or fatal cardiac failure.
  - Insufficient cement infiltration of the bone structure or postoperative events may affect the cement / bone interface, subsequently resulting in micro-motions between cement and bone. This could lead to the development of a layer of fibrous tissue between cement and bone, or even to bone resorption. Regular, long-term supervision is therefore recommended for all patients.
  - The hardening, or setting, of Vertecem V+ Cement Kit bone cement is an exothermic reaction, meaning it produces heat. The temperature may become high enough to damage bone or other tissue at or near the implant.
  - The liquid monomer and its fumes are aggressive and highly flammable. The operating room should be adequately ventilated to minimize the monomer fume concentration and care should be exercised to prevent exposure, as they may produce irritation of the respiratory tract, eyes and possibly liver. Highly concentrated monomer fumes may have an adverse reaction with soft contact lenses and therefore personnel wearing contact lenses should not be involved in the mixing of this product.

## Sterile device

The Vertecem V+ Cement Kit is supplied in sterile form. Prior to use, check the product expiration date and verify the integrity of the sterile packaging; if it is opened or otherwise damaged, do not use. If all packaging appears to be intact, remove product(s) from the package using aseptic technique.

The monomer fluid inside the ampoule is sterilized using an aseptic processing technique (ultra filtration).

The double blister and all its content is sterilized using ethylene oxide.

Do not re-sterilize

## Single-use device

Do not re-use

The Vertecem V+ Cement Kit is intended for single-use only and must not be re-used. Re-use or reprocessing (e.g. cleaning and re-sterilization) may compromise the structural integrity of the device and / or lead to device failure which may result in patient injury, illness or death.

Furthermore, reuse or reprocessing of single-use devices may create a risk of contamination e.g. due to the transmission of infectious material from one patient to another. This could result in injury or death of the patient or user.

A single Vertecem V+ Cement Kit unit may not be used to treat more than one patient.

## Precautions

- A thorough preoperative check-up of the patient must be carried out.
- Store this package between 0°C and 25 °C and protect it from light to prevent premature polymerization of the monomer fluid. Always check the condition of the monomer before use; if it shows any sign of thickening or premature polymerization, do not use.
- If stored below OR temperature, always allow the product to equilibrate to OR temperature prior to use.
- Make sure to follow preparation and handling instructions for the Vertecem V+ Cement Kit carefully.
- The monomer liquid is a powerful lipid solvent; it should not be allowed to come into direct contact with sensitive tissue or rubber or latex gloves. Wearing of a second pair of gloves and strictly adhering to the mixing instructions may reduce the risk of hypersensitivity reactions.

## Instructions for Use

- Always use the full amounts of monomer liquid and polymer powder provided in the kit, respectively, when mixing Vertecem V+ Cement Kit bone cement. Otherwise the behavior of the Vertecem V+ Cement Kit can no longer be guaranteed. Using only one of the components is not permitted.
- Ensure that the powder and liquid component are thoroughly mixed before starting cement transfer.
- The mixing device is designed for mixing of Vertecem V+ Cement Kit bone cement and transferring it to a compatible application system. It is not designed for the injection of mixed cement directly into a vertebra. Such use is strictly prohibited. For information on compatibility, consultation with a DePuy Synthes representative is recommended.
- The formation of large air pockets within the cement should be avoided when filling the application system.
- The physician should have specific training and experience in performing a vertebroplasty or kyphoplasty procedure. She or he should also be thoroughly familiar with the properties, handling characteristics and application of the Vertecem V+ Cement Kit, as well as with percutaneous cement delivery.
- The use of Vertecem V+ Cement Kit bone cement during pregnancy or lactation is not recommended. In considering use of Vertecem V+ Cement Kit bone cement in a pregnant patient, the responsible physician should weigh the benefits of the application of Vertecem V+ Cement Kit bone cement against the potential risks to the patient and the pregnancy.
- Long-term efficacy data for Vertecem V+ Cement Kit bone cement is limited. The responsible physician should weigh the benefits of using Vertecem V+ Cement Kit bone cement in painful non-osteoporotic, acute traumatic fractures against the potential risks, especially when treating younger patients.
- Exercise caution in cases involving significant fragmentation and collapse of the vertebral body (i.e. a fractured vertebral body height corresponding to less than 1/3 of the original height). Such cases may require treatment with a more technically complex procedure.
- Application of Vertecem V+ Cement Kit bone cement should be monitored using real time imaging procedures capable of delivering high quality images. Use appropriate imaging techniques to confirm correct needle placement, the absence of damage to surrounding structures and appropriate location and quantity of injected material.
- After last cement injection, the patient should remain immobile for 15 minutes to facilitate proper cement curing.
- If repeated access is required within a single vertebral body, the access needle should be closed with a trocar to maintain access and to prevent bone cement from leaking out of the vertebra through said needle.

**Warnings**

- Use of Vertecem V+ Cement Kit bone cement is not established for prophylactic augmentation of non-fractured vertebral bodies.
- Percutaneous application of polymethyl methacrylate should only be performed in medical settings where emergency decompressive surgery is available.
- During surgery Vertecem V+ Cement Kit bone cement may spread beyond the region of its intended use (i.e. outside of the vertebral body or into the circulatory system). If this occurs, injection should be halted immediately.
- Extravasation of bone cement can result in cement entering fracture cracks or venous circulation, but is not necessarily limited to said areas.
- Additives (such as antibiotics) are not to be mixed with the Vertecem V+ Cement Kit bone cement, as this will alter cement properties.

**Important physician information**

- As mentioned in the section “Side Effects”, adverse reactions affecting the cardiovascular system have been attributed to acrylic cement. Recent data indicate that the monomer undergoes rapid hydrolysis into methacrylic acid, and that a significant fraction of the circulating methacrylate is present in form of free acid rather than methyl ester. A correlation between changes in circulating concentrations of methyl methacrylate / methacrylic acid and changes in blood pressure has not been established.
- The physician is responsible for any complication or harmful consequences resulting from use of Vertecem V+ Cement Kit in an unapproved indication or in an unapproved manner, or from failure to follow the described operating technique or failure to observe safety instructions in the Instructions for Use.

**Combination of medical devices**

The mixed cement contains a solvent that can soften or etch certain polymeric materials. Therefore only use Vertecem V+ together with DePuy Synthes augmentation hardware with which it has been tested and validated (ex. Vertecem V+ Syringe Kit). For further information, refer to the corresponding device’s or system’s surgical technique. DePuy Synthes has not tested compatibility with devices provided by other manufacturers and assumes no liability in such instances. For information on compatibility with other devices or systems, consultation with a DePuy Synthes representative is recommended.

**Safe disposal**

Mixed Cement Kit bone cement should set before it is discarded in the hospital waste. For separate disposal of powder and liquid refer to local laws and regulations.

### **Magnetic Resonance environment**

Vertecem V+ is inherently electrically nonconductive, nonmetallic, and non-magnetic. Thus the material cannot be heated or act as an antenna when bringing the patient into the magnetic resonance imaging field during examinations. However, Vertecem V+ bone cement is used in combination with metal hardware. Therefore it is important to consider potential hazards related to the use of the metal hardware.

Potential hazards may include but are not limited to:

- Heating or migration of the device
- Artifacts on MR image

### **Operating instructions**

For full details please refer to the Vertecem V+ System surgical technique or that of any other product or system validated for use with the Vertecem V+ Cement Kit. For information on compatibility with other devices or systems consultation with a DePuy Synthes representative is recommended.


During preparation, mixing and transfer make sure to always handle the mixing device by gripping the blue parts. If the transparent part is used as grip the excess body heat provided by the users hand might result in a shortened application time.

1. Open the glass ampoule by breaking off its neck with the plastic cap. Place the opened ampoule in the ampoule holder in the Vertecem V+ Cement Kit inner blister or on a flat, sterile surface.
2. Hold the mixing device upright and make sure the blue handle is in its outmost position. Tap gently on its lid with your finger to ensure that no powder sticks to transportation lid or mixer walls. Remove the transportation lid from the mixing device and dispose of it.
3. Pour the full content of the ampoule into the mixer and close it tightly with the separate mixing and transfer lid. Make sure that both mixing lid and the sealing plug on top of it are securely tightened.
4. Start mixing the Vertecem V+ Cement Kit bone cement by pushing and pulling the blue handle from endpoint to endpoint for 20 seconds (1–2 strokes per second). Perform the first few mixing strokes slowly with an oscillating-rotating movement. Once properly mixed, the blue handle must be left in its outmost position.
5. Once cement has been mixed, remove the sealing plug and connect the mixing device to a compatible application system.


### **Vertecem V+ Cement Kit application time charts**

Note that cement handling and setting times are heavily dependent of temperature. Higher temperatures reduce the setting time and lower temperatures extend it. Handling time also depends on many other factors, including (without necessarily being limited to) syringe diameter, cannula diameter and length.

**Figure 1:** Temperature-time application chart for the Vertecem V+ Cement Kit using a 1 ml syringe and an 8Ga, 150 mm front-opening cannula.

**Figure 2:** Representation of the influence of syringe and cannula geometry on application time.

## Device-related storage and handling information

Temperature limitations: 0°C (32F) – 25°C (77F)

Keep dry

Keep away from sunlight

- Store products in their original protective packaging, and do not remove them from the packaging until immediately before use.
- Use products in the order they were received, i.e. the first-in, first-out principle.

## Additional device-specific information

Does not contain latex

Manufacturer

Highly flammable

Expiration date

Harmful

Do not use when packaging is damaged

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