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Instructions for Use

Please read these instructions for use, any enclosed information and relevant surgical techniques carefully before use. Ensure that you are familiar with the appropriate surgical technique.

Article number

07.702.0405

Article name

TRAUMACEM V+ Bone Cement, injectable

Description

TRAUMACEM V+ Injectable Bone Cement consists of sterile packed components to create radiopaque TRAUMACEM V+ Injectable Bone Cement. It is intended for treatment of osteoporotic fractures of the upper and lower extremities in combination with DePuy Synthes fixation devices intended for PMMA augmentation. The powder component is pre-packed in the mixer. The liquid component is stored in a glass ampoule. A transfer lid for the mixing and transfer of the bone cement in the application system is also included in the kit.

Material(s)

Containers for storage, transport and use:

Container	No. of units	Material(s)
Mixing device	1 unit	Polypropylene high density polyethylene terephthalate (contains powder component of bone cement)
Monomer ampoule	1 unit	Glass (contains liquid component of bone cement)
Mixing and transfer lid	1 unit	Polypropylene UHMW Polyethylene

Bone cement:

Component	Quantity	Constituent(s)
Polymer powder	13 g	<i>Polymethyl methacrylate/acrylate</i> <i>Zirconium dioxide</i> <i>Hydroxyapatite</i> <i>Benzoyl peroxide</i>
	44.6 %	
	40.0 %	
	15.0 %	
Monomer liquid	0.4 %	
	5 ml	<i>Methyl methacrylate (stabilized with 60 ppm Hydroquinone)</i> <i>N,N-dimethyl-p-toluidine</i>
	99.35 %	
0.65 %		

Intended use

The TRAUMACEM V+ Injectable Bone Cement is intended for the treatment of osteoporotic fractures in combination with internal fixation devices intended for PMMA augmentation.

Indications

The TRAUMACEM V+ Injectable Bone Cement is indicated for osteoporotic fractures of the upper and lower extremities in combined procedures with internal fixation devices intended for PMMA augmentation.

When using TRAUMACEM V+ Injectable Bone Cement in combination with DePuy Synthes fixation devices intended for PMMA augmentation, please refer to the surgical technique of the specific device.

Contraindications

- Risk for cement leakage outside the intended application site
- Intracapsular hip fractures
- PMMA bone cement is contraindicated in the presence of active or incompletely treated infection at the site where the bone cement is to be applied
- Patients with severe cardiac and/or pulmonary insufficiency
- Patients with known hypersensitivity or allergy to any of the components of cement (see Composition of Cement)
- Application to children or to women during pregnancy or lactation
- Arthroplasty procedures
- Acute traumatic fractures with good bone quality

Adverse Events

- Serious adverse events, some with fatal outcome, associated with the use of polymethyl methacrylate (PMMA), such as those for vertebroplasty and kyphoplasty, include (but are not necessarily limited to): myocardial infarction, cardiac arrest, cerebrovascular accident, pulmonary and/or cardiac embolism, anaphylaxis, and loosening or displacement of the implant. Although the majority of these adverse events present early within the post-operative period, there have been some reports of diagnoses beyond a year or more after the procedure.
- The most frequent adverse reactions reported in conjunction with the application of PMMA are: transitory drop in blood pressure, thrombophlebitis, hemorrhage and hematoma, superficial or deep wound infection, bursitis, short-term cardiac irregularities and heterotopic bone formation.
- Additional reported adverse events for PMMA bone cements include pyrexia due to an allergy to the bone cement, hematuria, dysuria, bladder fistula, nerve entrapment due to extrusion of the bone cement beyond the region of its intended application, adhesions and stricture of the ileum due to the heat released during polymerization, and transitory worsening of pain due to heat released during the polymerization. Deep wound infection is a serious complication and may require removal of the implant.

Deep wound infection may be latent and not manifest itself for several days or weeks postoperatively.

- Monitor patients carefully for any change in blood pressure during and immediately following the application of PMMA bone cement. Hypotensive reactions have occurred between 10 and 165 seconds following application of bone cement; they have lasted from 30 seconds to 5 or more minutes. Some have progressed to cardiac arrest. Patients should be monitored carefully for any change in blood pressure during and immediately following the application of bone cement.
- Polymethyl methacrylate may cause hypersensitivity among high-risk persons, which can result in an anaphylactic reaction.
- Avoid uncontrolled cement injection, as this may lead to extrusion of the bone cement beyond the site of its intended application and damage the surrounding tissues and may cause patient injury or illness.
- Inadequate fixation or unanticipated postoperative events may affect the cement-bone interface and lead to micromotion of cement against bone surface. A fibrous tissue layer may develop between the cement and the bone, and loosening of the implant may occur leading to implant failure. Long-term follow-up is advised for all patients on a regularly scheduled basis.
- Polymerization of PMMA bone cement is an exothermic reaction that occurs when the cement is hardening in situ. The released heat during this reaction may damage bone or other tissues surrounding the implant.
- The liquid monomer is highly volatile and flammable. The operating room should be adequately ventilated so as to minimize the concentration of monomer vapor. Care should be exercised to prevent exposure to the monomer vapors, which may produce irritation of the respiratory tract and eyes and possibly liver. Exposure to monomer vapors must be avoided, especially during pregnancy. Concentrated vapors of the liquid may have an adverse reaction with soft contact lenses. Personnel wearing contact lenses should not be involved with the mixing of this product.

Sterile device

The TRAUMACEM V+ Injectable Bone Cement is supplied in sterile form. Prior to use, check the product expiration date and verify the integrity of the sterile packaging; if it is opened or otherwise damaged, do not use. If all packaging appears to be intact, remove product(s) from the package using aseptic technique.

The expiry date is printed on the outer carton, outer blister, and patient adhesive labels. TRAUMACEM V+ Injectable Bone Cement must not be used after this date. The polymer powder must not be used if it shows yellow discoloration. The liquid monomer, the ampoule, the cement powder, the mixer, the transfer lid, the inner blister and the inside of the outer blister are all sterile.

The monomer fluid inside the ampoule is sterilized using an aseptic processing technique (ultra filtration).

The double blister and all its content is sterilized using ethylene oxide.

Do not resterilize

Single-use device

Do not re-use

The TRAUMACEM V+ Injectable Bone Cement is intended for single-use only and must not be re-used. Re-use or reprocessing (e.g. cleaning and re-sterilization) may compromise the structural integrity of the device and/or lead to device failure which may result in patient injury, illness, or death.

Reuse of single-use devices may create a risk of contamination. This could result in injury or death of the patient or user.

A single TRAUMACEM V+ Injectable Bone Cement unit may not be used to treat more than one patient.

Precautions

- Read instructions carefully prior to use. Always check the condition of the liquid before performing the procedure. Do not use the liquid component if the ampoule shows any sign of weakening or premature polymerization. Do not use the product after the expiry date. Do not use TRAUMACEM V+ Injectable Bone Cement if the center of the temperature indicator located on the packaging is black.
- If stored below OR temperature, always allow the product to equilibrate to OR temperature prior to use.
- Follow handling, mixing and preparation instructions carefully.
- Liquid methyl methacrylate is a powerful lipid solvent; it should not be allowed to come into direct contact with sensitive tissue. The liquid component should not be allowed to come into contact with rubber or latex gloves. Wearing of a second pair of gloves and strict adherence to the mixing instructions may diminish the risk of hypersensitivity reactions.
- Always use the full amounts of monomer liquid and polymer powder provided in the kit, respectively, when mixing TRAUMACEM V+ Injectable Bone Cement. Ensure that the powder and liquid component are thoroughly mixed before starting cement transfer.
- The mixing system is only appropriate for transferring mixed TRAUMACEM V+ Injectable Bone Cement into an adequate application system. The mixing system is not permitted for direct application of the bone cement.
- The formation of large air pockets within the cement should be avoided when filling the application system.

Instructions for Use

- Application of TRAUMACEM V+ Injectable Bone Cement should be monitored using real time imaging procedures capable of delivering high quality images. Use appropriate imaging techniques to confirm correct needle placement, the absence of damage to surrounding structures and appropriate location and quantity of injected material.
- At the completion of the cement injection process, the patient should remain immobile for 15 minutes to facilitate proper cement curing. Maintain the patient's position until the polymerization process of the cement is completed. Cement at body temperature (37°C) hardens faster than at room temperature.
- Deep wound infection is a serious complication and may require removal of the implant. Deep wound infection may be latent and not manifest itself for several days or weeks postoperatively.
- For cancer patients a careful logical analysis is needed prior to surgery to limit risks to the patient and optimize chances for a favorable outcome.
- Leaks can occur through fracture cracks.
- Prior to cement injection use a contrast medium to determine possible leaks so as not to damage the articular surface due to uncontrolled cement injection. Ensure to select an appropriate injection site and correct quantity of injected material.
- If TRAUMACEM V+ Injectable Bone Cement is found outside the bony tissue, injection of the material must be stopped immediately.

Warnings

- Extravasation of PMMA bone cement can result in cement entering fracture cracks or the affected joint, but is not necessarily limited to said areas if prescribed procedure is not followed.
- Additives (such as antibiotics) are not to be mixed with the TRAUMACEM V+ Injectable Bone Cement, as this will alter cement properties.
- Do not augment if X-ray contrast media leaks into the joint.
- In the event that there is a danger of cement leakage into the joint or fracture gap, stop injection immediately.
- In the event of PMMA bone cement leakage into an unintended area tissue damage, nerve or circulatory problems, and other serious adverse events may occur.
- Caution should be taken prior to and during implantation in the presence of tumor(s) at the augmentation site.

Important information for physician

- Adverse reactions affecting the cardiovascular system have been associated with use of PMMA and acrylic cement. Recent data indicate that the monomer undergoes rapid hydrolysis into methacrylic acid, and that a significant fraction of the circulating methacrylate is present in form of free acid rather than methylester. A correlation between changes in circulating concentrations of methyl methacrylate/methacrylic acid and changes in blood pressure has not been established.

- The physician is responsible for any complication or harmful consequences resulting from use of TRAUMACEM V+ Injectable Bone Cement in an unapproved indication or in an unapproved manner, or from failure to follow the described operating technique or failure to observe safety instructions in the instructions for use.
- The physician should, by specific training and experience, be thoroughly familiar with the properties, handling characteristics, and application of bone cements. Because the handling and curing characteristics of TRAUMACEM V+ Injectable Bone Cement vary with temperature, humidity, and mixing technique, they are best determined by the surgeon's actual experience. The physician is strongly recommended to carry out a trial run of the entire mixing, handling, and setting process prior to performing a surgical procedure with TRAUMACEM V+ Injectable Bone Cement.
- TRAUMACEM V+ Injectable Bone Cement must only be used by qualified physicians, with knowledge and experience of the specific operating techniques for this product. A thorough preoperative evaluation of the patient must be carried out prior to the operation. The surgeon should have specific training and experience in performing and augmentation procedure to be thoroughly familiar with the properties, handling characteristics and application of the product with percutaneous cement delivery.

Combination of medical devices

Only use TRAUMACEM V+ Injectable Bone Cement together with DePuy Synthes augmentation hardware with which it has been tested and validated (ex. TRAUMACEM V+ Syringe Kit – Sterile). The mixed cement contains a solvent that can soften or etch certain polymeric materials. For further information, refer to corresponding device's or system's surgical technique. DePuy Synthes has not tested compatibility with devices provided by other manufacturers and assumes no liability in such instances.

Magnetic Resonance environment

The medical device TRAUMACEM V+ Injectable Bone Cement is MR Safe. However, TRAUMACEM V+ Injectable Bone Cement is only used in combination with metal hardware; refer to the corresponding device's labeling for additional instructions or information essential to safe use in the MR environment.

Disposal

Mixed TRAUMACEM V+ Injectable Bone Cement should set before it is discarded in the hospital waste. For separate disposal of powder and liquid refer to local laws and regulations.

Special Operating Instructions

The blister pack containing all components should be removed from the carton in the non-sterile area. The portion of the outerblister which is to be peeled back is not sterile and should be opened by the circulating nurse. The

contained inner blister is passed to a member of the surgical team in the sterile area, using aseptic technique. Mix the TRAUMACEM V+ Injectable Bone Cement by using the mixing system and follow the instructions given. For mixing TRAUMACEM V+ Injectable Bone Cement a sterile working area is needed.

During preparation, mixing and transfer make sure to always handle the mixing device by gripping the blue parts. If the transparent part is used as grip the excess body heat provided by the users hand might result in a shortened application time.

1. Pull the handle of the mixer in the fully back position.
2. Open the glass ampoule by breaking the bottle neck just before steps 3–5.
3. Hold the mixing device upright by the junction between the cartridge and the blue part. Gently tap the top of the cement mixer three times in order to ensure no cement powder sticks to the top of the mixer cartridge and the mixer sterilization lid.
4. Remove the sterilization lid of the mixer and dispose it. Take care not to lean the handle on the table in order to avoid that the powder is spilled.
5. Pour all the monomer from the glass ampoule into the cement powder and close the mixer with the enclosed transfer lid tightly.
6. Hold the mixer by the blue part, mix TRAUMACEM V+ Injectable Bone Cement by moving the blue plunger back and forth from stop to stop for about 20 seconds (1–2 cycles/second). Perform the first mixing strokes slowly with oscillating-rotating movements.
7. When mixing is completed, pull the plunger in the fully back position, open the small translucent plug of the mixer and connect the cement mixer to the application system. The Transfer of the cement into the application system should be done immediately after mixing is complete.

For details please follow the instructions for use and the surgical technique of the device which is used for internal fixation.

TRAUMACEM V+ Injectable Bone Cement application time charts

Note that cement handling and setting times are heavily dependent on temperature. Higher temperatures reduce the setting time and lower temperatures extend it. At cement temperatures lower than 19°C (66°F), a longer waiting time is needed to reach an appropriate viscosity of the cement. Handling time also depends on many other factors, including (without necessarily being limited to) syringe diameter, cannula diameter and length.


Figure 1: Temperature-time application chart for the TRAUMACEM V+ Injectible Bone Cement using a 1 ml syringe and an 8Ga, 150 mm front-opening cannula.


Figure 2: Representation of the influence of syringe and cannula geometry on application time.

Device-related storage and handling information

Temperature limitations: 0°C (32°F) – 25°C (77°F)



Keep dry

Keep away from sunlight

- Products should be stored in a dry and clean environment and protected from direct sunlight.
- Store products in their original protective packaging, and do not remove them from the packaging until immediately before use.
- Use products in the order they were received, i.e. the first-in, first-out principle.
- The cement is designed for handling between 19°C (66°F) and 25°C (77°F) as standard. Before using TRAUMACEM V+ Injectable Bone Cement allow the product to reach the necessary temperature described in the time-temperature graph.
- Do not use TRAUMACEM V+ Injectable Bone Cement if the center of the temperature indicator located on the packaging is black.

Use	 <p>REF 07.702.040S TRAUMACEM™ V+ Bone Cement, injectable</p> <p>File: 130-0868-03</p>
Do not use	 <p>REF 07.702.040S TRAUMACEM™ V+ Bone Cement, injectable</p> <p>File: 130-0868-03</p>

Additional device-specific information

 LOT	Lot or batch number		Harmful
	Expiration date		Does not contain latex
 REF	Reference number	 0123	Notified body
	Manufacturer		Do not use if packaging is damaged
	Highly flammable		Caution, see instructions for use

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