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Instruction for Use

Please read these instructions for use, any enclosed information and relevant surgical techniques carefully before use. Ensure that you are familiar with the appropriate surgical technique.

The "TRAUMACEM™ V+ Syringe Kit, 4×1 mL, 2.3 mm Adapter" is a sterile packed, single use syringe kit used for the injection of DePuy Synthes PMMA-based cement into cancellous bone.

Since the "TRAUMACEM V+ Syringe Kit, 4×1 mL, 2.3 mm Adapter" is used in technically demanding procedures, the surgeon should be familiar with the use of the syringe kit as well as with the rest of the instrumentation and the procedure.

Product name and materials

"TRAUMACEM V+ Syringe Kit, 4×1 mL, 2.3 mm Adapter". One kit contains the following parts:

Product name	Number of units	Materials
1 mL syringe, blue plunger	4 units	PA12 (Polyamide) ABS (Acrylonitrile Butadiene Styrene) Silicone
Adapters	6 units	PEEK (Polyether ether ketone)
Stopcock	1 unit	PA 12 (Polyamide) POM (Polyoxymethylen)

Intended use

The "TRAUMACEM V+ Syringe Kit, 4x1 mL, 2.3 mm Adapter" is intended for the application of TRAUMACEM V+ PMMA-based bone cement with the intent to augment cancellous bone.

The product must only be used with DePuy Synthes-approved products. Refer to the corresponding directions regarding indications, contraindications, compatibility, use, precautions, warnings and side effects of the TRAUMACEM V+ PMMA-based bone cement and access solution used in conjunction with the "TRAUMACEM V+ Syringe Kit, 4x1 mL, 2.3 mm Adapter".

The 2.3 mm adapter is also intended for the application of a radiographic contrast agent/contrast agent. Refer to the corresponding directions regarding indications, contraindications, use, precautions, warnings and side effects of the radiographic contrast agent/contrast agent.

For information on compatibility with other devices or systems, consultation with a DePuy Synthes representative is recommended.

Side effects

There are no known side effects or interactions associated with the use of the "TRAUMACEM V+ Syringe Kit, 4x1 mL, 2.3 mm Adapter" together with DePuy Synthes PMMA-based bone cement.

Sterile device

Sterilized using ethylene oxide

Store instruments in their original protective packaging, and do not remove them from the packaging until immediately before use.

Prior to use, check the product expiration date and verify the integrity of the sterile packaging. If all packaging appear to be intact, remove product(s) from the package using aseptic technique. Do not use if the package is damaged or expiration date is passed.

Do not resterilize

Single-use device

Do not reuse

The "TRAUMACEM V+ Syringe Kit, 4x1 mL, 2.3 mm Adapter" is intended for single use only and must not be reused. Reuse or reprocessing (e.g. cleaning and/or resterilization) may compromise the structural integrity of the device and/or lead to device failure which may result in patient injury, illness or death. Furthermore, reuse or reprocessing of single-use devices may create a risk of contamination, e.g. due to the transmission of infectious material from one patient or person to another. This could result in injury or death of the patient or user.

Device-specific warnings

The purpose of the milliliter (mL) scale on the side of the "TRAUMACEM V+ Syringe Kit, 4x1 mL, 2.3 mm Adapter" syringes is to give physicians an indication of how much PMMA-based bone cement has been injected. This scale is for information purposes only.

Precautions

- The "TRAUMACEM V+ Syringe Kit, 4x1 mL, 2.3 mm Adapter" is used in technically demanding procedures. Therefore only physicians familiar with the proper technique and instrumentation used for the delivery and use of PMMA-based bone cement should use the "TRAUMACEM V+ Syringe Kit, 4x1 mL, 2.3 mm Adapter".
- The "TRAUMACEM V+ Syringe Kit, 4x1 mL, 2.3 mm Adapter" should only be used to deliver PMMA-based bone cement in conjunction with the use of fluoroscopy imaging equipment capable of delivering high-quality images.
- Ensure a good fit between the "TRAUMACEM V+ Syringe Kit, 4x1 mL, 2.3 mm Adapter" syringe and the "TRAUMACEM V+ Syringe Kit, 4x1 mL, 2.3 mm Adapter" stop-cock/used access solution, but make sure to be on axis and avoid using excessive force when coupling them. They are both made of plastic and could otherwise break.
- After filling a "TRAUMACEM V+ Syringe Kit, 4x1 mL, 2.3 mm Adapter" syringe with cement, grip it by its wings. Gripping it anywhere else could transfer heat from your hand to the cement, shortening the working time of the PMMA-based bone cement.
- Always fill all "TRAUMACEM V+ Syringe Kit, 4x1 mL, 2.3 mm Adapter" syringes with PMMA-based bone cement directly after mixing.
- Always make sure the PMMA-based bone cement has reached the preferred viscosity prior to injection.
- To halt the flow of cement out of the "TRAUMACEM V+ Syringe Kit, 4x1 mL, 2.3 mm Adapter" syringe at any given time, simply stop pushing the plunger.
- Make sure to always inject cement in a slow, controlled fashion.
- If cement flow is hindered at any time during injection, halt, investigate and correct the reason for the flow hindrance.
- Do not use tools or excessive force to operate the "TRAUMACEM V+ Syringe Kit, 4x1 mL, 2.3 mm Adapter".

- The syringes in the “TRAUMACEM V+ Syringe Kit, 4×1 mL, 2.3 mm Adapter” are solely intended to be used for PMMA-based bone cement augmentation and must not be used with radiographic contrast agent/contrast agent.
- For injection of radiographic contrast agent/contrast agent, only syringes with luer lock with an appropriate intended use should be used.
- Always consult the radiographic contrast agent/contrast agent instruction for use and the warnings therein prior to administration to ensure that the used radiographic contrast agent/contrast agent is approved for this application.

Combination of medical devices

The “TRAUMACEM V+ Syringe Kit, 4×1 mL, 2.3 mm Adapter” is tested and validated with the DePuy Synthes PMMA-based bone cement and DePuy Synthes augmentation hardware only.

DePuy Synthes has not tested compatibility with devices provided by other manufacturers and assumes no liability in such instances.

For information on compatibility with other devices or systems, consultation with a DePuy Synthes representative is recommended. Refer to corresponding device’s or system’s directions for use.

Operating instructions

1. Once cement has been mixed connect the stopcock to the mixer. Attach the side without funnel of the stopcock to the mixer. The handle in the initial position is turned 90° away from the mixer and the “off” sign is on the opposite side of the funnel. Ensure a tight fitting between stopcock and mixer but avoid breakage of the stopcock due to excessive torque when applying.

2. First the air has to be removed from the system. Gently turn the handle of the cement mixer clockwise. You will see the piston of the mixer advancing in the translucent cartridge and after a while cement will flow into the stopcock. As soon as the cement is visible in the stopcock close the stopcock by turning the handle ("off") toward the mixer (90°).
3. Attach a syringe to the stopcock (funnel side).
4. Open the stopcock by turning the handle (90° turn), back to its original position.
5. Carefully fill the syringe by slowly turning the handle of the cement mixer clockwise.
6. When the syringe is filled turn the handle of the stopcock again (90°) towards the mixer. The "off" sign is directed toward the mixer stopping the cement flow.
7. Disconnect the filled syringe and attach the next empty syringe. Repeat steps 3 to 7 for all syringes. Continue to fill the syringes in the same manner. Always fill all syringes at once.
8. To access adapters remove the cover lid. Screw the adapters on the filled syringes.
9. Follow the instructions for use of the bone cement to determine when to inject the bone cement best.
10. Attach the syringe with the adapter to the recess of the canulated screw which needs to be augmented – ensure a good interlock between recess and adapter.
11. Carefully inject the bone cement under the use of adequate imaging technique.

Troubleshooting

If the force needed to inject cement becomes too high, this indicates that the maximum working time for the combination of DePuy Synthes PMMA-based bone cement and access solution used has been reached.

Device-related storage and handling information

Keep dry

Keep away from sunlight

Additional device-specific information

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© Synthes GmbH 2017. All rights reserved. 60105404 SE_501031 AG 09/2017
Synthes GmbH