

---

# Instructions for Use

## TFN-ADVANCED Proximal Femoral Nailing System

This instruction for use is not intended for distribution in the USA.

Not all products are currently available in all markets.



### **Authorised Representative**

DePuy Ireland UC  
Loughbeg  
Ringaskiddy  
Co. Cork Ireland

# Instructions for Use

## TFNA, Ø 9 mm, long

| Right       | Left        | Length (mm) | Femoral Neck Angle |
|-------------|-------------|-------------|--------------------|
| 04.037.916S | 04.037.917S | 260         | 125°               |
| 04.037.918S | 04.037.919S | 280         | 125°               |
| 04.037.920S | 04.037.921S | 300         | 125°               |
| 04.037.922S | 04.037.923S | 320         | 125°               |
| 04.037.924S | 04.037.925S | 340         | 125°               |
| 04.037.926S | 04.037.927S | 360         | 125°               |
| 04.037.928S | 04.037.929S | 380         | 125°               |
| 04.037.930S | 04.037.931S | 400         | 125°               |
| 04.037.932S | 04.037.933S | 420         | 125°               |
| 04.037.934S | 04.037.935S | 440         | 125°               |
| 04.037.936S | 04.037.937S | 460         | 125°               |
| 04.037.938S | 04.037.939S | 480         | 125°               |
| 04.037.946S | 04.037.947S | 260         | 130°               |
| 04.037.948S | 04.037.949S | 280         | 130°               |
| 04.037.950S | 04.037.951S | 300         | 130°               |
| 04.037.952S | 04.037.953S | 320         | 130°               |
| 04.037.954S | 04.037.955S | 340         | 130°               |
| 04.037.956S | 04.037.957S | 360         | 130°               |
| 04.037.958S | 04.037.959S | 380         | 130°               |
| 04.037.960S | 04.037.961S | 400         | 130°               |
| 04.037.962S | 04.037.963S | 420         | 130°               |
| 04.037.964S | 04.037.965S | 440         | 130°               |
| 04.037.966S | 04.037.967S | 460         | 130°               |
| 04.037.968S | 04.037.969S | 480         | 130°               |

## TFNA, Ø 10 mm, long

| Right       | Left        | Length (mm) | Femoral Neck Angle |
|-------------|-------------|-------------|--------------------|
| 04.037.016S | 04.037.017S | 260         | 125°               |
| 04.037.018S | 04.037.019S | 280         | 125°               |
| 04.037.020S | 04.037.021S | 300         | 125°               |
| 04.037.022S | 04.037.023S | 320         | 125°               |
| 04.037.024S | 04.037.025S | 340         | 125°               |
| 04.037.026S | 04.037.027S | 360         | 125°               |
| 04.037.028S | 04.037.029S | 380         | 125°               |
| 04.037.030S | 04.037.031S | 400         | 125°               |
| 04.037.032S | 04.037.033S | 420         | 125°               |
| 04.037.034S | 04.037.035S | 440         | 125°               |
| 04.037.036S | 04.037.037S | 460         | 125°               |
| 04.037.038S | 04.037.039S | 480         | 125°               |
| 04.037.046S | 04.037.047S | 260         | 130°               |
| 04.037.048S | 04.037.049S | 280         | 130°               |
| 04.037.050S | 04.037.051S | 300         | 130°               |
| 04.037.052S | 04.037.053S | 320         | 130°               |
| 04.037.054S | 04.037.055S | 340         | 130°               |
| 04.037.056S | 04.037.057S | 360         | 130°               |
| 04.037.058S | 04.037.059S | 380         | 130°               |
| 04.037.060S | 04.037.061S | 400         | 130°               |
| 04.037.062S | 04.037.063S | 420         | 130°               |
| 04.037.064S | 04.037.065S | 440         | 130°               |
| 04.037.066S | 04.037.067S | 460         | 130°               |
| 04.037.068S | 04.037.069S | 480         | 130°               |

## TFNA, Ø 11 mm, long

| Right       | Left        | Length (mm) | Femoral Neck Angle |
|-------------|-------------|-------------|--------------------|
| 04.037.120S | 04.037.121S | 300         | 125°               |
| 04.037.122S | 04.037.123S | 320         | 125°               |
| 04.037.124S | 04.037.125S | 340         | 125°               |
| 04.037.126S | 04.037.127S | 360         | 125°               |
| 04.037.128S | 04.037.129S | 380         | 125°               |
| 04.037.130S | 04.037.131S | 400         | 125°               |
| 04.037.132S | 04.037.133S | 420         | 125°               |
| 04.037.134S | 04.037.135S | 440         | 125°               |
| 04.037.136S | 04.037.137S | 460         | 125°               |
| 04.037.138S | 04.037.139S | 480         | 125°               |
| 04.037.150S | 04.037.151S | 300         | 130°               |
| 04.037.152S | 04.037.153S | 320         | 130°               |
| 04.037.154S | 04.037.155S | 340         | 130°               |
| 04.037.156S | 04.037.157S | 360         | 130°               |
| 04.037.158S | 04.037.159S | 380         | 130°               |
| 04.037.160S | 04.037.161S | 400         | 130°               |
| 04.037.162S | 04.037.163S | 420         | 130°               |
| 04.037.164S | 04.037.165S | 440         | 130°               |
| 04.037.166S | 04.037.167S | 460         | 130°               |
| 04.037.168S | 04.037.169S | 480         | 130°               |
| 04.037.180S | 04.037.181S | 300         | 135°               |
| 04.037.182S | 04.037.183S | 320         | 135°               |
| 04.037.184S | 04.037.185S | 340         | 135°               |
| 04.037.186S | 04.037.187S | 360         | 135°               |
| 04.037.188S | 04.037.189S | 380         | 135°               |
| 04.037.190S | 04.037.191S | 400         | 135°               |
| 04.037.192S | 04.037.193S | 420         | 135°               |
| 04.037.194S | 04.037.195S | 440         | 135°               |
| 04.037.196S | 04.037.197S | 460         | 135°               |
| 04.037.198S | 04.037.199S | 480         | 135°               |

## TFNA, Ø 12 mm, long

| Right       | Left        | Length (mm) | Femoral Neck Angle |
|-------------|-------------|-------------|--------------------|
| 04.037.220S | 04.037.221S | 300         | 125°               |
| 04.037.222S | 04.037.223S | 320         | 125°               |
| 04.037.224S | 04.037.225S | 340         | 125°               |
| 04.037.226S | 04.037.227S | 360         | 125°               |
| 04.037.228S | 04.037.229S | 380         | 125°               |
| 04.037.230S | 04.037.231S | 400         | 125°               |
| 04.037.232S | 04.037.233S | 420         | 125°               |
| 04.037.234S | 04.037.235S | 440         | 125°               |
| 04.037.236S | 04.037.237S | 460         | 125°               |
| 04.037.238S | 04.037.239S | 480         | 125°               |
| 04.037.250S | 04.037.251S | 300         | 130°               |
| 04.037.252S | 04.037.253S | 320         | 130°               |
| 04.037.254S | 04.037.255S | 340         | 130°               |
| 04.037.256S | 04.037.257S | 360         | 130°               |
| 04.037.258S | 04.037.259S | 380         | 130°               |
| 04.037.260S | 04.037.261S | 400         | 130°               |
| 04.037.262S | 04.037.263S | 420         | 130°               |
| 04.037.264S | 04.037.265S | 440         | 130°               |
| 04.037.266S | 04.037.267S | 460         | 130°               |
| 04.037.268S | 04.037.269S | 480         | 130°               |

**TFNA, Ø 14 mm, long**

| Right       | Left        | Length (mm) | Femoral Neck Angle |
|-------------|-------------|-------------|--------------------|
| 04.037.450S | 04.037.451S | 300         | 130°               |
| 04.037.452S | 04.037.453S | 320         | 130°               |
| 04.037.454S | 04.037.455S | 340         | 130°               |
| 04.037.456S | 04.037.457S | 360         | 130°               |
| 04.037.458S | 04.037.459S | 380         | 130°               |
| 04.037.460S | 04.037.461S | 400         | 130°               |
| 04.037.462S | 04.037.463S | 420         | 130°               |
| 04.037.464S | 04.037.465S | 440         | 130°               |
| 04.037.466S | 04.037.467S | 460         | 130°               |
| 04.037.468S | 04.037.469S | 480         | 130°               |

**TFNA, short, length 170 mm**

| Right       | Dia. (mm) | Femoral Neck Angle |
|-------------|-----------|--------------------|
| 04.037.912S | 9         | 125°               |
| 04.037.942S | 9         | 130°               |
| 04.037.972S | 9         | 135°               |
| 04.037.012S | 10        | 125°               |
| 04.037.042S | 10        | 130°               |
| 04.037.072S | 10        | 135°               |
| 04.037.112S | 11        | 125°               |
| 04.037.142S | 11        | 130°               |
| 04.037.172S | 11        | 135°               |
| 04.037.212S | 12        | 125°               |
| 04.037.242S | 12        | 130°               |
| 04.037.272S | 12        | 135°               |

**TFNA, short, length 200 mm**

| Right       | Dia. (mm) | Femoral Neck Angle |
|-------------|-----------|--------------------|
| 04.037.913S | 9         | 125°               |
| 04.037.943S | 9         | 130°               |
| 04.037.973S | 9         | 135°               |
| 04.037.013S | 10        | 125°               |
| 04.037.043S | 10        | 130°               |
| 04.037.073S | 10        | 135°               |
| 04.037.113S | 11        | 125°               |
| 04.037.143S | 11        | 130°               |
| 04.037.173S | 11        | 135°               |
| 04.037.213S | 12        | 125°               |
| 04.037.243S | 12        | 130°               |
| 04.037.273S | 12        | 135°               |

**TFNA, short, length 235 mm**

| Right       | Left        | Length (mm) | Femoral Neck Angle |
|-------------|-------------|-------------|--------------------|
| 04.037.914S | 04.037.915S | 9           | 125°               |
| 04.037.944S | 04.037.945S | 9           | 130°               |
| 04.037.974S | 04.037.975S | 9           | 135°               |
| 04.037.014S | 04.037.015S | 10          | 125°               |
| 04.037.044S | 04.037.045S | 10          | 130°               |
| 04.037.074S | 04.037.075S | 10          | 135°               |
| 04.037.114S | 04.037.115S | 11          | 125°               |
| 04.037.144S | 04.037.145S | 11          | 130°               |
| 04.037.174S | 04.037.175S | 11          | 135°               |
| 04.037.214S | 04.037.215S | 12          | 125°               |
| 04.037.244S | 04.037.245S | 12          | 130°               |
| 04.037.274S | 04.037.275S | 12          | 135°               |

**TFNA Screws\***

|            | Length (mm) |            | Length (mm) |
|------------|-------------|------------|-------------|
| 04.038.070 | 70          | 04.038.105 | 105         |
| 04.038.075 | 75          | 04.038.110 | 110         |
| 04.038.080 | 80          | 04.038.115 | 115         |
| 04.038.085 | 85          | 04.038.120 | 120         |
| 04.038.090 | 90          | 04.038.125 | 125         |
| 04.038.095 | 95          | 04.038.130 | 130         |
| 04.038.100 | 100         |            |             |

**TFNA Helical Blades\***

|            | Length (mm) |            | Length (mm) |
|------------|-------------|------------|-------------|
| 04.038.270 | 70          | 04.038.305 | 105         |
| 04.038.275 | 75          | 04.038.310 | 110         |
| 04.038.280 | 80          | 04.038.315 | 115         |
| 04.038.285 | 85          | 04.038.320 | 120         |
| 04.038.290 | 90          | 04.038.325 | 125         |
| 04.038.295 | 95          | 04.038.330 | 130         |
| 04.038.300 | 100         |            |             |

**TFNA Helical Blades, perforated**

|             | Length (mm) |             | Length (mm) |
|-------------|-------------|-------------|-------------|
| 04.038.370S | 70          | 04.038.405S | 105         |
| 04.038.375S | 75          | 04.038.410S | 110         |
| 04.038.380S | 80          | 04.038.415S | 115         |
| 04.038.385S | 85          | 04.038.420S | 120         |
| 04.038.390S | 90          | 04.038.425S | 125         |
| 04.038.395S | 95          | 04.038.430S | 130         |
| 04.038.400S | 100         |             |             |

**TFNA Screws, perforated**

|             | Length (mm) |             | Length (mm) |
|-------------|-------------|-------------|-------------|
| 04.038.170S | 70          | 04.038.205S | 105         |
| 04.038.175S | 75          | 04.038.210S | 110         |
| 04.038.180S | 80          | 04.038.215S | 115         |
| 04.038.185S | 85          | 04.038.220S | 120         |
| 04.038.190S | 90          | 04.038.225S | 125         |
| 04.038.195S | 95          | 04.038.230S | 130         |
| 04.038.200S | 100         |             |             |

**Locking Screw for Medullary Nails, Ø 5 mm\***

|            | Length (mm) |            | Length (mm) |
|------------|-------------|------------|-------------|
| 04.045.026 | 26          | 04.045.066 | 66          |
| 04.045.028 | 28          | 04.045.068 | 68          |
| 04.045.030 | 30          | 04.045.070 | 70          |
| 04.045.032 | 32          | 04.045.072 | 72          |
| 04.045.034 | 34          | 04.045.074 | 74          |
| 04.045.036 | 36          | 04.045.076 | 76          |
| 04.045.038 | 38          | 04.045.078 | 78          |
| 04.045.040 | 40          | 04.045.080 | 80          |
| 04.045.042 | 42          | 04.045.082 | 82          |
| 04.045.044 | 44          | 04.045.084 | 84          |
| 04.045.046 | 46          | 04.045.086 | 86          |
| 04.045.048 | 48          | 04.045.088 | 88          |
| 04.045.050 | 50          | 04.045.090 | 90          |
| 04.045.052 | 52          | 04.045.095 | 95          |
| 04.045.054 | 54          | 04.045.100 | 100         |
| 04.045.056 | 56          | 04.045.105 | 105         |
| 04.045.058 | 58          | 04.045.110 | 110         |
| 04.045.060 | 60          | 04.045.115 | 115         |
| 04.045.062 | 62          | 04.045.120 | 120         |
| 04.045.064 | 64          |            |             |

**Locking Screw for Medullary Nails, Low Profile, Ø 5 mm\***

|            | Length (mm) |            | Length (mm) |
|------------|-------------|------------|-------------|
| 04.045.326 | 26          | 04.045.366 | 66          |
| 04.045.328 | 28          | 04.045.368 | 68          |
| 04.045.330 | 30          | 04.045.370 | 70          |
| 04.045.332 | 32          | 04.045.372 | 72          |
| 04.045.334 | 34          | 04.045.374 | 74          |
| 04.045.336 | 36          | 04.045.376 | 76          |
| 04.045.338 | 38          | 04.045.378 | 78          |
| 04.045.340 | 40          | 04.045.380 | 80          |
| 04.045.342 | 42          | 04.045.382 | 82          |
| 04.045.344 | 44          | 04.045.384 | 84          |
| 04.045.346 | 46          | 04.045.386 | 86          |
| 04.045.348 | 48          | 04.045.388 | 88          |
| 04.045.350 | 50          | 04.045.390 | 90          |
| 04.045.352 | 52          | 04.045.395 | 95          |
| 04.045.354 | 54          | 04.045.400 | 100         |
| 04.045.356 | 56          | 04.045.405 | 105         |
| 04.045.358 | 58          | 04.045.410 | 110         |
| 04.045.360 | 60          | 04.045.415 | 115         |
| 04.045.362 | 62          | 04.045.420 | 120         |
| 04.045.364 | 64          |            |             |

**End Caps**

|             | Length (mm) |
|-------------|-------------|
| 04.045.870S | 0           |
| 04.045.875S | 5           |
| 04.045.880S | 10          |
| 04.045.885S | 15          |

**Nut and Washers**

|             |                    |
|-------------|--------------------|
| 04.045.780S | Washer Ø 14/7      |
| 04.045.781S | Nut Ø 14           |
| 04.045.782S | Washer Ø 17.5/11.8 |

Alternatively, the TFNA implants can be applied using associated instrumentation and a set of the following compatible screw implants:

**Locking Screw Stardrive® Ø 5 mm\***

|            | Length (mm) |            | Length (mm) |
|------------|-------------|------------|-------------|
| 04.005.516 | 26          | 04.005.548 | 58          |
| 04.005.518 | 28          | 04.005.550 | 60          |
| 04.005.520 | 30          | 04.005.552 | 62          |
| 04.005.522 | 32          | 04.005.554 | 64          |
| 04.005.524 | 34          | 04.005.556 | 66          |
| 04.005.526 | 36          | 04.005.558 | 68          |
| 04.005.528 | 38          | 04.005.560 | 70          |
| 04.005.530 | 40          | 04.005.562 | 72          |
| 04.005.532 | 42          | 04.005.564 | 74          |
| 04.005.534 | 44          | 04.005.566 | 76          |
| 04.005.536 | 46          | 04.005.568 | 78          |
| 04.005.538 | 48          | 04.005.570 | 80          |
| 04.005.540 | 50          | 04.005.575 | 85          |
| 04.005.542 | 52          | 04.005.580 | 90          |
| 04.005.544 | 54          | 04.005.585 | 95          |
| 04.005.546 | 56          | 04.005.590 | 100         |

**End Caps**

|             | Length (mm) |
|-------------|-------------|
| 04.038.000S | 0           |
| 04.038.005S | 5           |
| 04.038.010S | 10          |
| 04.038.015S | 15          |

\* Available non-sterile or sterile packed. Add "S" to the catalogue number to order sterile products."

Products available non-sterile and sterile can be differentiated with the suffix "S" added to the article number for sterile products.

The screws are also available in sterile tube packaging (corresponding article number with suffix "TS").

Screw length designations are defined to reflect the readings on the length measurement tools and do not necessarily correspond to the actual total length of the screw.

**Introduction**

The TFN-ADVANCED™ Proximal Femoral Nailing System (TFNA) consists of cannulated femoral nails, helical blades or screws, end caps and locking screws. The TFNA Nail is anatomically contoured and tapers to a nominal diameter of 9, 10, 11, 12, or 14 mm. The proximal locking hole accommodates angles ranging from 125°–135°. TFNA Nails are available in short lengths (170–235 mm) and long nail lengths (260–480 mm), with the lengths 235 mm and above available in right or left versions. The TFNA accepts commercially available Synthes 4.9mm Locking Bolts and/or 5.0mm Locking Screws. This system is manufactured from titanium alloy and are provided in sterile and non-sterile packaging. TFNA also has the option for cement augmentation of the TFNA Blade (perforated or non-perforated) and TFNA Screw (perforated or non-perforated).

Important note for medical professionals and operating room staff: These instructions for use do not include all the information necessary for selection and use of a device. Please read the instructions for use and the Synthes brochure "Important Information" (SE\_023827) carefully before use. Ensure that you are familiar with the appropriate surgical procedure.

**Materials**

| Device(s)   | Material(s)                        | Standard(s)               |
|---|------------------------------------|---------------------------|
| Nails   | Ti-15Mo (TiMo)<br>Titanium Alloy   | ASTM F2066                |
| Nails (Locking Mechanism), End Caps, Head Elements (helical blades and screws), Locking Screws, Nut | Ti-6Al-7Nb (TAN)<br>Titanium Alloy | ISO 5832-11<br>ASTM F1295 |
| Nails (Locking Mechanism)   | 40Co-20Cr-16Fe-15Ni-7Mo (Elgiloy)  | ISO 5832-7<br>ASTM F1058  |
| Washers   | TiCP                               | ISO 5832-2<br>ASTM F67    |

### Intended Use

Proximal Femoral Nailing Implants – including the TFNA implants – are intended to be used for temporary fixation and stabilization of the proximal femur and the femoral shaft.

### Indications

TFNA SHORT (lengths 170 mm, 200 mm, 235 mm)

- Pertrochanteric fractures (31-A1 and 31-A2)
- Intertrochanteric fractures (31-A3)
- 235 mm nails are additionally indicated for high subtrochanteric fractures

TFNA LONG (lengths 260 mm–480 mm)

- Pertrochanteric fractures (31-A1 and 31-A2)
- Intertrochanteric fractures (31-A3)
- Fractures of the trochanteric area (31-A1/A2/A3) with diaphyseal extension
- Combined fractures of the trochanteric area (31-A1/A2/A3) and the femoral shaft (32-A/B/C)
- Pathological fractures, including prophylactic use
- Malunions
- Non-unions

TFNA AUGMENTATION

- For fractures in the proximal femur with poor bone quality and/or increased risk of fixation failure at the implant/bone interface.

For indications and contraindications of the "TRAUMACEM™ V+ Injectable Bone Cement", the "TRAUMACEM V+ Syringe Kit" and the "TRAUMACEM V+ Injection Cannula", please consult the corresponding "Instructions for Use".

### Contraindications

TFNA SHORT (lengths 170 mm, 200 mm, 235 mm)

- Femoral neck fractures (31-B)
- Femoral shaft fractures (32-A/B/C)

TFNA LONG (lengths 260 mm–480 mm)

- Femoral neck fractures (31-B)

TFNA AUGMENTATION

- Risk for intraarticular or vascular cement leakage
- Acute traumatic fractures with good bone quality

### Patient Target Group

The TFNA implants are recommended for use in skeletally mature patients.

### Intended User

This IFU alone does not provide sufficient background for direct use of the device or system. Instruction by a surgeon experienced in handling these devices is highly recommended.

This device is intended to be used by qualified health care professionals e.g. surgeons, physicians, operating room staff, and individuals involved in preparation of the device. All personnel handling the device should be fully aware of the IFU, the surgical procedures, if applicable, and/or the Synthes brochure "Important Information" (SE\_023827) as appropriate.

Implantation is to take place according to the instructions for use following the recommended surgical procedure. The surgeon is responsible for ensuring that the device is suitable for the pathology/condition indicated and that the operation is carried out properly.

### Expected Clinical Benefits

Expected clinical benefits of internal fixation devices such as TFNA implants when used according to instructions for use and recommended technique are:

- Stabilize bone segment and facilitate healing
- Restore anatomical alignment and extremity function

### Performance Characteristics of the Device

TFNA is designed to reduce the risk of post-operative complications associated with hip fractures by providing surgical options to enhance stability in poor bone, improved anatomical fit, and increased implant strength. It is also designed to provide a range of options to support surgical preferences and patient anatomies including choice of augmentable blade or screw head elements, various locking options, and a range of nail sizes.

### General Notes:

- Implants are designed for temporary fixation. Therefore, if bone consolidation is not sufficient the system may fail over time.
- There are many types of implant failures, including but not limited to implant breakage.
- There are several factors that can influence implant failure including fracture reduction, surgical technique, obesity, level of activity/weight-bearing, and non-union or delayed union. Surgeons should consider these factors in intra-operative care for bone consolidation. These failures can occur post-operatively and may require reoperation.
- The aim of post-operative care must be the promotion of bone consolidation.


### Potential Adverse Events, Undesirable Side Effects and Residual Risks

- Adverse tissue reaction, allergy/hypersensitivity reaction
- Bone damage including intra- and post-operative bone fracture, osteolysis, or bone necrosis
- Damage to vital organs or surrounding structures
- Embolism
- Infection
- Injury to user
- Malunion/Non-union
- Neuro-vascular damage
- Pain or discomfort
- Poor joint mechanics
- Soft tissue damage (including compartment syndrome)
- Symptoms resulting from implant migration, loosening, bending, cut-out, or breakage

### Sterile Device


**STERILE R** Sterilized using irradiation

Store sterile devices in their original protective packaging, and do not remove them from the packaging until immediately before use.

 Do not use when packaging is damaged.

Prior to use, check the product expiration date and verify the integrity of the sterile packaging. Do not use if the package is damaged or date of expiration has passed.

### Single-use Device

 Do not re-use.

Indicates a medical device that is intended for one use, or for use on a single patient during a single procedure.

Re-use or clinical reprocessing (e.g. cleaning and resterilization) may compromise the structural integrity of the device and/or lead to device failure which may result in patient injury, illness or death.

Furthermore, reuse or reprocessing of single-use devices may create a risk of contamination e.g. due to the transmission of infectious material from one patient to another. This could result in injury or death of the patient or user.

Contaminated implants must not be reprocessed. Any Synthes implant that has been contaminated by blood, tissue, and/or bodily fluids/matter should never be used again and should be handled according to hospital protocol. Even though they may appear undamaged, the implants may have small defects and internal stress patterns that may cause material fatigue.

## Warnings and Precautions

The general risks associated with surgery are not described in these instructions for use. For more information, please refer to the Synthes brochure "Important Information" (SE\_023827).

It is strongly advised that TFNA implants are implanted only by operating surgeons who are familiar with the general problems of trauma surgery and who are able to master the product-specific surgical procedures. Implantation is to take place with the instructions for the recommended surgical procedure. The surgeon is responsible for ensuring that the operation is carried out properly.

The manufacturer is not responsible for any complications arising from incorrect diagnosis, choice of incorrect implant, incorrectly combined implant components and/or operating techniques, the limitations of treatment methods, or inadequate asepsis.

## WARNINGS

- It is critical to ensure proper selection of the implant meets the needs of the patient anatomy and the presenting trauma.
- The TFNA Nail is not intended for full weight-bearing in patients with complex unstable fractures until sufficient bone consolidation is confirmed in the follow up X-rays.
- Conditions that place excessive stresses on bone and implant such as severe obesity or degenerative diseases should be considered. The decision whether to use these devices in patients with such conditions must be made by the physician taking into account the risks versus the benefits to the patients.
- Use of these devices is not recommended when there is systemic infection, infection localized to the site of the proposed implantation or when the patient has demonstrated allergy or foreign body sensitivity to any of the implant materials.
- Physician should consider patient bone quality to ensure it provides adequate fixation to promote healing.
- Compromised vascularity in the site of proposed implantation may prevent adequate healing and thus preclude the use of this or any orthopaedic implant.
- Physician should take into account an increase in medullary pressure occurring during medullary nailing or reaming. This releases varying amounts of bone marrow and fat into the venous blood system.
- Do not augment if X-ray contrast media leaks into the joint.
- The 6 mm minimum distance is recommended to reduce the risk of thermal injury to the adjacent cartilage tissue.
- In the event that there is danger of cement leakage into the joint, fracture gap or venous system, stop injection immediately.
- If the extravasated cement conforms to the architecture of the hip joint, it may not need to be removed. However, if it does not conform and is abrasive or damages the articular surface, then the extruded cement will need to be removed.
- To remove the cement, the treating physician has the option of either hip arthroscopy, arthroplasty, or open arthrotomy to remove the extruded pieces. The timing of the removal is at the discretion of the physician after appropriate evaluation of the patient.

## PRECAUTIONS

For additional precautions specific to a surgical step, please refer to the Special Operating Instructions.

## Combination of Medical Devices

Synthes has not tested compatibility with devices provided by other manufacturers and assumes no liability in such instances.

## Magnetic Resonance Environment

### Torque, Displacement and Image Artifacts according to ASTM F 2213-06, ASTM F 2052-14 and ASTM F2119-07

Non-clinical testing of worst case scenario in a 3 T MRI system did not reveal any relevant torque or displacement of the construct for an experimentally measured local spatial gradient of the magnetic field of 3.69 T/m. The largest image artifact extended approximately 169 mm from the construct when scanned using the Gradient Echo (GE). Testing was conducted on a 3 T MRI system.

### Radio-Frequency-(RF)-induced heating according to ASTM F2182-11a

Non-clinical electromagnetic and thermal testing of worst case scenario lead to peak temperature rise of 9.5 °C with an average temperature rise of 6.6 °C (1.5 T) and a peak temperature rise of 5.9 °C (3 T) under MRI Conditions using RF Coils (whole body averaged specific absorption rate [SAR] of 2 W/kg for 6 minutes [1.5 T] and for 15 minutes [3 T]).

## Precautions:

The above mentioned test relies on non-clinical testing. The actual temperature rise in the patient will depend on a variety of factors beyond the SAR and time of RF application. Thus, it is recommended to pay particular attention to the following points:

- It is recommended to thoroughly monitor patients undergoing MR scanning for perceived temperature and/or pain sensations.
- Patients with impaired thermoregulation or temperature sensation should be excluded from MR scanning procedures.
- Generally, it is recommended to use a MR system with low field strength in the presence of conductive implants. The employed specific absorption rate (SAR) should be reduced as far as possible.
- Using the ventilation system may further contribute to reduce temperature increase in the body.

## Treatment before Device is used

### Non-sterile Device:

Synthes products supplied in a non-sterile condition must be cleaned and steam-sterilized prior to surgical use. Prior to cleaning, remove all original packaging. Prior to steam-sterilization, place the product in an approved wrap or container. Follow the cleaning and sterilization instruction given by the Synthes brochure "Important Information" (SE\_023827).

### Sterile Device:

The devices are provided sterile. Remove products from the package in an aseptic manner.

Store sterile devices in their original protective packaging, and do not remove them from the packaging until immediately before use.

Prior to use, check the product expiration date and verify the integrity of the sterile packaging. Do not use if the package is damaged.

## Implant Removal

In case the physician decides to remove the implants, the following steps shall be followed:

1. Carefully dissect the soft tissues and visualize the end cap. Remove the end cap with a retaining Synthes screwdriver.
2. Carefully dissect the soft tissues and visualize the screw heads. In the case of screw head overgrowth or damaged recess, optional instruments are available for screw removal for example if required, a curette and a sharp hook to clear recess from tissue; an extractor shaft and a conical extraction screw to remove screws with damaged recess. Remove all locking screws.
3. Thread the extraction screw into the nail.
4. Disengage locking mechanism and remove helical blade or screw.
5. Remove the nail.


## Troubleshooting

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

## Clinical Processing of the Device


Detailed instructions for processing of implants and reprocessing of reusable devices, instrument trays and cases are described in the Synthes brochure "Important Information" (SE\_023827).


## Device-related Storage and Handling Information

 Caution, see instructions for use.

## Additional Device-specific Information

 Reference number

 Lot or batch number

 Legal manufacturer

 Expiration date

## Disposal

Any Synthes implant that has been contaminated by blood, tissue, and/or bodily fluids/matter should never be used again and should be handled according to hospital protocol.

Devices must be disposed of as a healthcare medical device in accordance with hospital procedures.

## Special Operating Instructions

**1. Position the patient** in the lateral decubitus or supine position on a fracture table or radiolucent operating table. Position the image intensifier to allow visualization of the proximal femur in both the AP and lateral planes.

### 2. Reduce fracture.

Precaution: Instruments and screws may have sharp edges or moving joints that may pinch or tear user's glove or skin.

Precaution: Handle devices with care and dispose worn bone cutting instruments in an approved sharps container.

### 3. Determine CCD angle.

### 4. Determine nail length and diameter.

Measure using the radiographic ruler.

Alternative: Nail length may also be determined by using a reaming rod.

### 5. Identify nail entry point.

Make a longitudinal incision proximal to the greater trochanter. In the AP view, the nail insertion point is on the tip or slightly lateral to the tip of the greater trochanter, in the curved extension of the medullary cavity. This represents a point, 5° lateral of the femoral shaft axis, measured from a point just below the lesser trochanter, as the ml angle of the nail is 5°.

### 6. Insert guide wire.

Confirm guide wire placement in both the AP and lateral planes. Insert to a depth of approximately 15 cm.

### 7. Open canal.

Guide the cannulated drill bit over the guide wire through the protection sleeve to the bone and drill to the stop.

Precaution: Guide wires are single-use items, do not re-use.

Option: Open canal with the hollow reamer.

Precaution: Monitor the drill depth under image intensifier throughout the procedure.

### 8. Option: Ream medullary canal.

Option: The reamer protection tube can be used to help protect the proximal metaphysis during reaming.

### 9. Assemble insertion instruments.

Match the geometry of the insertion handle to the nail. Pass the connecting screw through the insertion handle and into the nail. Secure the assembly with the ball hex screwdriver.

Precaution: Ensure that the connection between the nail and the insertion handle is tight (retighten if necessary).

Precaution: Do not attach the aiming arm to the insertion handle yet.

Precaution: If a 235 mm or longer nail is selected, reconfirm that the correct nail (right or left) is assembled.

### 10. Insert nail.

Under image intensification, verify fracture reduction and insert the nail as far as possible by hand. Use the insertion assembly to manipulate the nail across the fracture.

Option: To use a hammer, screw the driving cap onto the hybrid insertion handle. Monitor the tip of the nail using image intensification.

Precaution: Using light blows, the hammer can also be used with the hammer guide to back slap the nail if the nail has been slightly over inserted.

Precaution: Confirm that the nail is tightly connected to the insertion handle as hammering may loosen the connection.

### 11. Verify nail insertion depth and anteversion.

Verify nail insertion depth and position for the helical blade/screw. Adjust nail rotation.

### 12. Insert guide sleeve

Precaution: The distal tooth of the guide sleeve should rest on the lateral cortex. Do not overtighten on the cortex as this may affect the accuracy of the aiming assembly.

Precaution: The fatigue strength of the nail may be affected and may contribute to the potential for the nail to fracture if the nail is damaged during any step of the helical blade/screw reaming in addition to other factors such as fracture reduction, surgical technique, obesity, level of activity/weight-bearing, non-union, or delayed union.

### 13. Insert guide wire for helical blade/screw

Precaution: If the nail must be repositioned to improve guide wire placement, remove the guide sleeve assembly and adjust with the insertion handle. Make a new incision for insertion of the guide sleeve, if necessary. Do not pull on the guide sleeve or power tool to make this adjustment as this could affect the accuracy of the aiming.

Precaution: The fatigue strength of the nail may be affected and may contribute to the potential for the nail to fracture if the nail is damaged during any step of the helical blade/screw reaming in addition to other factors such as fracture reduction, surgical technique, obesity, level of activity/weight-bearing, non-union, or delayed union.

Precaution: Do not reuse guide wires as they may bend during initial use. If the guide is deformed during insertion, use a new guide and discard the deformed guide wire.

Precaution: Insert the guide wire for the blade or screw carefully to avoid penetration of the guide wire into the joint. Penetration of the articular surface is a contraindication for the augmentation of the blade or screw.

### 14. Measure helical blade/screw length.

### 15. Open lateral cortex for helical blade/screw insertion.

Precaution: Monitor the drill depth under image intensifier throughout the procedure.

### 16. Option A: Helical blade insertion.

Precaution: Image intensifier should be used during helical blade insertion to monitor positioning.

Precaution: Assure that the guide wire is in place while inserting the helical blade to prevent the cannulation from clogging, impeding an optional augmentation procedure.

### 17. Option B: Screw insertion.

Precaution: There is no stop on the tap, therefore monitoring insertion via the following methods is recommended:

– Monitor the depth under image intensifier

– Monitor the respective graduations of the instrument shaft in relation to the guide sleeve

Precaution: Image intensifier should be used during screw insertion to monitor positioning.

Precaution: Assure that the guide wire is in place while inserting the screw to prevent the cannulation from clogging, impeding an optional augmentation procedure.

### 18. Rotational locking.

Precaution: If the locking mechanism is not turned back 1/2 turn after initial tightening as described above, controlled collapse and compression of the fracture may not occur.

### 19. Interfragmentary compression (option).

Precaution: Caution should be taken when using the buttress/compression nut with the pin wrench to avoid over compression that could potentially cause the helical blade to lose purchase in the bone, especially in patients with poor bone quality.

### 20. Augmentation.

It is recommended to use 3 ml of cement for augmentation.

Precaution: The working time for TRAUMACEM V+ Injectable Bone Cement at room temperature (20 °C) is approximately 27 minutes. At body temperature (37 °C) the setting time is 15 minutes. After last cement injection, the patient should remain immobile for 15 minutes to facilitate proper cement curing.

Precaution: Use only radiographic contrast agents that are indicated for this application.

Precaution: Consult the manufacturer's directions on indications, contraindications, use, precautions, warnings and side effects of the radiographic contrast agent.

Precaution: Always use the full amounts of monomer liquid and polymer powder provided in the kit, respectively, when mixing TRAUMACEM V+ Injectable Bone Cement. Otherwise, the behavior of the TRAUMACEM V+ Injectable Bone Cement can no longer be guaranteed. Using only one of the components is not permitted.

Precaution: Ensure that the powder and liquid component are thoroughly mixed before starting cement transfer.

Precaution: Ensure a good fit between the syringe and the stop-cock/used access solution, but make sure to be on axis and avoid using excessive force when coupling them. They are both made of plastic and could otherwise break.

Precaution: Do not advance the cannula more than 5 mm over the selected head element length. This would result in injection of cement in front of the head element tip where no additional stability is achieved and the risk of penetration and cement leakage is increased.

## 21. Distal locking.

Precaution: Confirm that the nail is securely connected to the insertion handle, especially after hammering.

Read length from the calibrated drill bits or measure length using depth gauge for locking screws.

Precaution: Select adequate screw length to avoid protrusion of the screw tip and irritation of soft tissue.

5.0 mm locking screws can be connected to the screwdriver with the associated retention pin; this does not apply to the alternative locking screws Stardrive®.

Precaution: The screw must not be tightened with the power tool. Disengage the power tool from the screwdriver shaft before the screw is fully seated and use the manual handle to bring the screw to its final position.

Ensure drill bits, screws, nuts or washers do not interfere with other medical devices (e.g. knee prosthesis, nail, other screws) and/or critical anatomy (e.g. condylar notch, joint space).

Nut should only be used with the 5.0 mm Locking Screws for Medullary Nails.

## 22. Insert end cap.

**0 mm end cap:** Remove the connecting screw using the ball hexagonal screwdriver while leaving the insertion handle connected to the nail. Insert the 0 mm end cap through the insertion handle.

**5–15 mm end cap:** Remove the connecting screw and insertion handle using the hexagonal screwdriver. Insert the end cap.



Synthes GmbH  
Eimattstrasse 3  
4436 Oberdorf  
Switzerland  
Tel: +41 61 965 61 11  
[www.jnjmedicaldevices.com](http://www.jnjmedicaldevices.com)