eMax 2 Plus High Speed Electric System

Contents

General Information.....	2
System Assembly	6
System Operation.....	6
System Shutdown.....	8
Additional Information.....	8
Console Alerts, Codes, and Displays	8
Trouble Shooting	9
Technical Specifications (IEC-60601-1-2 3rd Edition)	10
Technical Specifications (IEC-60601-1-2 4th Edition)	14

General Information

System Description

High Speed Electric Drill System consisting of a console (SC2100, SC2101 or SC2102), handpiece (EMAX2PLUS), foot control (E-FP, E-FP-DIR or E-FP-DIR/IRR), country or region specific medically rated A/C Power Cord, attachment (various options) and dissection tool (various options). There is also an optional hand control (EMAX2-HC or E-HC) and irrigation accessories.

Indications for Use

Cutting and shaping bone including bones of the cranium and spine.

Cleaning and Sterilization

- Prior to first use the equipment must be processed as per the included Anspach Cleaning and Sterilization instructions.
- At the point of use the device must be cleaned as soon as possible after use to prevent drying of blood, tissue, other biological debris and contaminants on the device.
- Refer to the Anspach Cleaning and Sterilization instructions accompanying this device.


Features & Compatibility

		Consoles			Foot Control Features		Hand Controls	
		SC2100	SC2101	SC2102	Active Direction Control Switch	Active Irrigation Control Switch	EMAX2-HC	E-HC
Foot Controls	E-FP	•	•	•				
	E-FP-DIR	•	•	•	•			
	E-FP-DIR/IRR	•	•	•	•	•		
	EPLUS-FP	•	•	•	•	•		
	EPLUS-FP-NS	•	•	•				
Hand-piece	EMAX2PLUS	•	•	•			•	•
Console Features	Two Foot Control Ports *	•						
	Two Handpiece Ports *	•						
	Console Irrigation	•	•					


* If two foot controls or two handpieces are connected upon startup of console it will default to Foot Control 1 and/or Handpiece 1.

The following instructions are written such that the handpiece and attachments are held with distal end pointed away from user.

Symbols

	Anspach Product Icon		Consult Operating Instructions
	Located on the Foot Pedal switch this indicates a change state (on/off) for the console Irrigation System		Located on the Foot Pedal switch this indicates a change of rotational direction for the handpiece
	Direction of rotation		Keep Dry (A.K.A. Protect from Moisture)
	Direction of rotation for lock position		Indicates position, alignment or location
	Indicates setting, position or location		AC Power Source
	Indicates setting, position or location		Indicates position or location
	Indicates the attachment setting, refer to MA-DRIVER & MA-D20 section in this document for further details		Located on the Console this indicates motor or handpiece connector
	Indicates the drill tip exposure		Symbol with yellow background: CAUTION: refer to accompanying documentation

Symbols

	Located on the Console this indicates Hand Control activation		Located on the Console this indicates Foot Control connector
30 Sec ON / 30 Sec OFF	Duty Cycle of the DRIVER	RPM	Revolution per Minute
	Located on the Console this indicates Forward (Clockwise rotation when looking at handpiece from proximal end).		Located on the Console this indicates Reverse (Counter-Clockwise rotation when looking at the handpiece from proximal end)
	Indicates action of rotation and position for Secure (LOCK) and Release (UNLOCK)		Located on the Console this indicates a change of irrigation state for the Console Irrigation System
	Direction of rotation		Direction of rotation
	Direction of rotation for lock position		1. Console Control Touch pad: indicates a change state (on/off) for the Console Irrigation System 2. Irrigation Pump: Indicates irrigation tube direction towards the irrigation source
	Located on the Foot Pedal switch this indicates a change of rotational direction for the handpiece		Located on the Foot Pedal switch this indicates a change state (on/off) for the Console Irrigation System
	Located on the Console this indicates a change of rotational speed rate for the Handpiece		NO NOT DISPOSE OF IN HOUSEHOLD WASTE
	Temperature Limits		Sterile unless damaged or open
CE 2797	CE Mark (A.K.A. CE Mark [Notified body number], Conformité Européenne) Meaning: Device complies with applicable EEC Directives	Rx ONLY	United States Federal Law restricts this device to sale by or on the order of a physician or other licensed health-care provider
	Found on the Console ON/OFF		Located in the Console Irrigation Pump this indicates a possible finger pinch point
	Lock (Run)		Unlock (Load)
	Located on the Console Irrigation Pump these symbols are not used by the user since the pump is self-adjusting to the Irrigation Tube Diameter		Located on the Handpiece this indicates position of the knurl knob to the desired location
REF	Reference Number (A.K.A. Item Number, Catalogue Number, Part Number)	SN	Serial Number
LOT	Lot (A.K.A. Lot Number, Batch Number, Batch Code)		Transmitter Interference may occur in the vicinity of equipment marked with this symbol.
	Manufacturer		Type B Electrical Equipment
	Use By Date (A.K.A. Expiration Date, Expiry)		Product warranty is void if seal is damaged or removed.
EC REP	Authorized European Union Representative		Single Use Only (A.K.A. Do Not Reuse)
STERILE R	Sterilized Using Irradiation		CAUTION: Refer to accompanying documentation
	Date of Manufacture	(01)	Global Trade Item Number (GTIN)
(10)	Batch or Lot Number	(17)	Expiration Date (YYMMDD)
(21)	Serial Number	QTY	Quantity
GTIN	Global Trade Item Number		Indicates Environmental Friendly Use Period of 50 years in China.
	Indicates Environmental Friendly Use Period of 10 years in China.		Indicates Environmental Friendly Use Period in China, this product does not contain any hazardous substances exceeding concentration limits.
	Certified/Accredited Test House		Compliant with IEC 60601-1-2 Ed. 4.0 (2014-02) Electromagnetic compatibility.

Warnings

- **⚠ Always use continuous irrigation to prevent heat build-up. Irrigation is necessary for proper performance.**
- Surgeon is responsible for learning proper techniques in use of equipment; improper use may cause serious injury to user or patient or damage to system.
- Instrument operator and all operating room personnel must wear eye protection.
- Visually inspect for damage before using; do not use if damage is evident.
- Do not use if the product sterilization barrier or its packaging is compromised.
- Do not use, or discontinue use of powered equipment exhibiting excessive temperatures that can cause patient injury (necrosis) and/or user discomfort.
- Use of damaged or improperly maintained power equipment and/or misused powered equipment can result in excessive temperatures.
- Use caution to avoid cutting or tearing gloves while handling dissection tools.
- Dissection tools must be adequately retained within attachment to prevent distal migration, which may cause injury. Confirm attachment is proper size for dissection tool and that it is secure.
- Gently pull on dissection tool shaft to ensure it is fully seated and properly installed.
- Only cut visible areas unless an image intensifier is utilized.
- Delicate structures in proximity to dissection must be thoroughly protected to prevent injury.
- Maintain firm control of instrument at all times.
- Do not bend or use as a lever.
- Use a gentle tapping motion or side-to-side motion and let instrument do cutting.
- Do not use excessive force.
- Forceful side loading of dissection tool may cause fracture of dissection tool, which may cause injury.
- Dissection tools are disposable and intended for single patient use only. Do not re-sterilize and/or re-use dissection tools. Reuse and/or re-sterilization of these devices could result in reduced performance, material degradation and/or patient contamination
- Use standard protocol for disposal of sharp instruments.
- Continuous extreme cutting at or near stalling conditions will quickly overheat handpiece.
- Do not operate in an explosive flammable environment.
- Do not modify ground or power cord.
- Do not allow liquid into console.
- Do not use a craniotome in conjunction with Hand Control. Excessive rotational force on hand control may move it, thus shutting handpiece down.
- Use of accessories or cables other than those provided by DePuy Synthes Power Tools and specified for use with the eMax 2 Plus System may result in increased emissions or decreased immunity.
- The eMax 2 Plus should not be used adjacent to or stacked with other equipment and that if adjacent or stacked use is necessary, the eMax 2 Plus should be observed to verify normal operation in the configuration in which it will be used.
- Medical Electrical Equipment needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information provided in this accompanying documentation.
- The emissions characteristics of this equipment make it suitable for use in industrial areas and hospitals (CISPR 11 class A). This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instruction manual, may cause harmful interference to radio communications.
- Operation of this equipment in a residential area (for which CISPR 11 class B is normally required) this equipment might not offer adequate protection to radio frequency communication services.
- In the event of interference, the user might need to take mitigation measures such as relocating or reorienting the equipment.
- Portable and Mobile RF communications equipment can affect Medical Electrical Equipment.
- Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the EMAX2PLUS system. Otherwise, degradation of the performance of this equipment could result.
- Power source should comply with applicable IEC, CEC, and NEC requirements. Grounding reliability can only be achieved when this equipment is connected to a receptacle marked "HOSPITAL GRADE."
- Do not use in oxygen rich environment.
- No modification of this equipment is allowed.
- Do not modify. Modifications could result in loss of electrical safety.
- Dispose of items contaminated with body fluids with other biohazardous waste.
- At end of life recycle or dispose of device in accordance with local and national regulations.
- To avoid the risk of electric shock, this equipment must only be connected to a supply mains with protective earth. The use of Accessories, transducers, and cables other than those specified, with the exception of transducers and cables sold by the Manufacturer of this device as replacement parts for internal components, may result in increased Emissions or decreased Immunity of the eMax 2 Plus System.

Cautions

- United States Federal law restricts this device to sale by or on order of a physician or other licensed healthcare provider.
- Do not use accessories other than those provided by DePuy Synthes Power Tools and specified for use with Anspach Systems.
- To insure equipment operates as designed, read and follow manufacturer's instructions.
- Do not operate handpiece without an attachment and the corresponding dissection tool.
- Do not engage safety mechanism while handpiece is running; doing so makes safety mechanism inoperable.
- Only Anspach Dissection Tools should be used with eMax 2 Plus Systems.
- Use care to protect hose when handling, cleaning, and during system use.
- Damage to hose can cause leaking, rupture, or other related failures.
- Do not step on, set equipment on, pinch, kink, clamp, or otherwise occlude handpiece hose during use.
- Do not engage safety mechanism while handpiece is running; doing so makes safety mechanism inoperable.
- Do not use the eMax 2 Plus System with MS-OSC or MS-SAG; attempting to do so will result in damage to saws and possible damage to the eMax 2 Plus Handpiece.

Latex Information

Not made with natural rubber latex.

Warranty & Return Policies

Warranty and return policy is available upon request.

Warning: All Instrument System components returned for servicing or repair should be properly cleaned and sterilized as applicable prior to shipping.

Warning: Transmissible Spongiform Encephalopathies (TSE)

DePuy Synthes Power Tools will not authorize or accept the return of products that directly contact patients or are contaminated with a patient's body fluids who is suspected or confirmed with a Transmissible Spongiform Encephalopathies/Creutzfeldt-Jakob Disease (TSE/CJD) diagnosis. DePuy Synthes Power Tools recommends that all Anspach products used on a patient confirmed with a TSE/CJD diagnosis be incinerated. Anspach Dissection Tools used on a patient suspected of TSE/CJD diagnosis must be incinerated.

Contact your Sales Representative for replacement of product incinerated under this policy or for temporary equipment while original equipment is quarantined. Contact the DePuy Synthes Power Tools Customer Service Department regarding TSE/CJD contamination for additional information.

System Assembly

Warning

Prior to first use the equipment must be processed as per the included Anspach High Speed System Cleaning and Sterilization instructions.

Note: When following these instructions hold handpiece and attachments with distal end pointing away from user.

1. Plug hospital grade AC power cord into AC power inlet on back of console. Plug opposite end of power cord into standard, hospital grade grounded wall outlet. Connect to supply mains with protective earth only.
2. Insert foot control cord connector into Foot Control 1 connector port. (A second foot control cord may be connected into Foot Control 2 connector port on applicable consoles.) Foot control cord connector is keyed. Align foot control cord connector with connector port on front of console. Do not push foot control cord connector into console connector port when out of alignment.
3. Plug handpiece hose connector into Handpiece 1 connector port located on front of console. (A second handpiece may be connected into Handpiece 2 connector port on applicable consoles.) Handpiece hose connector is keyed. Align handpiece hose connector with connector port on front of console. Do not push handpiece hose connector into console connector port when out of alignment.
4. Activate console by depressing power switch located on back of console to "I" position. Appropriate face panel light emitting diodes (LED) illuminate and a beep will sound.
5. For operating instructions see Console Face Panel Operation and Foot Control Operation below.

Note: Handpiece is fully functional at this time.

Irrigation Assembly for Consoles with Optional Irrigation (Excludes SC2102)

1. Upon startup "Irrigation" is OFF, and LED is not illuminated.
2. Insert IV pole into bracket on back of console.
3. Hang irrigation bag onto IV pole.
4. Check integrity of irrigation tubing package and use by date.

Warning: Discard product if sterilization barrier or its packaging is compromised or if past its expiration date.

5. Remove tubing and hose clips from package.
6. Connect small diameter tubing to proximal end of attachment clip.
7. Secure clip to handpiece and attachment.
8. Route tubing along handpiece hose and secure to hose using hose clips.
9. Remove protective cap from bayonet and insert bayonet into irrigation bag.
10. Hold irrigation bag and tubing vertical and allow fluid to completely fill large diameter tubing.
11. Depress button located on right side of console to eject pump head.
12. Lift top of pump head and insert large diameter section of tubing into pump. Refer to diagram on pump housing to ensure tubing is routed properly (smaller diameter tubing flows towards Anspach Handpiece).

Note: Verify there are no air bubbles between irrigation bag and pump. If air bubbles are present, remove tubing from pump and allow fluid to completely fill large diameter tubing and reinstall to pump.

13. Check position of tube in pump to avoid damage to irrigation tubing and close pump head.
14. To activate, depress "Irrigation" button on face panel. LED illuminates.
15. To control flow of irrigation, depress arrows located to right of irrigation button on face panel.

Note: Irrigation can also be activated from certain optional Foot Controls (see Step 5 of Foot Control Operation).

16. Depress Irrigation Control Switch on top right of foot control. Upon activation, LED on face panel illuminates.

Caution: Tubing can disconnect from connectors without warning if occluded. Do not step on, set equipment on, pinch, kink, clamp, or otherwise occlude tubing during use. Do not operate with pump head open.

Warning: Irrigation tubing, attachment clip, and hose clip are SINGLE USE ONLY. Reuse and/or re-sterilization of these devices could result in reduced performance, material degradation and/or patient contamination. For proper operation only use at ambient temperatures between 65°F to 85°F (18°C to 29°C).

System Operation

Upon startup, console is set to operate handpiece in forward (clockwise when viewed from proximal end of handpiece) direction at maximum speed and display reads 80,000 RPM. System operation, including optional irrigation, can be controlled from console face panel. On applicable consoles LED's on face panel correspond to connector ports on connector panel. System operation can also be controlled by several available foot control options and an optional hand control.

Console Face Panel Operation

1. To increase or decrease speed of handpiece, depress blue arrows located on console face panel. Speed increases and decreases in 10,000 RPM increments.
2. To change direction, depress "R" or "F" arrow located on console face panel. Console beeps once, indicating handpiece direction has changed. Direction can only be changed when handpiece is not running. A series of three beeps indicates console is set to operate handpiece in reverse/counterclockwise direction. LED on console face panel located to left of "F" indicates forward direction. LED located to left of "R" indicates reverse direction.
3. (Optional) To activate irrigation, depress "Irrigation" button on face panel. LED illuminates. To control flow of irrigation, depress arrows located to right of irrigation button on face panel.

Foot Control Operation

Certain optional foot controls have direction and irrigation control switches.

1. See Step 2 of Console Face Panel Operation to set operational direction of handpiece.
2. (Optional) Depress "Foot" button on face panel to activate Foot Control 1, depress again to activate Foot Control 2. Only one foot control may be active at a time.

Note: Handpiece is fully functional at this time.

3. Press on foot control pedal to start handpiece. Increase pressure on pedal to increase speed of handpiece and release pressure on pedal to decrease speed of handpiece.
4. (Optional) To change direction of handpiece, depress Directional Control Switch located on top left of foot control for a minimum of one second. When direction changes, a single beep will sound and LED on console face panel indicates direction of rotation. Direction can only be changed when handpiece is not running.
5. (Optional) To activate irrigation, depress Irrigation Control Switch on top right of foot control for a minimum of one second and then step on center of foot control pedal. Upon activation, LED on face panel illuminates; irrigation pump and handpiece start.

Hand Control Assembly and Operation

Optional hand-control attachment for handpiece allows user to control operation with a lever attached to handpiece.

Warning: Do not expose handpiece to magnets (such as laying handpiece on magnetic drapes) when in hand control mode, as handpiece operation may occur.

Warning: It is recommended that a craniotome not be used in conjunction with Hand Control. Excessive rotational force on hand control may move it resulting in handpiece shutting down.

1. Upon console startup Hand Control is OFF. LED is not illuminated.
2. Rotate handpiece until small dimple at proximal end (hose end) of handpiece faces up. Slide hand control over distal end of handpiece such that notch on hand control slides over dimple at proximal end (hose end) of handpiece. Slide Hand Control onto handpiece until secure.
3. To activate Hand Control, depress "Hand" button on console face panel. This deactivates Foot Control; Hand Control LED on face panel illuminates.
4. To prevent inadvertent operation, ensure silver lever on Hand Control is fully inserted in "Lock" position.
5. Extend silver lever from "Lock" position for operation.
6. Apply pressure on Hand Control to operate handpiece.

Note: Handpiece is fully functional at this time.

7. To deactivate Hand Control, depress "Foot" button on console face panel. To remove hand control from handpiece, slide towards distal end of handpiece and remove.

System Operational Assessment

1. Confirm attachment and dissection tools are properly locked in place.
2. Depress foot pedal to assure you have properly functioning device prior to proceeding.
3. Confirm device functions as expected. If device does not function as expected refer to troubleshooting chart.

Prior to Use Inspection	
Straight Attachments	Visually inspect for any damage to the tube.
Angle Attachments	Visually inspect for bent or broken drive shaft and for any damage to the tube.
Craniotomes	Visually inspect for bent or broken foot.
Perforator Driver	Visually inspect for overall damage or missing components.
Foot Pedals	Visually inspect for damage to the electrical cord or connector. Visually inspect for damage or cracks to the housing or pedal.
EMAX2PLUS	Visually inspect for damage to the silicone hose or to the electrical connector. Connect to console and operate. The handpiece should operate smoothly. There is no requirement to operate with attachment or dissection tool.
SC2100	Visually inspect for damage or cracks to the housing.
SC2101	Visually inspect for damage to the electrical power cord.
SC2102	Power the system and inspect for illumination of the LEDs. If there is an irrigation pump present check for proper function.

NOTE: DRIVER Duty Cycle 30 Sec ON, 30 Sec OFF for eight cycles.

The recommendations for times of use for the DRIVER attachments have been determined under average load and worst case ambient air temperature of 85°F (29°C).

System Shutdown

1. There are no special procedures for system shutdown. Depress power switch on rear of console to shut down system.
2. To disconnect power cord from mains, grasp ridged portion of connector between thumb and forefinger of one hand. Gently pull connector away from console. Connector should remove easily. If not, ensure only ridged portion of connector is pulled. Do not pull on smooth part of connector, as connector will not disengage from console.

Additional Information

End of Life

Return all equipment for proper disposal.

Recommended Manufacturer Service Interval

It is recommended that the equipment be returned to DePuy Synthes Power Tools at a minimum every 9 months for full product inspection and service.

Console Alerts, Codes, and Displays

Console Indicators

Code	Description	Audible Signal **
E1, E8 *	System Fault	none
E2	Handpiece Lock Engaged	10 slow beeps
E3	Handpiece Stall Warning	10 fast beeps
E4, E5 *	Fault in Handpiece	none
E6, E7	Handpiece Overheat Warning	5 second beep
E9 *	Handpiece Not Recognized	none
(none)	Handpiece Set to Reverse	3 beeps

* If these codes are displayed, return handpiece and console to DePuy Synthes Power Tools for servicing.

** Note: Consoles with software version 3.0 and higher do not generate an "Audible Signal" for Codes E1 – E9. To determine software version press the irrigation rate buttons at the same time on the console touchpad. The console will display two digits. Example Only: "30" means version 3.0. On non-irrigation console models, the buttons can be found in the same location, although they are not marked.

Trouble Shooting

Trouble Shooting Chart		
Problem	Possible Cause	Possible Solution
Excessive handpiece noise.	1. Faulty internal component.	1. Return attachment to DePuy Synthes Power Tools for service.
Lack of power to console.	1. Plug is not fully inserted/power switch not turned on. 2. Outlet is not functional.	1. Verify plugs to wall and back of console are pushed in completely. 2. Verify another piece of electrical equipment can receive power from outlet.
Handpiece vibration or extremely hot.	1. Faulty internal component. 2. Internal motor overheats due to continuous extreme cutting at near stalling conditions. 3. Hose may be kinked. 4. Handpiece out of balance. 5. Handpiece hose damaged.	1. Return to DePuy Synthes Power Tools for service. 2. Decrease cutting force, allow motor to cool. 3. Unkink hose (verify it is not being pinched or clamped to table). 4. Return to DePuy Synthes Power Tools for service. 5. Return to DePuy Synthes Power Tools for service.
Excessive vibration of dissection tool.	1. Dissection tool may be bent. 2. Dissection tool may not be fully seated. 3. Improper attachment and dissection tool combination. 4. Possible attachment bearing damage.	1. Replace with new dissection tool. 2. Reassemble dissection tool and attachment. 3. Only use correct dissection tool with appropriate attachment. 4. Return to DePuy Synthes Power Tools for service.
Attachment is hot.	1. Debris lodged inside attachment. 2. Possible bearing damage.	1. Clean attachment using Anspach Attachment Cleaning Instructions. 2. Return attachment to DePuy Synthes Power Tools for service.
Craniotome attachment is bent.	1. Excessive force used in operation.	1. Do not use, replace attachment.
Inoperative Foot Control (LEDs are illuminated on face panel)	1. Selected foot control on face panel does not match foot control in use 2. System is set in hand control mode (LED light is illuminated) 3. Plug may not be fully inserted 4. Handpiece in "Safe" position or cutter is not rotating 5. Defective attachment	1. Press Foot Control button on face panel to select other Foot Control 2. Press Foot Control button on face panel to set it for Foot Control 1 or 2, depending which one is in use. 3. Verify plug for handpiece and foot control are fully inserted. Plug will latch in place once fully inserted. 4. Verify cutting burr is completely engaged and handpiece's knurled knob is in secure position. 5. Clean attachment or replace attachment.
Inoperative Hand Control	1. Defective internal component 2. Hand Control may not be fully installed. 3. Selected Hand Control on face panel not in use. 4. Internal handpiece switch activating Hand Control may be damaged.	1. Return hand control to DePuy Synthes Power Tools for service. 2. The circumferential groove on handpiece housing will be exposed when the Hand Control is fully installed. 3. Press Hand Control button on face panel to illuminate Hand Control LED 4. Return handpiece to DePuy Synthes Power Tools for service.

Technical Specifications for Devices Compliant with IEC-60601-1-2 3rd Edition

The device complies with the following standards:

IEC 60601-1:2005: + CORR. 1 (2006) + CORR. 2 (2007), EN 60601-1:2006 +A11:2011, ANSI/AAMI ES60601-1:2005, CAN/CSA-C22.2 No.60601-1:08, (3rd edition) Medical electrical equipment part 1: General requirements for basic safety and essential performance;

UL 60601-1:2003 (Revised 2006), IEC 601-1:1988 +A1:1991 + A2:1995; EN60601-1:1990 with A1 and A12:1993, A2:1995 and A13:1996; and CAN/CSA C22.2 No. 601.1-M90, Standard for Medical Electrical Equipment – Part 1 General Requirements for Safety

With regard to electrical shock, fire, mechanical hazards, this ETL classified device conforms to AAMI STD ES60601-1, UL STD 60601-1, and is certified to CSA STD C22.2 No.60601-1 and 601.1

This device complies with applicable EEC directives. Medical Device Directive 93/42/EEC as amended by 2007/47/EC. Machinery Directive 2006/46/EC.

Environmental Conditions	
Temperature	Operating: 18°-30°C (65°-85°F)
	Transportation and Storage: -29° – +50°C (-20° - +120°F)
Relative Humidity	Operating: 30% to 70%
	Transportation and Storage: 10% to 90%
Atmospheric Pressure	Operating: 0.7–1.06bar
	Transportation and Storage: Not applicable

Console (SC2100, SC2101, and SC2102) Specifications
Primary: 100-240VAC, 50/60 Hz, 250 VA
Class I: Protective Earth
Fluid Ingress Protection: IPX0
Type B: Applied Part, Continuous Operation

Foot Pedal (E-FP, E-FP-DIR, E-FP-DIR/IRR) Specifications
Size: 26.7 cm x 16.5 cm x 14.7 cm (10.5 in x 6.5 in x 8.5 in)
Weight: 2.0 kg (4.41 lbs)
Cord: 3.66 m (12 ft) in length
Fluid Ingress Protection: IPX8

Table 1: Emission

The Anspach High Speed Electric System is composed of the eMAX 2 Plus Handpiece model: EMAX2PLUS; Electric Consoles models SC2100, SC2101, & SC2102; Foot Control models E-FP, E-FP-DIR, E-FP-DIR/IRR, EPLUS-FP, EPLUSFPNS; Hand Controls EMAX2-HC and E-HC.

Guidance and manufacturer's declaration – electromagnetic emissions		
The Anspach High Speed Electric System is intended for use in the electromagnetic environment specified below. The customer or the user of the Anspach High Speed Electric System should assure that it is used in such an environment.		
Emission test	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group 1	The Anspach High Speed Electric System uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment. RF emissions.
RF emissions CISPR 11	Class A	The Anspach High Speed Electric System is suitable for use in all establishments, other than domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations / flicker emissions IEC 61000-3-3	Complies	

Table 2: Immunity (all devices)

Guidance and manufacturer's declaration – electromagnetic immunity			
The Anspach High Speed Electric System is intended for use in the electromagnetic environment specified below. The customer or the user of Anspach High Speed Electric System should assure that it is used in such an environment.			
Immunity test standard	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±6 kV contact ±8 kV air	±6 kV contact ±8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%
Electrical fast transient / burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/output lines	±2 kV for power supply lines ±1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1 kV line to line ±2 kV line to earth	±1 kV line to line ±2 kV line to earth	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply lines IEC 61000-4-11	<5% U_T (>95% dip in U_T) for 0.5 cycle) 40% U_T (60% dip in U_T) for 5 cycle) 70% U_T (30% dip in U_T) for 25 cycle) <5% U_T (>95% dip in U_T) for 5 sec)	<5% U_T (>95% dip in U_T) for 0.5 cycle) 40% U_T (60% dip in U_T) for 5 cycle) 70% U_T (30% dip in U_T) for 25 cycle) <5% U_T (>95% dip in U_T) for 5 sec)	Mains power quality should be that of a typical commercial or hospital environment.
Note: U_T is the a.c. mains voltage prior to application of the test level.			
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

Table 3: Immunity (not life-supporting devices)

Guidance and manufacturer's declaration – electromagnetic immunity			
The Anspach High Speed Electric System is intended for use in the electromagnetic environment specified below. The customer or the user of the Anspach High Speed Electric System should assure that it is used in such an environment.			
Immunity test standard	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	3 Vrms	Portable and mobile RF communications equipment should be used no closer to any part of the Anspach High Speed Electric System, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d = 1.2\sqrt{P}$ $d = 1.2\sqrt{P}$ 80 MHz to 800 MHz $d = 2.3\sqrt{P}$ 800 MHz to 2.7 GHz Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, ^a should be less than the compliance level in each frequency range. ^b Interference may occur in the vicinity of equipment marked with the following symbol: 
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	
Note 1: At 80 MHz and 800 MHz, the higher frequency range applies. Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			
a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Anspach High Speed Electric System is used exceeds the applicable RF compliance level above, the Anspach High Speed Electric System should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the device containing the Anspach High Speed Electric System.			
b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.			

Table 4: Recommended separation distances (not life-supporting devices)

Recommended separation distances between portable and mobile RF communications equipment and the Anspach High Speed Electric System			
The Anspach High Speed Electric System is intended for use in the electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the Anspach High Speed Electric System can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Anspach High Speed Electric System as recommended below, according to the maximum output power of the communication equipment.			
Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter		
	150 kHz to 80 MHz $d = 1.2\sqrt{P}$	80 MHz to 800 MHz $d = 1.2\sqrt{P}$	800 MHz to 2.5 GHz $d = 2.3\sqrt{P}$
0.01	12 cm	12 cm	23 cm
0.1	38 cm	38 cm	73 cm
1	1.2 m	1.2 m	2.3 m
10	3.8 m	3.8 m	7.3 m
100	12 m	12 m	23 m
For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.			
Note 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.			
Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			

Technical Specifications for Devices Compliant with IEC-60601-1-2 4th Edition

The device complies with the following standards:

AAMI ES60601-1 Issued: 2006/03/09 Medical Electrical Equipment - Part 1: General Requirements for Basic Safety and Essential Performance, Amendment 1 – 2012

CSA C22.2#60601-1 Issued: 2014/03/01 Ed: 3 Medical Electrical Equipment - Part 1: General Requirements for Basic Safety and Essential Performance

IEC 60601-1-6 Issued: 2013/10/29 Ed: 3 Medical Electrical Equipment - Part 1-6: General Requirements for Basic Safety and Essential Performance - Collateral Standard: Usability; Amendment 1

IEC 60601-1: 2005 + CORR. 1:2006 + CORR. 2:2007 + AM1:2012 (or IEC 60601-1: 2012 reprint) EN 60601-1: 2005 + A1:2013, IEC 60601-1 / EN 60601-1, Medical electrical equipment Part 1: General requirements for basic safety and essential performance (CB Scheme)

IEC 62366 Medical devices – Application of usability engineering to medical devices

IEC 62366: 2007 (First Edition) + A1: 2014 (CB Scheme)

IEC 60601-1-6:2010 (Third Edition) + A1:2013 for use in conjunction with IEC62366:2007 (First Edition) + A1:2014 and IEC 60601-1:2005 (Third Edition) + Corr.1 (2006) + Corr.2 (2007) + A1: 2012 or equivalent consolidated version IEC 60601-1:2012 (Ad.3.1), Medical electrical equipment Part 1-6 General requirements for safety - Collateral Standard: Usability

IEC 60601-1-2:2014 (Edition 4) Collateral Standard: Electromagnetic Compatibility- Medical Electrical Equipment.

This device complies with applicable EEC directives. Medical Device Directive 93/42/EEC as amended by 2007/47/EC. Machinery Directive 2006/46/EC.

Environmental Conditions	
Temperature	Operating: 18°-30°C (65°-85°F)
	Transportation and Storage: -29° – +50°C (-20° - +120°F)
Relative Humidity	Operating: 30% to 70%
	Transportation and Storage: 10% to 90%
Atmospheric Pressure	Operating: 0.7–1.06bar
	Transportation and Storage: Not applicable
Console (SC2100, SC2101, and SC2102) Specifications	
Primary: 100-240VAC, 50/60 Hz, 250 VA	
Class I: Protective Earth	
Fluid Ingress Protection: IPX0	
Type B: Applied Part, Continuous Operation	
Foot Pedal (E-FP, E-FP-DIR, E-FP-DIR/IRR) Specifications	
Size: 26.7 cm x 16.5 cm x 14.7 cm (10.5 in x 6.5 in x 8.5 in)	
Weight: 2.0 kg (4.41 lbs)	
Cord: 3.66 m (12 ft) in length	
Fluid Ingress Protection: IPX8	

Table 1: Emission

The Anspach High Speed Electric System is composed of the eMax 2 Plus Handpiece model: EMAX2PLUS; Electric Consoles models SC2100, SC2101, & SC2102; Foot Control models E-FP, E-FP-DIR, E-FP-DIR/IRR, EPLUS-FP, EPLUS-FP-NS; Hand Controls EMAX2-HC and E-HC.

Guidance and manufacturer's declaration – electromagnetic emissions		
The Anspach High Speed Electric System is intended for use in the electromagnetic environment specified below. The customer or the user of the Anspach High Speed Electric System should assure that it is used in such an environment.		
Emission test	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group 1	The Anspach High Speed Electric System uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment. RF emissions.
RF emissions CISPR 11	Class A	The Anspach High Speed Electric System is suitable for use in all establishments, other than domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations / flicker emissions IEC 61000-3-3	Complies	

Table 2: Immunity (all devices)

Guidance and manufacturer's declaration – electromagnetic immunity			
The Anspach High Speed Electric System is intended for use in the electromagnetic environment specified below. The customer or the user of Anspach High Speed Electric System should assure that it is used in such an environment.			
Immunity test standard	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV contact ±15 kV air	±8 kV contact ±15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%
Electrical fast transient / burst IEC 61000-4-4	±2 kV for AC lines ±1 kV for I/O lines 100 kHz Pulse Repetition Frequency	±2 kV for AC lines ±1 kV for I/O lines 100 kHz Pulse Repetition Frequency	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±0.5 kV line to line ±0.5, 1, 2 kV line to earth	±0.5, 1 kV line to line ±0.5, 1, 2 kV line to earth	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply lines IEC 61000-4-11	0% U _T (100% dip in U _T) for 0.5 periods at 0°, 45°, 90°, 135°, 180°, 225°, 270°, 315° 0% U _T (100% dip in U _T) for 1 period 70% U _T (30% dip in U _T) for 25 periods (50 Hz) 30 periods (60 Hz) Interruptions 0% U _T (100% dip in U _T) for 250 periods (50 Hz) 320 periods (60 Hz)	0% U _T (100% dip in U _T) for 0.5 periods at 0°, 45°, 90°, 135°, 180°, 225°, 270°, 315° 0% U _T (100% dip in U _T) for 1 period 70% U _T (30% dip in U _T) for 25 periods (50 Hz) 30 periods (60 Hz) Interruptions 0% U _T (100% dip in U _T) for 250 periods (50 Hz) 320 periods (60 Hz)	Mains power quality should be that of a typical commercial or hospital environment. If the user of the EMAX2PLUS system requires continued operation during power mains interruptions, it is recommended that the system be powered from an uninterruptible power supply or a battery.
Note: U _T is the AC mains voltage prior to application of the test level.			
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

Table 3: Immunity (not life-supporting devices)

Guidance and manufacturer's declaration – electromagnetic immunity			
The Anspach High Speed Electric System is intended for use in the electromagnetic environment specified below. The customer or the user of the Anspach High Speed Electric System should assure that it is used in such an environment.			
Immunity test standard	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz 6 Vrms in ISM bands 150 kHz to 80 MHz	3 Vrms 150 kHz to 80 MHz 6 Vrms in ISM bands 150 kHz to 80 MHz	Portable and mobile RF communications equipment should be used no closer to any part of the Anspach High Speed Electric System, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d = 1.2\sqrt{P}$ $d = 1.2\sqrt{P}$ 80 MHz to 800 MHz $d = 2.3\sqrt{P}$ 800 MHz to 2.7 GHz Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, ^a should be less than the compliance level in each frequency range. ^b Interference may occur in the vicinity of equipment marked with the following symbol: 
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz 80% at 1 kHz (AM modulation)	3 V/m 80 MHz to 2.5 GHz 80% at 1 kHz (AM modulation)	
Note 1: At 80 MHz and 800 MHz, the higher frequency range applies.			
Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			
a	Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Anspach High Speed Electric System is used exceeds the applicable RF compliance level above, the Anspach High Speed Electric System should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the device containing the Anspach High Speed Electric System.		
b	Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.		

EMAX2PLUS System Enclosure Port Immunity to RF wireless communications equipment (IEC 61000-4-3)						
Test frequency (MHz)	Band ^a (MHz)	Service ^a	Modulation ^b	Maximum power (W)	Distance (m)	IMMUNITY TEST LEVEL (V/m)
385	380–390	TETRA 400	Pulse modulation ^b 18 Hz	1.8	0.3	27
450	430–470	GMRS 460; FRS 460	FM: ±5 kHz deviation 1 kHz sine	2	0.3	28
710	704–787	LTE Band 13, 17	Pulse modulation ^b 217 Hz	0.2	0.3	9
745						
780						
810	800–960	GSM 800/900; TETRA 800; iDEN 820; CDMA 850; LTE Band 5	Pulse modulation ^b 18 Hz	2	0.3	28
870						
930						
1,720	1,700–1,990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; UMTS	Pulse modulation ^b 217 Hz	2	0.3	28
1,845						
1,970						
2,450	2,400–2,570	Bluetooth; WLAN 802.11b/g/n; RFID 2450; LTE Band 7	Pulse modulation ^b 217 Hz	2	0.3	28
5,240	5,100–5,800	WLAN 802.11a/n	Pulse modulation ^b 217 Hz	0.2	0.3	9
5,500						
5,785						
Note: If necessary to achieve the IMMUNITY TEST LEVEL, the distance between the transmitting antenna and the EMAX2PLUS System may be reduced to 1 m. The 1 m test distance is permitted by IEC 61000-4-3.						
a For some services, only the uplink frequencies are included.						
b The carrier shall be modulated using a 50% duty cycle square wave signal.						
c As an alternative to FM modulation, 50% pulse modulation at 18 Hz may be used because while it does not represent actual modulation, it would be worst case.						

Table 4: Recommended separation distances (not life-supporting devices)

Recommended separation distances between portable and mobile RF communications equipment and the Anspach High Speed Electric System			
The Anspach High Speed Electric System is intended for use in the electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the Anspach High Speed Electric System can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Anspach High Speed Electric System as recommended below, according to the maximum output power of the communication equipment.			
Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter		
	150 kHz to 80 MHz $d = 1.2\sqrt{P}$	80 MHz to 800 MHz $d = 1.2\sqrt{P}$	800 MHz to 2.5 GHz $d = 2.3\sqrt{P}$
0.01	12 cm	12 cm	23 cm
0.1	38 cm	38 cm	73 cm
1	1.2 m	1.2 m	2.3 m
10	3.8 m	3.8 m	7.3 m
100	12 m	12 m	23 m
For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.			
Note 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.			
Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			

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© DePuy Synthes 2019. All rights reserved.
The Anspach Effort, Inc.
Synthes GmbH