English

Otologic Curved Micro (OCM) Attachment and Otologic Curved Micro (OCM) Burr Support Sleeves**List of Devices**

OCM	OCM-15SD	OCM-06SD	OCM-1SD
OCM-05SD	OCM-1S	OCM-07SD	OCM-2SD
OCM-06SD	OCM-1SD	OCM-08SD	OCM-2SD
OCM-07SD	OCM-2SD	OCM-15SB	
OCM-08SD	OCM-2SD	OCM-15SD	
OCM-09SD	OCM-1SD	OCM-15SD	

Device Description

The Otologic Curved Micro (OCM) attachment connects to XMax® Pneumatic and eMax 2 and eMax 2 Plus Electric Systems and accepts OCM Burr Support Sleeves which are available in various cutting burs.

The OCM Burr Support Sleeve (OCS) consists of the cutter and an outer sleeve that protects and supports the rotating shaft of the cutting burr.

Side View

Indications for Use

When used with the ANSPACH® Systems, the OCM Attachment and OCM Burr Support Sleeves are intended for cutting and shaping bone primarily in otology procedures such as colectomies.

US Sterilization Instructions

Option	Exposure Temperature	Exposure Time	Minimum Dry Time
1. Pre-Vacuum (wrapped or unwrapped)	132°C	4 minutes	20 minutes
2. Gravity Air Displacement (wrapped or unwrapped)	132°C	15 minutes	15 minutes

Warning:

- To first use the OCM Attachment must be processed as per the included Cleaning and Sterilization instructions.
- It is not intended for gross bone removal. Heavy loading and / or lack of irrigation while performing bone dissection may result in burr fracture.
- Prior to use visually inspect the instruments. Do not use the instruments if damaged.
- To avoid overheating this device, employ a duty cycle of 10 seconds ON and 10 seconds OFF for an unlimited number of cycles.

Note: When following these instructions hold handpiece and attachments with distal end pointing away from user. OCIM Installation for XMax® and eMax Handpieces

OCM Installation for XMax® and eMax Handpieces

- Ensure handpiece is in "RUN" or "SECURE" position.
- With the OCM Attachment's drive shaft into distal end of handpiece, push and rotate to right until it stops. Listen for a click as handpiece knurl knobs move up to seat next to OCM Attachment.
- Gently twist back and forth to ensure OCM Attachment is locked into position.

OCM Burr Support Sleeve Assembly

- Line up arrow located on the distal end of OCM Attachment with flat end of OCM Burr, push until fully seated. Listen for a click. (Illustration 1)
- Gently pull OCM Burr to ensure it is locked into position.


Illustration 1 (Top View)

Curved Burr Support Sleeve Disassembly

- Remove OCM Burr from distal end of OCM Attachment by grasping OCM Burr at groove(s). (Illustration 1)
- Retract the OCM Burr Support Sleeve as bihazardous waste and dispose of accordingly.

OCM Burrs

Warning: Disposable. Single use only. Reuse and/or re-sterilization of dissection tools could result in reduced performance, material degradation and/or patient contamination.

OCM Disassembly for XMax® and eMax Handpieces

- Move handpiece knurl knob to "SAFE" position.
- Rotate OCM Attachment to left approximately one quarter turn, and remove from distal end of handpiece housing.

Manual Cleaning & Lubrication**Warnings:**

- Device must be cleaned as soon as possible after use to prevent drying of blood, tissue, other biological debris and contaminants on the device.
- National regulations must be observed.
- You must follow standard sterilization procedures.
- You must follow the manufacturer's recommendations for any detergents, disinfectants and washing equipment used.
- Do not reprocess OCIM Burrs. OCIM Burrs are disposable and SINGLE USE ONLY. Reuse and/or re-sterilization of dissection tools could result in reduced performance, material degradation and/or patient contamination.

Cautions:

- Do not use an alkaline detergent with a pH greater than 11. Alkaline detergents with higher pH attack greases and seals which can increase wear and cause device to malfunction.
- Device must not be cleaned using ultrasound equipment.
- Deionized (deionized or purified) water must be used during final rinse to prevent corrosion and spots.
- Avoid exposure to hypochlorite solutions (i.e. bleach), as these will promote corrosion.

Note: Use of alkaline detergents can cause device color to fade but this does not impair function.

Operating Room (OR) Personnel Instructions

- After completion of procedure, remove external debris by wiping attachments and handpiece.
- Dispose of open or used dissection tools.

Reprocessing Personnel Instructions

- Assemble all necessary supplies including Attachment Cleaning Brush (ACB), sponge, lint-free cloth and soft bristle brush to assist in removing gross soil.
- Remove OCM Attachment and OCM Burr from handpiece.

Rinse

- Rinse OCM Attachment under cold deionized (DI) or purified water (PURW) for a minimum of 2 minutes.
- Manipulate any moving parts such as release sleeves under running cold deionized (DI) or purified water (PURW) to loosen or remove gross soil.
- Use a sponge, soft lint-free cloth or soft bristle brush to assist in removing gross soil.

Caution: Do not use high pressure water or air. Do not rinse OCM Attachment with saline solution.

Clean

- Prepare cleaning solution using alkaline detergent (non-chlorinated) or neutral pH enzymatic cleaner or EnzymePro with deionized water for correct temperature, water quality (i.e. pH neutral) and product concentration.
- Fully immerse OCM Attachment in cleaning solution as per manufacturer's instructions and manually agitate for 2 minutes.
- Use a sponge, soft lint-free cloth or soft bristle brush to remove all visible soil.
- Gently use ACB to clean cannulations of OCM Attachment while immersed in cleaning solution.

Lubricate

Caution: This step is for the OCM attachment only. Do not lubricate handpieces.

Option 1: 05.001.071 Lubricant for Ansapach Systems, 110 ml

Warning: ACB is single use and should only be used to clean one set of attachment(s).

Rinse Thoroughly

- Rinse device thoroughly using running deionized (DI) or purified water (PURW) for a minimum of 2 minutes.
- Use a syringe or pipette filled with hot deionized (DI) or purified water (PURW) to remove any excess liquid.

End of Life

Return all equipment for proper disposal.

Recommended Manufacturer Inspection Interval

It is recommended that the equipment be returned to the DePuy Synthes Power Tools Customer Service Department regarding TSE/CJD contamination at a minimum every 12 months to perform a full product inspection.

Cautions: Do not use high pressure water or air. Do not rinse OCM Attachment with saline solution.

Inspect

- Visionally inspect OCM Attachment for visible soil. Repeat Steps 2 thru 5 if visible soil remains.

Warning: Transmissible Spongiform Encephalopathies (TSE)

DePuy Synthes Power Tools will not authorize or accept the return of products that are damaged or contaminated with a patient's body fluids who is suspected or confirmed to have a transmissible Spongiform Encephalopathy (TSE/CJD) disease. DePuy Synthes Power Tools recommends that all Ansapach products used on a patient confirmed with a TSE/CJD diagnosis be incinerated. Ansapach OCIM Burrs used on a patient suspected of TSE/CJD diagnosis must be incinerated.

Contact your Sales Representative for replacement of product incinerated under this policy or for temporary equipment while original equipment is quarantined. Contact the DePuy Synthes Power Tools Customer Service Department regarding TSE/CJD contamination for additional information.

Option 2: Non-Silicone Based Medical Lubricant

- Prepare a lubricating solution of instrument milk (non-silicone based medical lubricant) per the manufacturer's directions.
- Fully immerse OCM Attachment in lubricating solution as per manufacturer's instructions, at room temperature, and agitate for 15 seconds.
- Remove OCM Attachment and allow it to drain completely until no visible droplets are coming from it.

Caution: Do not rinse out instrument milk (lubricant). Do not apply mineral oil or other lubricants which may cause OCM Attachment to overheat.

Note: The clinical processing instructions provided have been validated by DePuy Synthes Power Tools for preparing a non-sterile medical device; this instruction is provided in accordance with ISO 17664 and ANSI/AAMI ST81.

Sterilization

Before sterilization:

- Ensure sterilization equipment is in proper working order as specified by manufacturer.
- Ensure equipment manufacturer's instructions are properly employed by trained and qualified personnel.

For patent information go to www.DePuy.com/patentmarketing.

© DePuy Synthes Power Tools, Division of DJO 2014.

Deutsch-German Otologisch gebogenes Mikro-Zubehör (OCM) und otologisch gebogene Mikro-Stützhülsen (OCM)**English****Micro Curved Micro (OCM) Attachment and Micro Curved Micro (OCM) Burr Support Sleeves****List of Devices**

OCM	OCM-15SD	OCM-06SD	OCM-1SD
OCM-05SD	OCM-1S	OCM-07SD	OCM-2SD
OCM-06SD	OCM-1SD	OCM-08SD	OCM-2SD
OCM-07SD	OCM-2SD	OCM-15SB	
OCM-08SD	OCM-2SD	OCM-15SD	
OCM-09SD	OCM-1SD	OCM-15SD	

The following instructions apply to the devices listed below.

Assure actual cycle employed has been properly sterilized for the device(s)/load configuration being processed and appropriate sterilization indicator device are included for each process and cycle.

Warning: Device must be cleaned prior to sterilization using a validated cleaning method from instructions provided with device in accordance with ISO 17664 and ANSI/AAMI ST81.

1. Packaging for Sterilization:

- When wrapping the sterilization container a lid must be used.
- In US double wrap basket or instruments in PTFE covered sterilization wrap using standard wrap techniques described in ANSI/AAMI ST79 (for example: Simultaneous double wrapping or sequential wrapping).
- Outside US double wrap basket or instruments in accordance with local procedures, using standard wrap and wrapping techniques such as those described in ANSI/AAMI ST79.

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