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**Óleo de manutenção Synthes, 40 ml,  
para Electric Pen Drive e Air Pen Drive**

Instruments and implants  
approved by the AO Foundation

 **DePuy Synthes**

**POWER TOOLS**

COMPANIES OF *Johnson & Johnson*

## 05.001.095

### Synthes Maintenance Oil, 40 ml, for Electric Pen Drive and Air Pen Drive

#### Intended Use

Synthes Maintenance Oil is intended to lubricate surgical instruments, its components and for their maintenance.

#### General Information

Technically, Synthes Maintenance Oil may support the mechanical properties of surgical instruments and its components, may help to reduce friction and corrosion and therefore heat generation and may help to improve the longevity of the instruments. The oil film is not impairing effective sterilization. Synthes Maintenance Oil is a silicone free product.

#### Safety Instructions

The Synthes Maintenance Oil is only to be used for lubricating Synthes Power Tool Systems (Electric Pen Drive and Air Pen Drive including its components) and instruments after careful consultation of the instructions for use of the corresponding surgical instrument or power tool. Do not use Synthes Maintenance Oil for Anspach Power Tools.

Synthes Maintenance Oil is just to be applied prior to sterilization.


Synthes Maintenance Oil is designed for use by the reprocessing department or the Synthes service center and production personnel.

Do not use Synthes Maintenance Oil after the indicated expiry date.

To learn more about hazards, first aid measures, handling and storage, read the attached Material & Safety Data Sheet (MSDS).

#### Application

1. Check expiry date and remove lid.
2. Remove the safety seal of the dispenser with scissors.
3. For best handling and precise application we recommend the use of a dispense tip (Nordson, green,  $d=0.84$  mm).
4. Consult the instructions for use of the corresponding power tool system or surgical instrument for proper use.

### Precautions:

- Store in a cool, dry and dark place
- Keep away from direct sunlight
- Read, follow and keep the instructions for use

### Explanation of Symbols

Caution

Consult the Instructions for Use before operating the device

Product is delivered in an unsterile condition

Keep away from direct sunlight

Indicates the expiry date

05.001.095

**Synthes Maintenance Oil, 40 ml,  
for Electric Pen Drive and Air Pen Drive**

**Material Safety Data Sheet (MSDS) SE\_492001 AB**

Synthes Maintenance Oil, 40 ml

**1. Identification of the substance/mixture  
and of the company/ undertaking**

**1.1 Product identifier**

Substance name / Trade name:

**Synthes Maintenance Oil**

**1.2 Relevant identified usage of the substance or mixture  
and usage that is not recommended**

Indication: Usage as lubricant

Contraindication: None

**1.3 Details of the supplier of the safety data sheet**

**Manufacturer/Supplier**

Synthes GmbH

**Address / PO Box**

Eimattstr. 3

**National code / Postal code / City**

CH-4436 Oberdorf

**Contact for technical information**

msds@synthes.com

**Telephone / Telefax**

+41 61 965 61 11

+41 61 965 66 00

**1.4 Emergency telephone number**

+41 44 251 51 51 TOX-Center Zurich

**2. Hazards identification**

**2.1 Classification of the substance or mixture according to  
Regulation (EC) No. 1272/2008, Annex VII**

The product is not required to be labeled according to the relevant rules and standards.

## **Classification according to Directive 67/548/EEC or Directive 1999/45/EC**

The product is not required to be labeled according to the relevant rules and standards.

### **3. Composition/Information on ingredients**

#### **3.1 Mixtures**

Preparation of glycerin-fatty acid

### **4. First Aid Measures**

#### **4.1 Description of the First Aid Measures**

##### **Following inhalation**

Remove victim to fresh air. In case of breathing arrest: artificial respiration.

##### **Following skin contact**

Remove / Take off all contaminated clothing. Rinse skin with soap and water, then rinse with plenty of water.

##### **Following eye contact**

Rinse with plenty of water. If eye irritation persists, consult a doctor and take along this safety data sheet.

### **5. Fire fighting measures**

#### **5.1 Extinguishing media**

Suitable extinguishing media:

Dry powder extinguisher with ABC- or BC-powder, carbon dioxide extinguisher, foam extinguisher

Unsuitable extinguishing media: Water jet

#### **5.2 Advice for fire fighters**

A fire may release toxic gases (carbon monoxide, carbon dioxide).

Wear respiratory protection when fighting the fire.

05.001.095

## Synthes Maintenance Oil, 40 ml, for Electric Pen Drive and Air Pen Drive

### 6. Accidental release measures

#### 6.1 Personal precautions, protective equipment and actions to be taken in case of emergency

Keep area well ventilated. Remove residuals of liquids with cleaning cloth or oil absorbing agent. Danger of slipping on spilled product.

Avoid fire or naked lights.

Do not smoke.

### 7. Handling and storage

#### 7.1 Precautions for safe handling

Store dry at room temperature. Storage area should be well ventilated.

Local exhaust ventilation at the workplace is necessary.

#### 7.2 Conditions for safe storage, including any incompatibilities

##### Information and storage conditions

Do not store in escape routes, work rooms and in their close vicinity.

The substance should not be stored together with substances which are liable to cause dangerous chemical reactions (e.g. strong acids and alkalis, oxidizing substances).

Do not eat, drink or smoke if handling the product.

### 8. Exposure controls / personal protection

#### 8.1 Parameters to be monitored

##### 8.1.1 Limits for exposure at the workplace and/or biological limits

##### Maximum allowable concentration (MAC) Germany, SUVA 1903 CH

Comply with the limit values for cooling lubricants, 10 mg/m<sup>3</sup>  
(BiA-Testing 6750)

## 8.2 Limiting and monitoring exposure

### 8.2.1 Individual protection measures – personal protection

The general rules handling chemical products should be observed.

#### Eye / Face protection

Wear eye protection if handling the product.

#### Skin protection

In case of long, excessive or repeated skin contact wear protective gloves (e.g. made of butyl rubber suitable for applications up to 8 h).

#### Work hygiene

Avoid skin contact with the product.

Avoid inhaling aerosoles.

Remove contaminated clothing.

## 9. Physical and chemical properties

### 9.1 Details on the basic physical and chemical properties

Appearance:

– State of aggregation: Liquid

– Colour: Colourless to slightly yellow

Odour: Characteristic

pH-Value: Not applicable

Initial boiling point and boiling range: > 200 °C (392 °F)

Flashing point: > 230 °C (446 °F)

Flammability (solid, gaseous): Flammable

Vapour pressure: < 0.1 mbar at 20 °C, DIN 51519

Solubility(ies): Insoluble

Ignition temperature: > 300 °C (572 °F), DIN 51794

Viscosity: 15 mm<sup>2</sup>/s at 40 °C (104 °F), DIN EN ISO 3104

Explosive properties: Gas-air mixture is potentially explosive

05.001.095

## Synthes Maintenance Oil, 40 ml, for Electric Pen Drive and Air Pen Drive

### 10. Stability and reactivity

#### 10.1 Thermal decomposition

> 220 °C / 428 °F

#### 10.2 Hazardous products of decomposition

Unknown if stored and handled appropriately. In case of improper use (e.g. high temperature) or fire harmful vapours may be released to the environment.

#### 10.3 Hazardous reactions

With strong oxidants (e.g. liquid or condensed oxygen), with very strong acids and very strong bases, high temperatures, naked flames and other ignition sources.

### 11. Toxicological information

#### 11.1 Details on toxicological effects

##### Details on the base oil

Oral toxicity LD<sub>50</sub> > 20 000 mg/kg (mouse)

##### Irritant effect on the skin

Not sensitizing

### 12. Ecological information

The agent is classified as a substance not hazardous to water, without WHC, ID-Number 760. It is biodegradable. Degradable according to ISO 10708 > 60 % after 28 days.

### 13. Disposal considerations

Completely emptied cans may be disposed of with regular household garbage. Not completely emptied cans are subject to country-specific regulations concerning the disposal of special wastes.

Special waste code: 13 02 07, for oily absorbent material:  
15 02 02

### 14. Transport information

The product is not classified as hazardous material according to the transport regulations:

ADR/RID, SDR/RSD, IMDG/IMO, IATA/ICAO

## 15. Regulatory information

### 15.1 Safety, health and environmental provisions / specific legal provisions on the substance or mixture

The product is not subject to classification according to the relevant regulations. It does not contain any substances listed in the Swiss VOC Regulation or the 31. BImSchV.

## 16. Other information

### Legend

<b>ADR</b>	European Agreement on the International Transport of Dangerous Goods by Road
<b>BImSchV</b>	German Federal Emission Protection Ordinance
<b>CAS</b>	Chemical Abstracts Service
<b>ChemV</b>	Swiss Chemicals Regulation
<b>DIN</b>	German Institute for Standardization
<b>EC</b>	Effective Concentration
<b>EC</b>	European Community
<b>EN</b>	European Standard
<b>GHS</b>	Globally Harmonized System
<b>IATA-DGR</b>	International Air Transport Association – Dangerous Goods Regulations
<b>IBC Code</b>	International Code for Building and Equipment of ships for the carriage in bulk of dangerous chemicals
<b>ICAO-TI</b>	International Civil Aviation Organization – Technical Instructions
<b>IMDG Code</b>	International Maritime Code for Dangerous Goods
<b>ISO</b>	International Organization for Standardization
<b>IUCLID</b>	International Uniform Chemical Information Database
<b>LC</b>	Lethal Concentration
<b>LD</b>	Lethal Dose
<b>log Kow</b>	Distribution Coefficient between octanol and water
<b>MARPOL</b>	Maritime Pollution Convention
<b>OECD</b>	Organisation for Economic Co-operation and Development
<b>PBT</b>	Persistent, bioaccumulative, toxic
<b>RID</b>	Regulation on the International Carriage of Dangerous Goods by Rail
<b>TRGS</b>	Technical Rules for Hazardous Substances
<b>UN</b>	United Nations
<b>VOC</b>	Volatile Organic Compounds
<b>vPvB</b>	very persistent and very bioaccumulative
<b>VwVwS</b>	German Regulation on Substances Hazardous to Water
<b>WHC</b>	Water Hazard Class

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60099549 SE\_484988 AA 04/2013 40100164 © Synthes, Inc. or its affiliates
Subject to modifications Synthes is a trademark of Synthes, Inc. or its affiliates