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**05.001.098**

**Synthes Maintenance Spray, 400 ml,  
for Electric Pen Drive and Air Pen Drive**

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**Spray de manutenção Synthes, 400 ml,  
para Electric Pen Drive e Air Pen Drive**

Instruments and implants  
approved by the AO Foundation

 **DePuy Synthes**

**POWER TOOLS**

COMPANIES OF *Johnson & Johnson*

**05.001.098**

## **Synthes Maintenance Spray, 400 ml**

### **Intended Use**

The Synthes Maintenance Spray is intended to lubricate surgical instruments, its components and for their maintenance.

### **General Information**

Technically, the Synthes Maintenance Spray may support the mechanical properties of surgical instruments and its components, may help to reduce friction and corrosion and therefore heat generation and may help to improve the longevity of the instruments. The oil film is not impairing effective sterilization. Synthes Maintenance Spray is a silicone free product. Do not use Synthes Maintenance Spray for Anspach Power Tools.

### **Safety Instructions**

The Synthes Maintenance Spray is only to be used for lubricating Synthes Power Tool Systems (Electric Pen Drive and Air Pen Drive including its components) and instruments after careful consultation of the instructions for use of the corresponding surgical instrument or power tool.

Synthes Maintenance Spray is just to be applied prior to sterilization.


The Synthes Maintenance Spray is designed for use by the reprocessing department or the Synthes service center and production personnel.

Do not use Synthes Maintenance Spray after the indicated expiry date.

To learn more about hazards, first aid measures, handling and storage, read the attached Material & Safety Data Sheet (MSDS).

## Application

1. Check expiry date and remove lid
2. Attach the specific adapter for the instrument to lubricate:
  - Electric Pen Drive (EPD with adapter 05.001.101)
  - Air Pen Drive (APD, with adapter 05.001.092)
  - Attachments for EPD and APD (with adapter 05.001.102)
  - Screwdriver 90° (with adapter 05.505.002)
3. Consult the instruction for use of the corresponding power tool system or surgical instrument for proper use.

**Precautions:**

- Product contains extremely flammable gas under pressure; may explode if heated
- Store in a cool, dry and dark place
- Keep away from direct sunlight
- Read, follow and keep the instruction for use
- Wet cleaning under mechanical load may dissolve the etiquette of the Synthes Maintenance Spray

**Explanation of Symbols**

Caution

Consult the Instructions for Use before operating the device

Product is delivered in an unsterile condition

Keep away from direct sunlight

Indicates the expiry date

**Material Safety Data Sheet (MSDS), SE\_492023 AB**

Synthes Maintenance Spray, 400 ml

**1. Identification of the substance/mixture and of the company/undertaking****1.1 Product identifier**

Substance name / Trade name:

**Synthes Maintenance Spray****1.2 Relevant identified uses of the substance or mixture and uses advised against**

Indication: Usage as lubricant

Contraindication: None

**1.3 Details of the supplier of the safety data sheet****Manufacturer/Supplier**

Synthes GmbH

**Address / PO Box**

Eimattstr. 3

**National code / postal code (ZIP) / City**

CH-4436 Oberdorf

**Contact for technical information**

msds@synthes.com

**Telephone/Telefax**

+41 61 965 61 11

+41 61 965 66 00

**1.4 Emergency telephone number**

+41 44 251 51 51 TOX-Center Zurich

**2. Hazards identification****2.1 Classification of the substance or mixture according to Regulation (EC) No. 1272/2008, Annex VII**

Flame, R phrases: R 12, Extremely flammable, R 18, In use, may form flammable/explosive vapour-air mixture S phrases: S 3/9, Keep container in a cool, well-ventilated place, S 16, Keep away from sources of ignition, No smoking, S 33, Take precautionary measures against static discharges, S 51, Use only in well-ventilated areas

05.001.098

## Synthes Maintenance Spray, 400 ml

### Classification according to Directive 67/548/EEC or Directive 1999/45/EC

Flammable Aerosol (H222)

Gases under pressure, Category compressed gas, H280

## 2.2 Label Elements

### Labelling according to Regulation (EC) No 1272/2008 (substances) / Directive 1999/45/EC (mixtures)

#### Pictogram

**Danger:** Contains gas under pressure; may explode if heated

#### Precautionary statements (P-phrases)

Keep out of reach of children. (P102)

Do not spray on an open flame or other ignition source. (P211)

Pressurized container: Do not pierce or burn, even after use. (P251)

Protect from sunlight. Do not expose to temperatures exceeding 50 °C / 122 °F. (P410+P412)

## 3. Composition/Information on ingredients

### 3.1 Mixtures

Preparation of glycerin-fatty acid-ester and compressed gas (propane/butane)

Substance name: **Butane**

EC-No.: 203-448-7 / CAS-No.: 106-97-8

Contingent: 10–25 %

Substance name: **Propane**

EC-No.: 200-827-9 / CAS-No.: 74-98-6

Contingent: 10–25 %

Substance name: **Isobutane**

EC-No.: 200-857-2 / CAS-No.: 75-28-5

Contingent: 2–10 %

Classification according to regulation (EC) Nr. 1272/2008:  
Flammable gases, Category 1; H220  
Gases under pressure H280

(Full text of quoted hazard notes: see section 16)

## **4. First aid measures**

### **4.1 Description of first aid measures**

#### **Following inhalation**

Move victim to fresh air. In case of breathing arrest: artificial respiration.

#### **Following skin contact**

Remove / Take off all contaminated clothing. Rinse skin with soap and water, then rinse with plenty of water.

#### **Following eye contact**

Rinse with plenty of water. If eye irritation persists, consult a physician and take along this safety data sheet.

### **4.2 Most important symptoms and effects, both acute and delayed**

N-Butane, Propane and Isobutane have no irritating effects as gas; if released from aerosol cylinders, cold injuries may occur. Systemic effects are to be expected with very high gas concentrations.

## **5. Firefighting measures**

### **5.1 Extinguishing media**

Suitable extinguishing media:

Dry powder extinguisher with ABC- or BC-powder, carbon dioxide extinguisher, foam extinguisher, water spray

### **5.2 Advice for firefighters**

In case of fire eliminate cylinder from the danger zone or cool with firefighting water. Danger of explosion! Wear respiratory protection when fighting the fire.

## **6. Accidental release measures**

### **6.1 Personal precautions, protective equipment and emergency procedures**

Keep away from open flames, hot surfaces and sources of ignition

## 7. Handling and storage

### 7.1 Precautions for safe handling

Store dry at room temperature. Storage area should be well ventilated.

Local exhaust ventilation at the workplace is necessary.

### 7.2 Conditions for safe storage, including any incompatibilities

#### Information on storage conditions

Do not store in escape routes, work rooms and in their close vicinity.

The substance should not be stored with substances which are liable to cause dangerous chemical reactions (e.g. strong acids and alkalis, oxidizing substances).

**Storage class:** 2 Aerosols

## 8. Exposure controls / personal protection

### 8.1 Parameters to be monitored

#### 8.1.1 Limits for exposure at the workplace and/or biological limits

##### Maximum allowable concentration (MAC) Germany, SUVA 1903 CH

Substance name: **Butane; CAS-No.: 106-97-8:**

1000 ml/m<sup>3</sup>, 2400 mg/m<sup>3</sup>, peak limit: excess factor 4, duration 15 min, mean value; 4 times per shift; interval 1 h category II – Substances with a resorptive effect

Substance name: **Propane; CAS-No.: 74-98-6:**

1000 ml/m<sup>3</sup>, 1800 mg/m<sup>3</sup>, peak limit: excess factor 4, duration 15 min, mean value; 4 times per shift; interval 1 h category II – Substances with a resorptive effect

Substance name: **Isobutane; CAS-No.: 75-28-5:**

1000 ml/m<sup>3</sup>, 2400 mg/m<sup>3</sup>, peak limit: excess factor 4, duration 15 min, mean value; 4 times per shift; interval 1 h category II – Substances with a resorptive effect

## 8.2 Limiting and monitoring exposure

### 8.2.1 Individual protecting measures / personal protective equipment

#### Eye/Face protection

Wear eye protection if handling the product.

#### Skin protection

Not necessary if used appropriately.

In case of long, excessive or repeated skin contact wear protective gloves (e.g. made of butyl rubber suitable for applications up to 8 h).

#### Work hygiene

Avoid skin contact with the liquid phase, frostbites.

Avoid inhaling the gas.

Remove gas contaminated clothing and air them far from ignition sources.

## 9. Physical and chemical properties

### 9.1 Details on the basic physical and chemical properties

Appearance:

– Aggregate state: Aerosol

– Colour: Colourless

Odour: Odourless

pH-value: Not applicable

Initial boiling point / boiling range:

Compressed gas  $-44\text{ }^{\circ}\text{C}$  /  $-47.2\text{ }^{\circ}\text{F}$

Flashing point:

Agent:  $>210\text{ }^{\circ}\text{C}$  /  $410\text{ }^{\circ}\text{F}$ , DIN EN ISO 2592

Compressed gas:  $-97\text{ }^{\circ}\text{C}$  /  $-142.6\text{ }^{\circ}\text{F}$

Flammability (solid, gaseous): flammable

Upper/lower explosion threshold: 12.0 Vol. % / 1.5 Vol. %

Vapour pressure:

Agent:  $<0.1\text{ mbar}$  at  $20\text{ }^{\circ}\text{C}$  /  $68\text{ }^{\circ}\text{F}$ , DIN 51519

Compressed gas: 3.5 bar at  $20\text{ }^{\circ}\text{C}$  /  $68\text{ }^{\circ}\text{F}$

Solubility(ies): Insoluble

Spontaneous ignition temperature:

Agent:  $>300\text{ }^{\circ}\text{C}$  /  $572\text{ }^{\circ}\text{F}$ , DIN 51794

Compressed gas:  $-97\text{ }^{\circ}\text{C}$  /  $-142.6\text{ }^{\circ}\text{F}$

Viscosity: Agent:  $15\text{ mm}^2/\text{s}$  at  $40\text{ }^{\circ}\text{C}$  /  $104\text{ }^{\circ}\text{F}$ , DIN EN ISO 3104

Explosive properties: Gas-air mixture is potentially explosive

**10. Stability and reactivity**

**10.1 Thermal decomposition:** >220 °C / 428 °F

**10.2 Hazardous decomposition products**

Unknown if stored and handled appropriately. In case of improper use (e.g. high temperature) or fire, harmful vapours may be released to the environment.

**10.3 Hazardous reactions**

With strong oxidants (e.g. liquid or condensed oxygen), with very strong acids and bases, high temperatures, naked flames and other ignition sources.

**11. Toxicological information**

**11.1 Details on toxicological effects**

**Details on the base oil**

Oral toxicity LD<sub>50</sub> >20 000 mg/kg (mouse)

**Irritant effect on the skin**

Not sensitizing

The lubricant contained in the aerosol cans is registered at NSF International (NSF/H1-release).

**12. Ecological information**

The agent contained in the aerosol cans is not classified as a substance hazardous to water, without WHC, ID-Number 760. It is bio-degradable.

**13. Disposal considerations**

Completely emptied cans may be disposed of with regular household garbage. Not completely emptied cans underlie country-specific regulations concerning the disposal of special wastes. Special waste code: 16 05 04

## 14. Transport information

### 14.1 Land transport (ADR / RID / GGVS / GGVE)

Class(es): 2

UN-No.: 1950

Packing group: –

Hazard label(s): 2.1

Proper shipping name: Aerosol dispenser, flammable

Limited Quantity (LQ): 1 liter

Transport class: 2

Tunnel restriction code: D

### 14.2 Sea transport (IMDG-Code / GGVSee)

Class(es): 2

UN-No.: 1950

Hazard label(s): 2.1

Proper shipping name: Aerosols, flammable

Packing group: –

Marine Pollutant: No

EMS-Nr.: F-D, S-U

### 14.3 Air transport (ICAO / IATA-DGR)

Class(es): 2

UN / ID-No.: UN1950

Label: 2.1

Proper shipping name: Aerosols, flammable

Packing group: –

### 14.4 Environmental perils

Label for substances which are harmful to the environment

ADR / RID / IMDG-Code / ICAO-TI / IATA-DGR: no

Marine Pollutant: no

## 15. Regulatory information

### 15.1 Safety, health and environmental regulations / legislations specific for the substance or mixture

The safety regulations on the aerosol dispenser are valid.

Labeling of H- and P-statements see chapter 2.

The product contains up to 20 % volatile organic compounds according to the Swiss positive list Annex I and VOC according to TRGS 220, 31 respectively. BimSchV (approx. 50 g per can).

The nominal pressure (test pressure) of the aerosol dispenser is 15 bar.

The optimum pressure is approx. 18 bar.

## 16. Other information

### Legend

<b>ADR</b>	European Agreement on the International Transport of Dangerous Goods by Road
<b>BImSchV</b>	German Federal Emission Protection Ordinance
<b>CAS</b>	Chemical Abstracts Service
<b>ChemV</b>	Swiss Chemicals Regulation
<b>DIN</b>	German Institute for Standardization
<b>EC</b>	Effective Concentration
<b>EC</b>	European Community
<b>EN</b>	European Standard
<b>GHS</b>	Globally Harmonized System
<b>IATA-DGR</b>	International Air Transport Association – Dangerous Goods Regulations
<b>IBC-Code</b>	International Code for Building and Equipment of ships for the Carriage in bulk of dangerous chemicals
<b>ICAO-TI</b>	International Civil Aviation Organization – Technical Instructions
<b>IMDG-Code</b>	International Maritime Code for Dangerous Goods
<b>ISO</b>	International Organization for Standardization
<b>IUCLID</b>	International Uniform Chemical Information Database
<b>LC</b>	Lethal Concentration
<b>LD</b>	Lethal Dose
<b>log Kow</b>	Distribution Coefficient between octanol and water
<b>MARPOL</b>	Maritime Pollution Convention
<b>OECD</b>	Organisation for Economic Co-operation and Development
<b>PBT</b>	Persistent, bioaccumulative, toxic
<b>RID</b>	Regulation on the International Carriage of Dangerous Goods by Rail
<b>TRGS</b>	Technical Rules for Hazardous Substances
<b>UN</b>	United Nations
<b>VOC</b>	Volatile Organic Compounds
<b>vPvB</b>	very persistent and very bioaccumulative
<b>VwVwS</b>	German Regulation on Substances Hazardous to Water
<b>WHC</b>	Water Hazard Class

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60099550 SE\_484989 AA 04/2013 40100163 © Synthes, Inc. or its affiliates
Subject to modifications Synthes is a trademark of Synthes, Inc. or its affiliates