Cleaning and Sterilization

Warning: These instructions apply to the devices listed below. Do not use these instructions to replace any devices that are not listed.

SC200	CRANI-A	XMXX	B-GREEN-R
SC200U	CRANI-A-01	XMXX-H	B-TURO-X
SC210	CRANI-A-R	BLACKMAX-N	B-GREY-ADG
SC210L	CRANI-A-R	BLACKMAX-N-LP	B-TURO-ODA
SC2102	CRANI-L-R	AL-III-30	B-QD8
EMAX2PLUS	CRANI-P	AL-III-DISS-20	B-QD8-S
EMAX2	QD8	AL-III-DISS-30	B-QD11
E-FP	QD8-S	AL-III-RRM-MAT	B-QD11-S
E-FP-DR	QD11	AUTOLUBE-III	B-QD14
E-FP-DR/IRR	QD11-S	AUTOLUBE-III-20	B-QD14-S
EPLUS-FP	QD14	AUTOLUBE-III-DISS	CSR60B
EPLUS-FP-NS	QD14-S	AUTOLUBE-III-DISS	BM-COUPLER
EMAX2HC	CSR60	AUTOLUBE-III-FG7	B-ORANGE-45
EHC	QD8	B-SAGITTAL	B-BLACK
SHORT	QD8	AUTOLUBE-III-MAT	BM-DRIVER
SHORT-HD	COUPLER	AUTOLUBE-III-NK	B-MR-A
MEDIUM	ORANGE-45	B-SILVER	B-MR-M
MEDIUM-HD	ORANGE-90	B-TURQUOISE	B-MR-L
LONG	ORANGE-90	B-SAGITTAL	B-BLACK
LONG-01	SA-JACOBS	B-BLUES	BM-SAGITTAL
LONG-HD	SA-JLATCH	B-TURQ-L	IMR-MOTOR
LONG-S	S-SAW	B-BLACK	IMR-FOOT-CTRL
XL-HD	R-SAW	B-RED	IMR-SHORT
MIAT6	O-SAW	B-GREY/MIAT6	BM-SHORT-HD
		B-GOLD	IMR-CRANI-A
		B-GOLD	IMR-OL

Mechanical Cleaning
For mechanical/automated cleaning instructions see information supplied with appropriate manual and sterilization basket for instruments being cleaned.

Manual Cleaning and Maintenance

- At the point of use the device must be cleaned as soon as possible after use to prevent drying of blood, tissue, other biological debris and contaminants on the device.
- National regulators must be observed.
- You must follow standard hospital procedures.
- You must follow manufacturer’s recommendations for any detergents, disinfectants and washing equipment used.
- Do not reprocess dissection tools. Dissection tools are disposable and Single Use Only. Reprocessing dissection tools may result in reduced performance, material degradation and/or patient contamination.
- Use personal protective equipment.

Caution:

- Rx Only: United States Federal law restricts this device to sale by or on order of a physician or other licensed healthcare provider.
- Do not use any alkaline detergent with a pH greater than 11. Alkaline detergents with higher pH attack greases and seals which can increase wear and cause damage to metalization.
- The use of alkaline detergents can cause color to fade but does not impair function.
- Do not use ultrasonic equipment.
- Do not use hypochlorite solution (i.e. bleach).
- Do not rinse with saline solution.
- Handpieces must not be immersed in cleaning solution but may be sprayed with it.
- Fully demineralized (deionized, distilled, or purified) water must be used during final rinse to prevent corrosion and scale.
- Prior to use, and while cleaning, visually inspect the instruments). Do not use the instrument(s) if damaged.

Preparation

- Assemble all necessary supplies including Attachment Cleaning Brush (ACB), sponge, lint-free cloth and soft bristle brush.
- Remove attachment and dissection tool from handpiece.

- Rinse
 - Rinse handpiece with distal end pointing down under running cold deionized, distilled or purified water for a minimum of 2 minutes.
 - Use a sponge, soft lint-free cloth, or soft bristle brush to assist in removing gross soil.
- Caution: Do not immerse. Do not use high pressure water or air. Do not allow running water to enter distal end of handpiece or console connector (end of cable).
- Clean
 - Prepare cleaning solution using alkaline detergent or neutral pH enzymatic detergent per manufacturer’s directions for correct temperature, water quality (i.e. pH, hardness), and product concentration/dilution.
 - With distal end pointing down spray and wipe using cleaning solution for a minimum of 2 minutes.
 - Use a sponge, soft lint-free cloth, or soft bristle brush to remove all visible soil.
 - Hold handpiece with distal end pointing down, and insert ACB wetted with cleaning solution into distal end of handpiece.

Warning: ACB is single use and should only be used to clean one handpiece and associated attachments.

Caution: Do not immerse. Do not allow running water to enter distal end of handpiece or console connector. Do not insert ACB beyond last stride. Do not force ACB into handpiece.

- Rinse thoroughly
 - With distal end pointing down rinse device thoroughly using running hot deionized, distilled or purified water for a minimum of 2 minutes.
 - Hold handpiece with distal end pointing down and flush inside of distal end using a syringe or pipette filled with hot deionized, distilled or purified water.
- Caution: Do not use high pressure water or air. Do not immerse. Do not allow running water to enter distal end of handpiece or console connector.

- Clean
 - Prepare cleaning solution using alkaline detergent or neutral pH enzymatic detergent per manufacturer’s directions for correct temperature, water quality (i.e. pH, hardness), and product concentration/dilution.
 - With distal end pointing down spray and wipe using cleaning solution for a minimum of 2 minutes.
 - Use a sponge, soft lint-free cloth, or soft bristle brush to remove all visible soil.
 - Hold handpiece with distal end pointing down, and insert ACB wetted with cleaning solution into distal end of handpiece.

Warning: The following sterilization methods and cycles have been validated in accordance with ISO 17665. Do not use any sterilization method or cycle that has not been validated in accordance with ISO 17665.

Option	Exposure Temperature	Exposure Time	Minimum Dry Time
1. Pre-Vacuum (wrapped or unwrapped)	132°C	4 minutes	30 minutes
2. Gravity Air Displacement (wrapped or unwrapped)	132°C	15 minutes	30 minutes

Option	Exposure Temperature	Exposure Time	Minimum Dry Time
1. Pre-Vacuum (wrapped or unwrapped)	134-138°C	4 minutes	30 minutes
2. Pre-Vacuum (wrapped or unwrapped)	132-138°C	3-18 minutes	30 minutes
3. Gravity Air Displacement (wrapped or unwrapped)	132-138°C	15-18 minutes	30 minutes

Caution: Do not use high pressure water or air. Do not rinse attachment with saline solution.

Preparation

- Assemble all necessary supplies including Attachment Cleaning Brush (ACB), sponge, lint-free cloth and soft bristle brush.
- Remove attachment and dissection tool from handpiece.

- Rinse
 - Rinse attachment under running cold deionized, distilled or purified water for a minimum of 2 minutes.
 - Manipulate any moving parts under running cold deionized, distilled or purified water to loosen or remove gross soil.
 - Use a sponge, soft lint-free cloth or soft bristle brush to assist in removing gross soil.

Caution: Do not use high pressure water or air. Do not rinse attachment with saline solution.

- Clean
 - Prepare cleaning solution using alkaline detergent or neutral pH enzymatic detergent per manufacturer’s directions for correct temperature, water quality (i.e. pH, hardness) and product concentration/dilution.
 - Fully immerse attachment in cleaning solution in a suitable container and agitate for a minimum of 2 minutes.
 - Use a sponge, soft lint-free cloth or soft bristle brush to remove all visible soil.
 - Gently use a sponge, soft lint-free cloth or soft bristle brush to remove all visible soil.
 - Manipulate any moving parts such as release sleeves.

Caution: Do not use high pressure water or air. Do not rinse attachment with saline solution.

Warning: ACB is single use and should only be used to clean one set of attachments. Do not use ACB on any other instrument.

Warning: All Instrument System components returned for servicing or repair should be properly cleaned and sterilized as applicable prior to shipping.

Warning: Transmissible Spongiform Encephalopathies (TSE)
DePuY Synthes Power Tools will not authorize or accept the return of products that directly contact patients or are contaminated with a patient’s body fluids who is suspected of, or confirmed to have, a transmissible spongiform encephalopathy (Creutzfeldt-Jakob Disease (TSE/CJD) diagnosis. DePuY Synthes Power Tools recommends that all Anspach products used on a patient confirmed with a TSE/CJD diagnosis be segregated. Anspach dissection tools used on a patient suspected of TSE/CJD diagnosis must be incinerated.

Contact your Sales Representative for replacement of product incriminated under this policy or for temporary equipment while original equipment is quarantined. Contact the DePuY Synthes Power Tools Customer Service Department regarding TSE/CJD contamination or additional information.

Warning: This step is for attachments only. Do not lubricate handpieces.

Warning: 1. 05.001.078, Lubricant for Anspach Systems_110 ml

Caution: Do not use Gravity Air Displacement Steam Autoclave Sterilization Process. Only for use with attachments listed below.

- For the following attachments apply one (1) spray pump spray to the proximal end of the attachment.

CRANI-A	LONG	ADG	SHORT
CRANI-A-01	LONG-01	ODA	SHORT-HD
CRANI-A-R	LONG-HD	MEDIUM	XL-HD
CRANI-L	LONG-S	MEDIUM-HD	
CRANI-R	CRANI-P	MIAT6	

- For the following attachments apply one (1) spray pump spray to the distal end of the attachment and one (1) to the proximal end of the attachment.

CSR60	QD11	SA-JACOBS	O-SAW
DRIVER	QD14	SA-JLATCH	ORANGE-45
QD8	QD11-S	R-SAW	ORANGE-90
QD8-S	QD14-S	R-SAW	ORANGE-90

- Applying excess amounts of lubricant will cause the lubricant to drip from the attachment. Clean off any excess lubricant.

Option 2: Non-Silicone Based Medical Lubricant
Note: Applicable to all attachments in [List of Devices](#) section at beginning of this document.

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