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05.001.078

Lubrificante para sistemas Anspach®, 110 ml

 **DePuy Synthes**
POWER TOOLS
COMPANIES OF *JOHNSON & JOHNSON*

Instruments and implants
approved by the AO Foundation

 **ANSPACH®**

Intended Use

The lubricant for Anspach systems is intended to lubricate Anspach system attachments during reprocessing, as part of their maintenance.

It is not intended for the lubrication of Anspach BlackMax system or iMRI system attachments.

General Information

Lubrication is MANDATORY to ensure the function of the device as intended by the manufacturer.

The lubricant does not impair effective sterilization. The lubricant for Anspach systems is a silicone-free product and is not made with natural rubber latex.

Safety Instructions

The lubricant for Anspach systems is only to be used for lubricating Anspach power tool system attachments after careful consultation of the instructions for use of the corresponding Anspach power tool system. Use this product only as indicated.

The lubricant for Anspach systems is to be applied just prior to sterilization.

The lubricant for Anspach systems is designed for use by the reprocessing department.

Do not use the lubricant for Anspach systems after the indicated expiry date.

To learn more about hazards, first aid measures, handling and storage, read the attached Material Safety Data Sheet (MSDS).

Application

Do not use lubricant for Anspach systems after the indicated expiry date.

1. Remove cover.
2. Remove protecting cap and install the spray pump provided in the package.
3. Remove the spray pump protective cap and attach the dispenser nozzle provided in the package by pushing it down on the spray pump head.
4. Consult the directions for use of the corresponding Anspach power tool system for proper use.

Precautions:

- Keep away from direct sunlight
- Read, follow and keep the instructions for use
- Wet cleaning under mechanical load may dissolve the etiquette of the lubricant for Anspach systems

Explanation of Symbols

Caution

Consult the Instructions for Use before operating the device

non sterile

Product is delivered in a nonsterile condition

Keep away from direct sunlight

Indicates the expiry date

Material Safety Data Sheet (MSDS), SE_502933 AA

Lubricant for Anspach Systems

1. Identification of the substance/mixture and of the company/undertaking

1.1 Product identifier

Substance name/Trade name:

Lubricant for Anspach Systems

1.2 Relevant identified uses of the substance or mixture

Lubricant for medical instruments

1.3 Details of the supplier of the safety data sheet **Manufacturer/Supplier**

Synthes GmbH

Address/PO Box

Eimattstrasse 3

National code/postal code (ZIP)/City

CH-4436 Oberdorf

Contact for technical information

msds@synthes.com

Telephone/Fax

+41 61 965 61 11

+41 61 965 66 00

1.4 Emergency telephone number

+41 44 251 51 51 Tox - Center Zurich

2. Hazards identification

2.1 Classification of the substance or mixture according to Regulation (EC) No. 1272/2008, Annex VII

Asp. Tox. 1; H304

Classification according to Directive 67/548/EEC or Directive 1999/45/EC:

Substance not classified according to Directive 67/548/EEC.

2.2 Label Elements

Labelling according to Regulation (EC) No 1272/2008 (substances) / Directive 1999/45/EC (mixtures)

68037-01-4 (1-Decene, homopolymer, hydrogenated)

Hazard symbol

GHS08

Signal word = Danger

Hazard statements

H304 May be fatal if swallowed and enters airways

Precautionary statements

P301 + P310 IF SWALLOWED: contact POISON CENTRE or doctor immediately

P331 DO NOT induce vomiting.

P405 Store locked up. Keep away from unauthorized persons and children.

P410 Keep away from sunlight.

2.3 Other hazards

PBT assessment

The product is not classified as PBT.

3. Composition/Information on ingredients

3.1 Substances

Main ingredient of the substance:

Substance name: 1-Decene, homopolymer, hydrogenated

EC-No.: –/CAS-No.: 500-183-1 ; 68037-01-4

Contingent: >95%

(See section 16 for the full text of the hazard statements declared above)

4. First aid measures

4.1 Description of first aid measures

Remove contaminated clothing immediately.
If health problems occur, consult a physician.

Following inhalation:

Move to fresh air in case of accidental inhalation of fumes from overheating or combustion. Use oxygen or artificial respiration if needed.

Following skin contact:

Following contact with the skin wash off with soap and water.

Following eye contact:

Rinse immediately with plenty of water, including under the eyelids, for at least 15 minutes.

If eye irritation persists, consult a specialist.

Following ingestion:

Drink 1 to 2 glasses of water. DO NOT induce vomiting.

If complaints persist, visit a doctor.

4.2 Most important symptoms and effects, both acute and delayed

Symptoms:

Cases of local necrosis due to delayed onset of pain and tissue damage a few hours after ingestion have been documented.

Effects:

Following ingestion and subsequent vomiting, the product may be aspirated into the lungs and cause chemical pneumonia or asphyxiation.

Indication of any immediate medical attention and special treatment needed:

No data available

5. Fire-fighting measures

5.1 Extinguishing media

Suitable extinguishing media: water, foam, dry chemical, carbon dioxide

Unsuitable extinguishing media: Direct water jet

5.2 Special hazards arising from the chemical or mixture

The following may be released in case of fire: carbon monoxide, carbon dioxide

5.3 Advice for firefighters

Cool endangered containers with water. Self-contained breathing apparatus. Wear protective clothing.

6. Accidental release measures

6.1 Personal precautions, protective equipment and emergency procedures

Refer to protective measures listed in sections 7 and 8. Particular risk of slipping on the leaked / spilled product.

6.2 Environmental precautions

Do not release into the environment.

Do not allow material to enter surface waters or sewers.

6.3 Methods and material for containment and cleaning up

Pick up with inert absorbent material (dry sand, diatomaceous earth; universal binder)

Shovel into suitable containers for disposal.

7. Handling and storage

7.1 Precautions for safe handling

- No special handling instructions required.
- Ensure adequate ventilation.
- Keep away from heat and sources of ignition.
- To avoid thermal decomposition, do not overheat (>200°C).
- Clean and dry equipment before starting work.

7.2 Conditions for safe storage, including any incompatibilities

Details of storage conditions

- No special storage conditions required.
- Keep away from heat and sources of ignition (>200°C).
- Keep in properly labeled containers.
- Keep away from combustible materials.
- Refer to protective measures listed in sections 7 and 8.

Requirements for storage areas and containers

Suitable materials: plastic, glass, metals

8. Exposure controls/personal protection

8.1 Control parameters

No parameters for use available

8.2 Exposure controls:

No data available

Personal protective equipment:

Respiratory protection:	No special precautions required
Eye/face protection:	No special precautions required
Hand protection:	No special precautions required, in case of extensive contact wear protective gloves.

9. Physical and chemical properties

9.1 Information on basic physical and chemical properties

Odour:	mild
Melting point/freezing point:	<-20°C
Initial boiling point/boiling range:	>300°C
Flash point:	>155°C
Pour point:	-66°C
Vapour pressure:	0.013 kPa 20°C
Relative density:	0.83 at 20°C
Solubility(ies):	In water; sparingly soluble
Decomposition temperature:	No data available
Kinematic viscosity [mm ² /s]:	10 cSt at 40°C

10. Stability and reactivity

10.1 Reactivity

No data available

10.2 Chemical stability

No data available

10.3 Conditions to avoid

High temperatures

10.4 Incompatible materials

Oxidizing agents

10.5 Hazardous decomposition products

None known.

11. Toxicological information

11.1 Information on toxicological effects

Acute oral toxicity:

LD₅₀, rat, >5000 mg/kg, method OECD 401

LD₅₀, rat, >5000 mg/kg (polytetrafluoroethylene)

Swallowing can cause the following symptoms: nausea, vomiting, diarrhoea

Acute inhalation toxicity:

LC₅₀, 4 h, rat, >5.2 mg/l, method: OECD 403

Acute dermal toxicity:

LD₅₀, hare, >2000 mg/kg, method: OECD 402

Skin corrosion/irritation:

Rabbits, no skin irritation;

Rabbits, no skin irritation, note: 14 days

Symptoms: erythema

Severe eye injury/irritation

Rabbits, no eye irritation

Eye contact may cause irritation/erythema

Sensitization of respiratory system/skin:

Guinea pigs, does not cause sensitisation in laboratory animals.

Note: skin

Germ cell mutagenicity:

Not mutagenic in the Ames-Test.

Chromosome aberration test in vitro, negative.

Reproductive toxicity

No data available.

Other information:

Following ingestion and subsequent vomiting, the product may be aspirated into the lungs and cause chemical pneumonia or asphyxiation

Eye contact with the product may cause slight irritation.

12. Ecological information

12.1 Toxicity

Fish, *Oncorhynchus mykiss*, LC₅₀, 96 h, >1000 mg/l remarks: based on the available data.

Crustaceans, *Daphnia magna*, EC₅₀, 48 h, >1000 mg/l remarks: based on the available data.

Daphnia magna, 21 days; 125 mg/l, remarks: based on the available data.

Algae: 72 h; EL₅₀; 1000 mg/l; remarks: based on the available data.

Bacteria: no data available.

12.2 Persistence and degradability

47.7 after 28 days

12.3 Bioaccumulative potential

No data available

12.4 Mobility in soil

No data available

12.5 Result of PBT and vPvB assessment

No data available

12.6 Other adverse effects

Ecological injuries are not known or expected during correct use. Do not release into the environment in an uncontrolled manner.

13. Disposal considerations

13.1 Waste treatment methods

Management of contaminated packaging

Completely emptied packaging must be removed and disposed of in accordance with the relevant regulations. Incompletely emptied packaging must be disposed of by the regional disposer.

Special precautions

Can be incinerated, when in compliance with local regulations.

14. Transport information

The product is not declared as a dangerous good.

15. Regulatory information

15.1 Safety, health and environmental regulations/legislation specific for the substance or mixture

EU regulations

Annex I, parts 1 + 2: not stated. See section 10 for details of any decomposition products.

National regulations

Water Hazard Class

Technical Rules for Hazardous Substances (TRGS)

Organic Solvents Ordinance (31st BImSchV)

German Hazardous Incident Ordinance (12th BImSchV)

Technical Instructions on Air Quality Control

16. Other information

References and data sources

Regulations

Dangerous Preparations Directive (1999/45/EC), last amended by Regulation (EC) No. 1907/2006.

Dangerous Substances Directive (67/548/EEC), last amended by Directive 2009/2/EC.

REACH Regulation (EC) No. 1907/2006, last amended by Regulation (EU) No. 253/2011.

Regulation (EC) No. 1272/2008, last amended by Regulation (EC) No. 790/2009

Legend

ADR	Accord européen relatif au transport international des marchandises dangereuses par route (European Agreement on the International Transport of Dangerous Goods by Road)
BImSchV	Bundes-Immissionsschutzverordnungen (German Federal Emission Protection Ordinance)
CAS	Chemical Abstracts Service
ChemV	Chemikalien Verordnung der Schweiz (Swiss Chemicals Regulation)
DIN	Norm des Deutschen Instituts für Normung (German Institute for Standardization)
EC	Effective Concentration
EC	European Community
EN	Europäische Norm (European Standard)
GHS	Globally Harmonized System
IATA-DGR	International Air Transport Association-Dangerous Goods Regulations
IBC Code	International Code for Building and Equipment of ships for the Carriage in bulk of dangerous chemicals
ICAO-TI	International Civil Aviation Organization - Technical Instructions
IMDG-Code	International Maritime Code for Dangerous Goods
ISO	International Organization for Standardization
IUCLID	International Uniform Chemical Information Database
LC	Lethal Concentration
LD	Lethal Dose
log Kow	Partition coefficient between octanol and water
MARPOL	Maritime Pollution Convention = Convention for the Prevention of Pollution from Ships
OECD	Organisation for Economic Cooperation and Development
PBT	Persistent, bioaccumulative, toxic
RID	Regulation on the International Carriage of Dangerous Goods by Rail
TRGS	Technische Regeln für Gefahrstoffe (Technical Rules for Hazardous Substances)
UN	United Nations
VOC	Volatile Organic Compounds
vPvB	very persistent and very bioaccumulative
VwVwS	German Regulation on Substances Hazardous to Water
WHC	Water Hazard Class

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60109948 SE_512018 AB 05/2014 © Synthes, Inc. or its affiliates
Subject to modifications Synthes is a trademark of Synthes, Inc. or its