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05.001.098

**Synthes Maintenance Spray, 400 ml,
for Electric Pen Drive and Air Pen Drive**

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**Spray de manutenção Synthes, 400 ml,
para Electric Pen Drive e Air Pen Drive**

Instruments and implants
approved by the AO Foundation

 **DePuy Synthes**

POWER TOOLS

COMPANIES OF *Johnson & Johnson*

05.001.098

Synthes Maintenance Spray, 400 ml

Intended Use

The Synthes Maintenance Spray is intended to lubricate surgical instruments, its components and for their maintenance.

General Information

Technically, the Synthes Maintenance Spray may support the mechanical properties of surgical instruments and its components, may help to reduce friction and corrosion and therefore heat generation and may help to improve the longevity of the instruments. The oil film is not impairing effective sterilization. Synthes Maintenance Spray is a silicone free product. Do not use Synthes Maintenance Spray for Anspach Power Tools.

Safety Instructions

The Synthes Maintenance Spray is only to be used for lubricating Synthes Power Tool Systems (Electric Pen Drive and Air Pen Drive including its components) and instruments after careful consultation of the instructions for use of the corresponding surgical instrument or power tool.

Synthes Maintenance Spray is just to be applied prior to sterilization.


The Synthes Maintenance Spray is designed for use by the reprocessing department or the Synthes service center and production personnel.

Do not use Synthes Maintenance Spray after the indicated expiry date.

To learn more about hazards, first aid measures, handling and storage, read the attached Material & Safety Data Sheet (MSDS).

Application

1. Check expiry date and remove lid
2. Attach the specific adapter for the instrument to lubricate:
 - Electric Pen Drive (EPD with adapter 05.001.101)
 - Air Pen Drive (APD, with adapter 05.001.092)
 - Attachments for EPD and APD (with adapter 05.001.102)
 - Screwdriver 90° (with adapter 05.505.002)
3. Consult the instruction for use of the corresponding power tool system or surgical instrument for proper use.

Precautions:

- Product contains extremely flammable gas under pressure; may explode if heated
- Store in a cool, dry and dark place
- Keep away from direct sunlight
- Read, follow and keep the instruction for use
- Wet cleaning under mechanical load may dissolve the etiquette of the Synthes Maintenance Spray

Explanation of Symbols

Caution

Consult the Instructions for Use before operating the device

Product is delivered in an unsterile condition

Keep away from direct sunlight

Indicates the expiry date

Material Safety Data Sheet (MSDS), SE_492023 AB

Synthes Maintenance Spray, 400 ml

1. Identification of the substance/mixture and of the company/undertaking**1.1 Product identifier**

Substance name / Trade name:

Synthes Maintenance Spray**1.2 Relevant identified uses of the substance or mixture and uses advised against**

Indication: Usage as lubricant

Contraindication: None

1.3 Details of the supplier of the safety data sheet**Manufacturer/Supplier**

Synthes GmbH

Address / PO Box

Eimattstr. 3

National code / postal code (ZIP) / City

CH-4436 Oberdorf

Contact for technical information

msds@synthes.com

Telephone/Telefax

+41 61 965 61 11

+41 61 965 66 00

1.4 Emergency telephone number

+41 44 251 51 51 TOX-Center Zurich

2. Hazards identification**2.1 Classification of the substance or mixture according to Regulation (EC) No. 1272/2008, Annex VII**

Flame, R phrases: R 12, Extremely flammable, R 18, In use, may form flammable/explosive vapour-air mixture S phrases: S 3/9, Keep container in a cool, well-ventilated place, S 16, Keep away from sources of ignition, No smoking, S 33, Take precautionary measures against static discharges, S 51, Use only in well-ventilated areas

05.001.098

Synthes Maintenance Spray, 400 ml

Classification according to Directive 67/548/EEC or Directive 1999/45/EC

Flammable Aerosol (H222)

Gases under pressure, Category compressed gas, H280

2.2 Label Elements

Labelling according to Regulation (EC) No 1272/2008 (substances) / Directive 1999/45/EC (mixtures)

Pictogram

Danger: Contains gas under pressure; may explode if heated

Precautionary statements (P-phrases)

Keep out of reach of children. (P102)

Do not spray on an open flame or other ignition source. (P211)

Pressurized container: Do not pierce or burn, even after use. (P251)

Protect from sunlight. Do not expose to temperatures exceeding 50 °C / 122 °F. (P410+P412)

3. Composition/Information on ingredients

3.1 Mixtures

Preparation of glycerin-fatty acid-ester and compressed gas (propane/butane)

Substance name: **Butane**

EC-No.: 203-448-7 / CAS-No.: 106-97-8

Contingent: 10–25 %

Substance name: **Propane**

EC-No.: 200-827-9 / CAS-No.: 74-98-6

Contingent: 10–25 %

Substance name: **Isobutane**

EC-No.: 200-857-2 / CAS-No.: 75-28-5

Contingent: 2–10 %

Classification according to regulation (EC) Nr. 1272/2008:
Flammable gases, Category 1; H220
Gases under pressure H280

(Full text of quoted hazard notes: see section 16)

4. First aid measures

4.1 Description of first aid measures

Following inhalation

Move victim to fresh air. In case of breathing arrest: artificial respiration.

Following skin contact

Remove / Take off all contaminated clothing. Rinse skin with soap and water, then rinse with plenty of water.

Following eye contact

Rinse with plenty of water. If eye irritation persists, consult a physician and take along this safety data sheet.

4.2 Most important symptoms and effects, both acute and delayed

N-Butane, Propane and Isobutane have no irritating effects as gas; if released from aerosol cylinders, cold injuries may occur. Systemic effects are to be expected with very high gas concentrations.

5. Firefighting measures

5.1 Extinguishing media

Suitable extinguishing media:

Dry powder extinguisher with ABC- or BC-powder, carbon dioxide extinguisher, foam extinguisher, water spray

5.2 Advice for firefighters

In case of fire eliminate cylinder from the danger zone or cool with firefighting water. Danger of explosion! Wear respiratory protection when fighting the fire.

6. Accidental release measures

6.1 Personal precautions, protective equipment and emergency procedures

Keep away from open flames, hot surfaces and sources of ignition

7. Handling and storage

7.1 Precautions for safe handling

Store dry at room temperature. Storage area should be well ventilated.

Local exhaust ventilation at the workplace is necessary.

7.2 Conditions for safe storage, including any incompatibilities

Information on storage conditions

Do not store in escape routes, work rooms and in their close vicinity.

The substance should not be stored with substances which are liable to cause dangerous chemical reactions (e.g. strong acids and alkalis, oxidizing substances).

Storage class: 2 Aerosols

8. Exposure controls / personal protection

8.1 Parameters to be monitored

8.1.1 Limits for exposure at the workplace and/or biological limits

Maximum allowable concentration (MAC) Germany, SUVA 1903 CH

Substance name: **Butane; CAS-No.: 106-97-8:**

1000 ml/m³, 2400 mg/m³, peak limit: excess factor 4, duration 15 min, mean value; 4 times per shift; interval 1 h category II – Substances with a resorptive effect

Substance name: **Propane; CAS-No.: 74-98-6:**

1000 ml/m³, 1800 mg/m³, peak limit: excess factor 4, duration 15 min, mean value; 4 times per shift; interval 1 h category II – Substances with a resorptive effect

Substance name: **Isobutane; CAS-No.: 75-28-5:**

1000 ml/m³, 2400 mg/m³, peak limit: excess factor 4, duration 15 min, mean value; 4 times per shift; interval 1 h category II – Substances with a resorptive effect

8.2 Limiting and monitoring exposure

8.2.1 Individual protecting measures / personal protective equipment

Eye/Face protection

Wear eye protection if handling the product.

Skin protection

Not necessary if used appropriately.

In case of long, excessive or repeated skin contact wear protective gloves (e.g. made of butyl rubber suitable for applications up to 8 h).

Work hygiene

Avoid skin contact with the liquid phase, frostbites.

Avoid inhaling the gas.

Remove gas contaminated clothing and air them far from ignition sources.

9. Physical and chemical properties

9.1 Details on the basic physical and chemical properties

Appearance:

– Aggregate state: Aerosol

– Colour: Colourless

Odour: Odourless

pH-value: Not applicable

Initial boiling point / boiling range:

Compressed gas $-44\text{ }^{\circ}\text{C}$ / $-47.2\text{ }^{\circ}\text{F}$

Flashing point:

Agent: $>210\text{ }^{\circ}\text{C}$ / $410\text{ }^{\circ}\text{F}$, DIN EN ISO 2592

Compressed gas: $-97\text{ }^{\circ}\text{C}$ / $-142.6\text{ }^{\circ}\text{F}$

Flammability (solid, gaseous): flammable

Upper/lower explosion threshold: 12.0 Vol. % / 1.5 Vol. %

Vapour pressure:

Agent: $<0.1\text{ mbar}$ at $20\text{ }^{\circ}\text{C}$ / $68\text{ }^{\circ}\text{F}$, DIN 51519

Compressed gas: 3.5 bar at $20\text{ }^{\circ}\text{C}$ / $68\text{ }^{\circ}\text{F}$

Solubility(ies): Insoluble

Spontaneous ignition temperature:

Agent: $>300\text{ }^{\circ}\text{C}$ / $572\text{ }^{\circ}\text{F}$, DIN 51794

Compressed gas: $-97\text{ }^{\circ}\text{C}$ / $-142.6\text{ }^{\circ}\text{F}$

Viscosity: Agent: $15\text{ mm}^2/\text{s}$ at $40\text{ }^{\circ}\text{C}$ / $104\text{ }^{\circ}\text{F}$, DIN EN ISO 3104

Explosive properties: Gas-air mixture is potentially explosive

10. Stability and reactivity

10.1 Thermal decomposition: >220 °C / 428 °F

10.2 Hazardous decomposition products

Unknown if stored and handled appropriately. In case of improper use (e.g. high temperature) or fire, harmful vapours may be released to the environment.

10.3 Hazardous reactions

With strong oxidants (e.g. liquid or condensed oxygen), with very strong acids and bases, high temperatures, naked flames and other ignition sources.

11. Toxicological information

11.1 Details on toxicological effects

Details on the base oil

Oral toxicity LD₅₀ >20 000 mg/kg (mouse)

Irritant effect on the skin

Not sensitizing

The lubricant contained in the aerosol cans is registered at NSF International (NSF/H1-release).

12. Ecological information

The agent contained in the aerosol cans is not classified as a substance hazardous to water, without WHC, ID-Number 760. It is bio-degradable.

13. Disposal considerations

Completely emptied cans may be disposed of with regular household garbage. Not completely emptied cans underlie country-specific regulations concerning the disposal of special wastes. Special waste code: 16 05 04

14. Transport information

14.1 Land transport (ADR / RID / GGVS / GGVE)

Class(es): 2

UN-No.: 1950

Packing group: –

Hazard label(s): 2.1

Proper shipping name: Aerosol dispenser, flammable

Limited Quantity (LQ): 1 liter

Transport class: 2

Tunnel restriction code: D

14.2 Sea transport (IMDG-Code / GGVSee)

Class(es): 2

UN-No.: 1950

Hazard label(s): 2.1

Proper shipping name: Aerosols, flammable

Packing group: –

Marine Pollutant: No

EMS-Nr.: F-D, S-U

14.3 Air transport (ICAO / IATA-DGR)

Class(es): 2

UN / ID-No.: UN1950

Label: 2.1

Proper shipping name: Aerosols, flammable

Packing group: –

14.4 Environmental perils

Label for substances which are harmful to the environment

ADR / RID / IMDG-Code / ICAO-TI / IATA-DGR: no

Marine Pollutant: no

15. Regulatory information

15.1 Safety, health and environmental regulations / legislations specific for the substance or mixture

The safety regulations on the aerosol dispenser are valid.

Labeling of H- and P-statements see chapter 2.

The product contains up to 20 % volatile organic compounds according to the Swiss positive list Annex I and VOC according to TRGS 220, 31 respectively. BimSchV (approx. 50 g per can).

The nominal pressure (test pressure) of the aerosol dispenser is 15 bar.

The optimum pressure is approx. 18 bar.

16. Other information

Legend

ADR	European Agreement on the International Transport of Dangerous Goods by Road
BImSchV	German Federal Emission Protection Ordinance
CAS	Chemical Abstracts Service
ChemV	Swiss Chemicals Regulation
DIN	German Institute for Standardization
EC	Effective Concentration
EC	European Community
EN	European Standard
GHS	Globally Harmonized System
IATA-DGR	International Air Transport Association – Dangerous Goods Regulations
IBC-Code	International Code for Building and Equipment of ships for the Carriage in bulk of dangerous chemicals
ICAO-TI	International Civil Aviation Organization – Technical Instructions
IMDG-Code	International Maritime Code for Dangerous Goods
ISO	International Organization for Standardization
IUCLID	International Uniform Chemical Information Database
LC	Lethal Concentration
LD	Lethal Dose
log Kow	Distribution Coefficient between octanol and water
MARPOL	Maritime Pollution Convention
OECD	Organisation for Economic Co-operation and Development
PBT	Persistent, bioaccumulative, toxic
RID	Regulation on the International Carriage of Dangerous Goods by Rail
TRGS	Technical Rules for Hazardous Substances
UN	United Nations
VOC	Volatile Organic Compounds
vPvB	very persistent and very bioaccumulative
VwVwS	German Regulation on Substances Hazardous to Water
WHC	Water Hazard Class

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60099550 SE_484989 AA 04/2013 40100163 © Synthes, Inc. or its affiliates
Subject to modifications Synthes is a trademark of Synthes, Inc. or its affiliates