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519.970
Olio speciale Synthes, 40 ml

519.970
Óleo especial Synthes, 40 ml

Instruments and implants
approved by the AO Foundation

 **DePuy Synthes**
POWER TOOLS

COMPANIES OF *Johnson & Johnson*

519.970

Synthes Special Oil, 40 ml

Intended Use

The Synthes Special Oil is intended to lubricate surgical instruments, its components and for their maintenance.

General Information

Technically, the Synthes Special Oil may support the mechanical properties of surgical instruments and its components, may help to reduce friction and corrosion and therefore heat generation and may help to improve the longevity of the instruments. The oil film is not impairing effective sterilization. Synthes Special Oil is a silicone free product.

Safety Instructions

The Synthes Special Oil is only to be used for lubricating Synthes Power Tool Systems and instruments after careful consultation of the instructions for use of the corresponding surgical instrument or power tool. Do not use this product for Electric Pen Drive, Air Pen Drive and Anspach Power Tools.

Synthes Special Oil is to be applied just prior to sterilization.

The Synthes Special Oil is designed for use by the reprocessing department or the Synthes service center and production personnel.

Do not use Synthes Special Oil after the indicated expiry date.


To learn more about hazards, first aid measures, handling and storage, read the attached Material & Safety Data Sheet (MSDS).

Application

1. Check expiry date and remove lid.
2. Remove the safety seal of the dispenser with scissors.
3. For best handling and precise application we recommend the use of a dispense tip (Nordson, green, $d=0.84$ mm).
4. Consult the instruction for use of the corresponding power tool system or surgical instrument for proper use.

Precautions:

- Store in a cool, dry and dark place
- Keep away from direct sunlight
- Read, follow and keep the instruction for use

Explanation of Symbols

Caution

Consult the Instructions for Use before operating the device

Product is delivered in an unsterile condition

Keep away from direct sunlight

Indicates the expiry date

519.970

Synthes Special Oil, 40 ml

Material Safety Data Sheet (MSDS), SE_492010 AB

Synthes Special Oil

1. Identification of the substance/mixture and of the company/undertaking

1.1 Product identifier

Substance name / Trade name:

Synthes Special Oil

1.2 Relevant identified uses of the substance or mixture and uses advised against

Indication: Usage as lubricant for medical devices

Contraindication: None

1.3 Details of the supplier of the safety data sheet Manufacturer/Supplier

Synthes GmbH

Address / PO Box

Eimattstr. 3

National code / postal code (ZIP) / City

CH-4436 Oberdorf

Contact for technical information

msds@synthes.com

Telephone / Telefax

+41 61 965 61 11

+41 61 965 66 00

1.4 Emergency telephone number:

+41 44 251 51 51 TOX-Center Zurich

2. Hazards identification

2.1 Classification of the substance or mixture according to Regulation (EC) No. 1272/2008, Annex VII

Not classified as hazardous according to the European regulation (EC) 1272/2008.

Classification according to Directive 67/548/EEC or Directive 1999/45/EC

Substance not classified according to Directive 67/548/EEC.

2.2 Label Elements

Labelling according to Regulation (EC) No 1272/2008 (substances) / Directive 1999/45/EC (mixtures)

none

2.3 Other hazards

The product is biologically inert.

Not hazardous in normal conditions of handling and use.

Ecological injuries are not known or expected under normal use.

Thermal decomposition at $>290^{\circ}\text{C}$ can lead to release of toxic and corrosive gases.

3. Composition/Information of ingredients

3.1 Substances

Main ingredient of the substance:

Substance name:

1-Propene, 1,1,2,3,3,3-hexafluoro-, oxidized, polymd.

EC-No.: – / CAS-No.: 69991-67-9

Contingent: –

(Full text of quoted hazard notes: see section 16)

4. First aid measures

4.1 Description of first aid measures

Remove contaminated clothing immediately.

If health problems occur, consult a physician.

Following inhalation

Move to fresh air in case of accidental inhalation of fumes from overheating or combustion.

Use oxygen or artificial respiration if needed.

Following skin contact

Wash off with soap and water.

If symptoms persist, call a physician.

Following eye contact

Rinse immediately with plenty of water, also under the eyelids, for at least 15 minutes.

If eye irritation persists, consult a specialist.

Following ingestion

Drink 1 or 2 glasses of water.

Do NOT induce vomiting.

If symptoms persist, call a physician.

4.2 Most important symptoms and effects, both acute and delayed**Inhalation**

No known effect

Skin contact

Symptoms: Redness

Eye contact

Contact with eyes may cause irritation, redness

Ingestion

Ingestion may provoke the following symptoms:

Nausea, Vomiting, Diarrhoea

5. Firefighting measures**5.1 Extinguishing media**

Suitable extinguishing media:

Water, powder, foam, dry chemical, carbon dioxide (CO₂)

Unsuitable extinguishing media:

None

5.2 Special hazards arising from the chemical or mixture

The product is not flammable and not explosive.

In case of fire at temperatures above 290 °C, hazardous decomposition products may be produced such as:

Gaseous hydrogen fluoride (HF), fluorophosgene.

5.3 Firefighting measures

Wear self-contained breathing apparatus and protective suit.

When intervention in close proximity wear acid-resistant oversuit.

Evacuate personnel to safe areas.

Approach from upwind.

Protect intervention team with a water spray as they approach the fire.

Keep containers and surroundings cool with water spray.

Keep product and empty container away from heat and sources of ignition.

6. Accidental release measures

6.1 Personal precautions, protective equipment and emergency procedures

Sweep up to prevent slipping hazard.

Prevent further leakage or spillage if safe to do so.

Keep away from open flames, hot surfaces and sources of ignition.

Refer to protective measures listed in sections 7 and 8.

6.2 Environmental precautions

Avoid the uncontrolled release of the substance, particularly into surface, ground or waste water and into soils.

6.3 Methods and material for containment and cleaning up

Pick up with suitable absorbent material, e.g. earth, sand or saw-dust, and shovel into suitable containers for disposal.

7. Handling and storage

7.1 Precautions for safe handling

– Store dry at room temperature.

– Ensure adequate ventilation.

– Keep away from sources of heat and ignition.

– To avoid thermal decomposition, do not overheat.

– Clean and dry equipment before starting work.

7.2 Conditions for storage, including incompatibilities

Details of storage conditions

– Store dry at room temperature.

– Keep away from heat and sources of ignition.

– Keep in properly labeled containers.

– Keep away from combustible materials.

– Keep away from incompatible products.

– Refer to protective measures listed in sections 6 and 8.

Requirements for storage areas and containers

Store dry at room temperature. Store in tight-sealing containers made of e.g. plastic, glass or metal.

8. Exposure controls / personal protection

8.1 Parameter to be monitored

Remark: Threshold limit value of by-products from thermal decomposition decomposition $>290^{\circ}\text{C}$.

8.1.1 Exposure Limit Values in the workplace and/or biological limit values

Workplace Limit Values (WLV) in Germany

Substance name: **Hydrogen fluoride anhydrous**
Specification: EU. Indicative Exposure and Directives relating to the protection of risks related to work exposure to chemical, physical, and biological agents / 12 2009.

Value: Time weighted average: 1.8 ppm / 1.5 mg/m³; Remarks: Indicative

Peak limit: Short term exposure limit: 3 ppm / 2.5 mg/m³

Teratogenic: –

Survey method: –

Substance name: **Hydrogen fluoride anhydrous**
Specification: Germany. DFG MAK List (advisory OELs). Commission for the Investigation of Health Hazards of Chemical Compounds in the Work Area (DFG 2008)
Remarks: Category I: Substances with a Workplace Exposure Limit for their local effect or for substances with a sensitizing effect in the airways.

Value: WEL: 1 ppm / 0.83 mg/m³

Peak limit: –

Teratogenic: –

Survey method: –

Substance name: **Hydrogen fluoride anhydrous**
Specification: Germany. TRGS 900, Threshold Limit Values in air in the workplace 07 / 2009.
Remarks: as F, Alveolar dust fraction.
Remarks: Can be absorbed through skin.

Value: WLV 1 ppm
WLV 0.83 mg/m³

Peak limit: –

Teratogenic: –

Survey method: –

Substance name: **Hydrogen fluoride anhydrous**
Specification: US. ACGIH Threshold Limit Values 2009.
Remarks: As F

Value: 0.5 ppm

Peak limit: 2 ppm

Teratogenic: –

Survey method: –

Substance name: **Carbonyl difluoride**
Specification: EU. Indicative Exposure and Directives relating to the protection of risks related to work exposure to chemical, physical, and biological agents / 12 2009.

Value: Time-weighted average: 2.5 mg/m³,
Remarks: indicative.

Peak limit: –

Teratogenic: –

Survey method: –

Substance name: **Carbonyl difluoride**
Specification: Germany. DFG MAK List (advisory OELs).
Commission for the Investigation of Health Hazards of Chemical Compounds in the Work Area (DFG 2008)

Remarks: Entered as F, Alveolar dust fraction. Remarks: As F, Alveolar dust fraction, can be absorbed through skin. Remarks: As F, Alveolar dust fraction, Category II: substances with a resorptive effect. Can be absorbed through skin.

Value: WEL 1 mg/m³

Peak limit: –

Teratogenic: –

Survey method: –

Substance name: **Carbonyl difluoride**
 Specification: Germany. TRGS 900, Threshold Limit Values in air in the workplace 07 / 2009.
 Remarks: as F, Alveolar dust fraction.
 Can be absorbed through skin

Value: WLV: 1 mg/m³

Peak limit: –

Teratogenic: –

Survey method: –

Substance name: **Carbonyl difluoride**
 Specification: US. ACGIH Threshold Limit Values 2009
 Value: Time-weighted average: 2 ppm
 Peak limit: Short-term exposure limit: 5 ppm
 Teratogenic: –
 Survey method: –

8.2.2 Environmental exposure controls

Dispose of rinse water in accordance with local and national regulations.

9. Physical and chemical properties

9.1 Details on the basic physical and chemical properties

Appearance

– Aggregate state: Viscous liquid

– Colour: Colourless

Odour: Odourless

Odour threshold: –

pH-value: –

Melting point/
freeze point: Not available

Initial boiling point/
boiling range: No data available

Flashing point: The product is not flammable

Evaporation rate: –

Flammability (solid,
gaseous): The product is not flammable.

Upper/lower flammability or explosion
threshold: Not explosive

Vapor pressure:	0.0000001–0.01 hPa / Temperature: 20 °C
Vapor density:	No data available
Relative density:	1.85–1.93 g/ml
Solubility(ies):	In water: insoluble In fluorinated solvents: soluble
Partition coefficient n-octanol/water:	–
Spontaneous ignition temperature:	Not flammable, not oxidising
Decomposition temperature:	>290 °C
Viscosity:	40 mm ² /s; Temperature 40 °C
Explosive properties:	Not explosive
Oxidizing properties:	–

10. Stability and reactivity

10.1 Reactivity

No dangerous reaction known under conditions of normal use.

10.2 Chemical stability

Powdered metals (Al, Mg) promote decomposition and lower the decomposition temperature.

10.3 Conditions to avoid

To avoid thermal decomposition, do not overheat.

Keep away from flames and sparks.

The product can decompose at temperatures above 290 °C

10.4 Incompatible materials

Flammable materials, combustible material, metals promote decomposition and lower decomposition temperature, Lewis acids (Friedel-Crafts) above 100 °C, Aluminium and magnesium in powder form above 200 °C.

10.5 Hazardous decomposition products

Gaseous hydrogen fluoride (HF), fluorophosgene at temperatures >290 °C

11. Toxicological information

11.1 Details on toxicological effects

Acute oral toxicity

LD50, rat, >15000 mg/kg

Ingestion may provoke the following symptoms:

Nausea, Vomiting, Diarrhoea

Acute dermal toxicity

LD50, rat, >5000 mg/kg

Skin corrosion / irritation

Rabbit, no skin irritation

Rabbit, no skin irritation, remarks: 14 days

Symptoms: Redness

Eye damage / irritation

Rabbit, no eye irritation

Contact with eyes may cause irritation, redness

Sensitization of respiratory system / skin

Guinea pig, did not cause sensitization on laboratory animals.

Remarks: Derma

Germ cell-mutagenicity

Not mutagenic in Ames Test

Chromosome aberration test in vitro, negative

Risk of aspiration

No known effect.

Other information

Description of possible hazardous to health effects is based on experience and / or toxicological characteristics of individual components.

The product is biologically inert.

Thermal decomposition can lead to release of toxic and corrosive gases.

Exposure to decomposition products causes severe irritation of eyes, skin and mucous membranes.

12. Ecological information**12.1 Toxicity**

Fishes, *Brachydanio rerio*, LC50, 96 h, >360 mg/l,
remarks: saturated aqueous solution

Crustaceans, *Daphnia magna*, EC50, >360 mg/l
remarks: saturated aqueous solution

12.2 Persistence and degradability

No data available

12.3 Bioaccumulative potential

No data available

12.4 Mobility in soil

No data available

12.5 Result of PBT and vPvB assessment

No data available

12.6 Other adverse effects

Ecological injuries are not known or expected during correct use.

13. Disposal considerations

13.1 Waste treatment methods

Management of contaminated packaging

Empty containers can be landfilled, when in accordance with the local regulations.

Special precautions

Can be incinerated, when in compliance with local regulations.

The incinerator must be equipped with a system for the neutralisation or recovery of HF.

In accordance with local and national regulations.

14. Transport information

The product is not declared as dangerous good.

15. Regulatory information

15.1 Safety, health and environmental regulations/legislations specific for the substance or mixture

EU regulations

Not classified according to Directive 67/548/EEC and ChemG SR 813.1 (ChemV SR 813.11).

National regulations

Observe national guidelines.

German Administrative Regulation on Water-endangering Substances – VwVwS

WGK=1 (slightly water endangering) – WGK Id number: 6526

Further relevant regulations

European Waste Catalogue, hazardous waste, waste codes should

519.970

Synthes Special Oil, 40 ml

be assigned by the user based on the application for which the product was used.

16. Other information

References and data origins

Regulations

Preparation regulation (1999/45/EC), last changed by regulation (EC) No. 1907/2006

Substance guideline (67/548/EWG), last changed by guideline 2009/2/EC

REACH regulation (EC) No. 1907/2006, last changed by regulation (EU) No. 253/2011

Regulation (EC) No. 1272/2008, last changed by regulation (EC) No. 790/2009

Swiss Chemicals Law: ChemG SR 813.1 and ChemV SR 813.11

Legend

ADR	European Agreement on the International Transport of Dangerous Goods by Road
BImSchV	German Federal Emission Protection Ordinance
CAS	Chemical Abstracts Service
ChemV	Swiss Chemicals Regulation
DIN	German Institute for Standardization
EC	Effective Concentration
EC	European Community
EN	European Standard
GHS	Globally Harmonized System
IATA-DGR	International Air Transport Association – Dangerous Goods Regulations
IBC-Code	International Code for Building and Equipment of ships for the Carriage in bulk of dangerous chemicals
ICAO-TI	International Civil Aviation Organization – Technical Instructions
IMDG-Code	International Maritime Code for Dangerous Goods
ISO	International Organization for Standardization
IUCLID	International Uniform Chemical Information Database
LC	Lethal Concentration
LD	Lethal Dose
log Kow	Distribution Coefficient between octanol and water
MARPOL	Maritime Pollution Convention
OECD	Organisation for Economic Co-operation and Development
PBT	Persistent, bioaccumulative, toxic
RID	Regulation on the International Carriage of Dangerous Goods by Rail

TRGS	Technical Rules for Hazardous Substances
UN	U nited N ations
VOC	V olatile O rganic C ompounds
vPvB	v ery p ersistent and v ery b ioaccumulative
VwVwS	German Regulation on Substances Hazardous to Water
WHC	W ater H azard C lass

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60099544 SE_484729 AA 04/2013 40100174 © Synthes, Inc. or its affiliates
Subject to modifications Synthes is a trademark of Synthes, Inc. or its affiliates