Instructions for Use

04.630.130S FACET WEDGE, Small

04.630.131S FACET WEDGE, Medium

04.630.132S FACET WEDGE, Large

04.630.135.02S FACET WEDGE Screw, \varnothing 3.0 mm,
L 12 mm

PART OF THE JOHNSON-JOHNSON FAMILY OF COMPANIES

Instructions for Use

04.630.130S FACET WEDGE, Small

04.630.131S FACET WEDGE, Medium

04.630.132S FACET WEDGE, Large

**04.630.135.02S FACET WEDGE Screw, Ø 3.0 mm,
L 12 mm**

Introduction

Prior to use:

Please read these instructions for use, the "Important Information" brochure and the corresponding FACET WEDGE surgical technique carefully before use. Ensure that you are familiar with the appropriate surgical technique.

Material

TAN (Ti-6Al-7Nb) ISO 5832-11

Intended use

FACET WEDGE is intended for the fixation of the spine as an aid to fusion through immobilization of the facet joints, with or without bone graft, at one or two levels, from L1 to S1. FACET WEDGE can be inserted minimal invasively and should

only be used to augment other fusion and stabilization techniques.

Indications

- Degenerative disc disease
- Degenerative facet joint disease (isolated facet based symptomatic back pain)
- Pseudarthrosis post anterior fusion with intact instrumentation

Contraindications

- Unilateral application, except in combination with pedicle screw fixation on the contralateral side
- Compromised facets due to decompression techniques
- Spondylolisthesis
- Fracture or other instabilities of the posterior elements
- Tumor
- Acute or chronic systemic or localized spinal infections

Potential adverse events

- Nonunion (Pseudarthrosis) or delayed union
- Bending, fracture, loosening, or migration of the implant
- Decrease in bone density due to stress shielding
- Sensitivity or foreign body reaction

English

- Postoperative pain or discomfort
- Nerve, soft tissue, or blood vessel damage
- Fracture of bony structures
- Nerve root or spinal cord impingement
- Dural leak
- Infection
- Necrosis of bone
- Hemorrhage

Sterile device

STERILE R Sterilized using irradiation

Store implants in their original protective packaging, and do not remove them from the packaging until immediately before use.

Prior to use, check the product expiration date and verify the integrity of the sterile packaging. Do not use, if the package is damaged.

Do not resterilize

Single use device

Do not reuse

Products intended for single-use must not be re-used.

Reuse or reprocessing (e.g. cleaning and resterilization) may compromise the structural integrity of the device and/or lead to device failure which may result in patient injury, illness or death.

Furthermore, reuse or reprocessing of single use devices may create a risk of contamination e.g. due to the transmission of infectious material from one patient to another. This could result in injury or death of the patient or user.

Contaminated implants must not be reprocessed. Any Synthes implant that has been contaminated by blood, tissue, and/or bodily fluids/matter should never be used again and should be handled according to hospital protocol. Even though they may appear undamaged, the implants may have small defects and internal stress patterns that may cause material fatigue.

Warnings

It is strongly advised that FACET WEDGE is implanted only by

English

operating surgeons who are familiar with the general problems of spinal surgery and who are able to master the product-specific surgical techniques. Implantation is to take place with the instructions for the recommended surgical procedure. The surgeon is responsible for ensuring that the operation is carried out properly.

For additional information, please refer to the FACET WEDGE surgical technique.

The manufacturer is not responsible for any complications arising from incorrect diagnosis, choice of incorrect implant, incorrectly combined implant components and/or operating techniques, the limitations of treatment methods, or inadequate asepsis.

Precautions

- Prepare the facet joints to ensure the correct placement of the FACET WEDGE implant
- The selection of a too large implant might lead to closing of the facet joint on the contralateral side
- The general risks associated with surgery are not described in these instructions for use. For more information, please refer to the enclosed Important Information

Combination of medical devices

Synthes has not tested compatibility with devices provided by other manufacturers and assumes no liability in such instances.

Magnetic Resonance environment

MRI safety

Non-clinical testing of the worst-case scenario has demonstrated that the implants of the FACET WEDGE system are MR conditional. These articles can be scanned safely under the following conditions:

- Static magnetic field of 1.5 Tesla and 3.0 Tesla.
- Spatial gradient field of 300 mT/cm (3000 Gauss/cm).
- Maximum whole body averaged specific absorption rate (SAR) of 2 W/kg for 15 minutes of scanning.

Based on non-clinical testing, the FACET WEDGE implant will produce a temperature rise not greater than 3.5 °C at a maximum whole body averaged specific absorption rate (SAR) of 2 W/kg, as assessed by calorimetry for 15 minutes of MR scanning in a 1.5 Tesla and 3.0 Tesla MR scanner.

MR Imaging quality may be compromised if the area of interest is in the exact same area or relatively close to the position of the FACET WEDGE device.

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© Synthes GmbH 2017. All rights reserved. 60095640 SE_463713 AF 02/2017
Synthes GmbH