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---

# **chronOS Strip** **Bone Void Filler**

Instruments and Implants  
approved by the AO Foundation

---

**DePuy Synthes**

PART OF THE *JOHNSON & JOHNSON* FAMILY OF COMPANIES

Version SE 460279 AD  
Date 2018 – 06

<b>English</b>	<b>3</b>
<b>Deutsch</b>	<b>9</b>
<b>Espanol</b>	<b>15</b>
<b>Français</b>	<b>21</b>
<b>Italiano</b>	<b>27</b>
<b>Português</b>	<b>33</b>

This instruction for use is not intended for distribution in the USA.

### chronOS Strip Bone Void Filler

Please read these instructions for use, the Synthes "Important Information" and the corresponding surgical technique carefully before use. Ensure that you are familiar with the appropriate surgical technique. Sufficient experience in the field of the use of biomaterials is highly recommended before using chronOS Strip Bone Void Filler Implants.

#### Material(s)

$\beta$ -tricalcium phosphate [Ca<sub>3</sub>(PO<sub>4</sub>)<sub>2</sub>]  
Poly(L-lactide-co- $\epsilon$ -caprolactone)

#### Standard(s):

ASTM-F 1088  
N/A

#### Description of the Device

chronOS Strip Bone Void Filler is a synthetic, porous, resorbable and biocompatible bone void filler made from chronOS Granules with a defined, uniform pore structure imbedded in a matrix of poly(L-lactide-co- $\epsilon$ -caprolactone). The flexible composite provides contourability and malleability of the device to the bony implant site. chronOS Strip Bone Void Filler is initially radio-opaque. chronOS Strip Bone Void Filler is resorbed in vital bone and replaced by autologous bone. chronOS Strip Bone Void Filler is provided in a sterile perfusion pouch. chronOS Strip Bone Void Filler can be perfused with the patient's bone marrow and/or blood.

#### Intended Use

chronOS Strip Bone Void Filler is a bone replacement material of  $\beta$ -tricalcium phosphate imbedded in a matrix of poly(L-lactide-co- $\epsilon$ -caprolactone). In its function as temporary bone substitute, it serves to fill and bridge aseptic bone defects.

#### Indications

chronOS Strip Bone Void Filler should be used as bone void fillers or augmentation material in zones requiring cancellous rather than cortical bone. This includes the filling of bone defects after trauma, reconstruction or correction in non-load bearing spinal indications only.

The porous structure of chronOS Strip Bone Void Filler acts as a matrix for the ingrowth of bone. chronOS Strip Bone Void Filler implants must always be applied by endosteal or subperiosteal implantation, i.e. by direct contact with the vital bone.

The implant must completely fill the bone defect (press-fit). However, it is essential to avoid overfilling to ensure tension-free wound closure.

chronOS Strip Bone Void Filler may be combined with autogenous bone marrow and/or blood. The size of the chronOS Strip Bone Void Filler implant used depends on the size of the bone defect.

As a rule, chronOS Strip Bone Void Filler is resorbed within 6 to 18 months depending on patient conditions and converted into autologous bone.

#### Contraindications

chronOS Strip Bone Void Filler should not be used in the following circumstances:

- Acute and chronic infections in the operation area, e.g. inflammation, bacterial bone diseases (posttraumatic or chronic osteomyelitis) and soft-tissue infections
- Untreated malignant myeloma, Burkitt's lymphoma and other lymphomas
- Defects and fractures in the region of an open epiphysis
- Osteoporosis
- Areas where there is significant vascular impairment proximal to the graft site
- Severe instability or deformation at the extraction point (harvesting site)
- Cranioplasty

chronOS Strip Bone Void Filler is not indicated for use in load bearing and unstable indications unless used in conjunction with appropriate osteosynthesis fixation systems or if the cortical bone can bear the full load.

**Side Effects/Adverse Events**

- General complications caused by invasive surgery
  - Nonunion or delayed union, which may lead to a failure of the implant
  - Pain, discomfort, abnormal sensation, or palpability due to the presence of the device
  - Increased fibrous tissue response around fracture site and/or the implant
- Apart from these possible adverse effects there is also the risk of complications associated with any surgical procedure featuring bony defects such as, necrosis of bone, infection, nerve damage and pain which may not be related to the implant. In general, good tissue response of  $\beta$ -tricalcium phosphate implants in bone is supported by experimental and clinical data. Nevertheless, the following complications are possible:
- Device fragment displacement that could potentially result in ectopic bone growth
  - Neurovascular injuries caused by surgical trauma
  - Foreign body reactions
  - Allergic reactions
  - Inflammatory reactions
  - Infections that can lead to failure of the procedure

**Sterile Device**

Sterilized using gamma-irradiation

- chronOS Strip Bone Void Filler is supplied sterile. Store implants in their original protective packaging, and do not remove them from the packaging until immediately before use.
- Prior to use, check the product expiration date and verify the integrity of the sterile packaging.
- Do not use, if the package is damaged or if the sterile barrier has been compromised. When removing the implant from its packaging, strictly observe the instruction concerning asepsis.

Do not re-sterilize

- This product is intended for single use and must not be re-used. Re-use or reprocessing (e.g. cleaning and/or re-sterilization) may compromise the sterility or integrity of the device which may lead to infection, device failure and/or contamination.

**Single Use Device**

Do not re-use

- This product is intended for single use and must not be re-used. Re-use or reprocessing (e.g. cleaning and/or re-sterilization) may compromise the sterility or integrity of the device which may lead to infection, device failure and/or contamination.
- The contents of an opened or damaged outer protective packaging may no longer be used and must be destroyed.

- Once the packaging has been opened, the content must be used immediately. Any unused parts, fragments or modified parts must be discarded in an approved biohazard waste container.
- Do not use after the expiry date.

### **Precautions/Restrictions**

The use of chronOS Strip Bone Void Filler has a more restricted indication in:

- Severe, endocrine-induced bone diseases (e.g. hyperparathyroidism)
- Current therapy with steroids and with drugs, which intervene in calcium metabolism (e.g. calcitonin)
- Severe, poorly controlled diabetes (diabetes mellitus) with bad wound healing tendencies
- Immunosuppressive therapy
- Poor bone quality

The use of chronOS Strip Bone Void Filler is not known in:

- Cardiovascular diseases
- Pregnancy and lactation period
- Filling of bone defects in congenital deformations
- Immunosuppressive or radiation therapy
- Compromised immune system, compromised wound healing
- When using osteosynthesis fixation systems, the use of Synthes implants is recommended.

The effect of mixing chronOS Strip Bone Void Filler with medicinal products or other synthetic bone void fillers is not known.

As there are no known clinical data available for chronOS Strip Bone Void Filler in spinal fusion procedures spanning more than 2 levels, its use in these cases is not recommended. The treating physician should weigh the benefits of applying chronOS Strip Bone Void Filler in more than 2 levels against the potential risk to the patient.

### **Warnings**

- Do not use chronOS Strip Bone Void Filler as a stand-alone product unless the cortical bone can bear the full load. Maximum mechanical stability in filling the defect with chronOS Strip Bone Void Filler is a prerequisite for a good bony incorporation. If there is instability, use appropriate internal fixation to assure rigid stabilization in all planes.
- chronOS Strip Bone Void Filler is radio-opaque until replaced with new bone. Radio-opacity may mask underlying pathological conditions. Radio-opacity may also make it difficult to radiographically assess the ingrowth of new bone.
- The Perfusion Pouch is intended for the impregnation of chronOS Strip Bone Void Filler with bone marrow or blood.
- Examine the package for damages before using the sterile bone substitute, as they might impair sterility. This applies to the aluminium peel-off package. In removing the implant from its package, strictly observe the instructions concerning aseptic procedures.
- Do not attempt to re-sterilize the unused contents of an opened package, but dispose it. Re-sterilization of chronOS Strip Bone Void Filler can result in product not being sterile, and/or not meeting performance specifications and/or alters material properties.
- Do not use chronOS Strip Bone Void Filler after the expiry date printed on the packaging.
- Do not use chronOS Strip Bone Void Filler instead of ODL Spacers with the Arch Laminoplasty System.
- Do not fill into artificial compartments (e.g. filling of interbody fusion devices).

### **Combination of the Medical Device with a Medicinal Product**

The effect of mixing chronOS Strip Bone Void Filler Implants with medicinal products is not known.

## **Interactions**

No negative interactions with bone marrow and autologous blood as of yet have been reported.

## **Magnetic Resonance (MR) Environment**

chronOS Strip Bone Void Filler Implants are of non-metallic origin. This material is inherently diamagnetic and cannot be heated up or act as an antenna either by bringing the patients into the MRI magnet or during the MRI examinations. However, chronOS Strip Bone Void Filler Implants may be used in combination with metal hardware. Unless stated otherwise, devices have not been evaluated for safety and compatibility within the MR environment. Please note that there are potential hazards which include:

- Heating or migration of the device
- Artifacts on MR image

## **Handling / Special Operating Instructions**

Keep the sterile chronOS Strip Bone Void Filler in its protective packaging until ready to use. Check that the packaging has not been damaged before opening as this might impair sterility. In removing the implant from its packaging, strictly observe the instructions concerning asepsis.

chronOS Strip Bone Void Filler must always be applied by endosteal or subperiosteal implantation, i.e. by direct contact with the vital bone. The surgical procedure depends on the localization, nature and extent of the bone defect. To prepare the implant site, remove all inflamed necrotic tissue and bone particles.

The appropriate implant size is determined by the size of the bone defect and the vascularity. Keep the regions of endosteal vessel and nerve cords clear to avoid pressure sores.

chronOS Strip Bone Void Filler may be contoured, molded, cut, twisted, or otherwise manipulated by the surgeon intraoperatively to fit the defect as dictated by surgical preference.

In case of using bone graft harvesting set, we do not recommend to fill the voids with chronOS Strip Bone Void Filler.

## **Perfusion Pouch**

chronOS Strip Bone Void Filler is provided in a perfusion pouch with a luer fitting to allow attachment of a syringe. The perfusion pouch facilitates homogeneous saturation of chronOS Strip Bone Void Filler implants with the patient's bone marrow and/or blood.

## **Steps for Using the Perfusion Pouch**

1. Aspirate bone marrow using a standard bone marrow aspiration system or fill a syringe with the patient's blood.
2. Unscrew protective cap on the perfusion pouch's port.
3. Attach the syringe containing autologous bone marrow or blood to the chronOS Strip Bone Void Filler perfusion pouch using the port with the luer fitting.
4. Inject the perfusion pouch with medium of choice following the recommended volume ranges to ensure adequate wetting of the chronOS Strip Bone Void Filler surface.


chronOS Strip Bone Void Filler Dimension	Perfusion Volume Range
50 mm × 25 mm × 3 mm	4–10 ml
100 mm × 25 mm × 3 mm	7–20 ml
50 mm × 25 mm × 6 mm	7–10 ml
100 mm × 25 mm × 6 mm	15–20 ml
47 mm × 18 mm × 3 mm (2 implants)	5–20 ml

**Note:** Do not exceed the recommended upper volume limit as this may over-pressurize the perfusion pouch.

- Perfuse chronOS Strip Bone Void Filler by gently pumping the syringe plunger 3–5 times or until the implant is visibly wetted on all surfaces.

**Note:** Do not pull the syringe plunger back to its maximum stroke as this may cause the plunger to disengage from the syringe barrel, resulting in the loss of syringe contents.

- Remove the syringe to relieve any pressure from the perfusion pouch.
- When ready to use the implant, remove it by grasping the implant through the perfusion pouch (to secure it) while tearing the designated notch.



- Implant chronOS Strip Bone Void Filler immediately or place in a sterile bowl for later use during the same surgical procedure.

**Note:** Remove the tabs of chronOS Strip Bone Void Filler precast prior to implantation.

### Device Related Storage and Handling Information

Keep the sterile bone substitute in its protective package in a dry place below 25°C.

### Additional Device-Specific Information

**REF**

Reference number

**LOT**

Lot or batch number

Manufacturer

Expiration date

**0123**

Notified body

Caution, see instructions for use

Do not use if package is damaged

### Additional Instructions for Use

For additional information, please contact your local Synthes representative.

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© Synthes GmbH 2018. All rights reserved. 60090701 SE\_460279 AD 06/2018 70100238
Synthes GmbH