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XRL Medium

Introduction

Please read these instructions for use, the enclosed Important Information and the corresponding technique guide carefully before use. Ensure that you are familiar with the appropriate surgical technique.

The modular XRL implant is made up of five components:

- 1× XRL Medium Central Body
- 2× XRL Medium Endplates
- 2× XRL Medium Endplate Screws

The XRL implant is assembled before implantation and consists of five components (2 endplates, 2 endplate screws and 1 central body). A cranial and caudal endplate are attached with endplate screws onto the XRL central body. The XRL central body must never be used without endplates.

Smaller heights are only available as integrated implants. These implants do not require assembly of endplates and set screws.

Material(s)

Material	Standard
PEEK (Polyetherether-ketone)	ASTM F 2026
TAN (Ti-6Al-7Nb)	ISO 5832-11
Tantalum	ISO 13782, ASTM F 560

Intended use

XRL serves as a replacement for injured vertebral bodies to stabilize the anterior thoracic and lumbar spine (Th3-L5). The approach can be anterior, anterolateral, lateral or posterolateral.

Depending on the pathological situation, a defect height of 22 mm to 142 mm can be corrected, allowing XRL to be used for one and two level vertebral body defects.

Indications

- Fractures with destruction of the anterior column in the thoracic and lumbar spine
- Post traumatic malalignment
- Replacement of vertebral bodies following tumor resection in the thoracic and lumbar spine
- Reconstruction of the anterior column following an infection

Contraindications

- For sole anterior treatment of advanced osteoporosis. In this case an additional dorsal treatment should be carried out in addition to the anterior treatment
- Diffuse spinal tumors
- In the absence of intact neighboring segments

Note

- As with all vertebral body replacement systems, XRL is not suitable as a stand-alone device, and must be used together with an internal fixation system (such as DePuy Synthes Supplemental Fixation) to absorb compression forces, tensile forces and torsional moments.
- Additional internal fixation system must not be removed before solid bony fusion occurs.

Potential risks

As with all major surgical procedures, risks, side effects and adverse events can occur. While many possible reactions may occur, some of the most common may include:

Problems resulting from anesthesia and patient positioning (e.g. nausea, vomiting, dental injuries, neurological impairments, etc.), thrombosis, embolism, infection, excessive bleeding, iatrogenic neural and vascular injury, damage to soft tissues incl. swelling, abnormal scar formation, functional impairment of the musculoskeletal system, Complex regional pain syndrome (CRPS), allergy/hypersensitivity reactions, side effects associated with implant or hardware loosening, malunion, non-union, ongoing pain; damage to adjacent bones (e.g. subsidence), discs (e.g. adjacent level degeneration), or soft tissue, dural tear or spinal fluid leak; spinal cord compression and/or contusion, partial displacement of the graft, vertebral angulation.

Sterile device

Sterilized using irradiation

Store devices in their original protective packaging, and do not remove them from the packaging until immediately before use.

Prior to use, check the product expiration date and verify the integrity of the sterile packaging. Do not use if the package is damaged.

Do not resterilize

Single use device

Products intended for single-use must not be re-used.

Reuse or reprocessing (e.g. cleaning and resterilization) may compromise the structural integrity of the device and / or lead to device failure which may result in patient injury, illness or death.

Furthermore, reuse or reprocessing of single use devices may create a risk of contamination e.g. due to the transmission of infectious material from one patient to another. This could result in injury or death of the patient or user.

Contaminated implants must not be reprocessed. Any Synthes implant that has been contaminated by blood, tissue, and/or bodily fluids/matter should never be used again and should be handled according to hospital protocol. Even though they may appear undamaged, the implants may have small defects and internal stress patterns that may cause material fatigue.

Warnings

- In addition to risks associated with any surgical procedure, adverse effects when using this instrumentation system can include:
 - device component fracture
 - loss of fixation
 - non-union
 - inflammation*
 - fracture of the vertebrae*
 - neurological injury*
 - vascular or visceral injury*
- The components of this device are manufactured from PEEK, titanium alloy and tantalum. For metallurgical, mechanical, and functional reasons, devices other than the provided components (e.g., supplemental fixation) should not contact the Synthes XRL construct if made from materials not listed here. Also, mixing titanium or titanium alloy with stainless steel implant components is not recommended for metallurgical, mechanical and functional reasons.
- Excessive strain by patient prior to sufficient bone bridge formation can lead to implant failure.
- Implants are not suitable to replace normal body structures or bear the body's weight in the cases of e.g. bad bone quality, infection, tumor or other disease of the adjacent vertebral bodies or for the case of a pseudoarthrosis (non-union).
- The performance of this product may be negatively influenced by patients whose activity, mental capacity, mental illness, alcoholism, drug abuse, occupation or lifestyle may interfere with their ability to follow postoperative restrictions, thereby placing undue stresses on the implant during bony healing.
- The manufacturer is not responsible for any complications arising from incorrect diagnosis, choice of incorrect implant, incorrectly combined implant components and/or operating techniques, the limitations of treatment methods, or inadequate asepsis.
- The implant components applied (name, article number, lot number) must be documented in each patient's record.
- Do not implant when an active systemic infection is present.

*These risks can also be associated with general risks of surgery.

Precautions

- The implantation of vertebral body replacement devices and supplemental internal fixation are technically demanding procedures that present a risk of serious injury to the patient. Therefore, only experienced spinal surgeons with specific training in the use of this system should perform these procedures. Preoperative planning and patient anatomy should be considered when selecting implant size and supplemental internal fixation.
- Do not use XRL for hypersensitive patients as this can lead to inflammation, resulting in medical treatment or re-operation and removal of implant.

Magnetic Resonance environment

MR Conditional

Non-clinical testing of the worst-case scenario has demonstrated that the implants of the XRL Medium system are MR conditional. These articles can be scanned safely under the following conditions:

- Static magnetic field of 1.5 Tesla and 3 Tesla.
- Spatial gradient field of 300 mT/cm (3000 Gauss/cm).
- Maximum whole body averaged specific absorption rate (SAR) of 4 W/kg for 15 minutes of scanning.

Based on non-clinical testing, the XRL Medium implant will produce a temperature rise not greater than 4.1 °C at a maximum whole body averaged specific absorption rate (SAR) of 4 W/kg, as assessed by calorimetry for 15 minutes of MR scanning in a 1.5 Tesla and 3 Tesla MR scanner.

MR Imaging quality may be compromised if the area of interest is in the exact same area or relatively close to the position of the XRL Medium device.

Special operating instructions

Selection

The selected endplate size should provide the widest support of the neighboring vertebral bodies. Make sure that the endplates of the

implant do not project in any direction beyond the endplates of the neighboring vertebral bodies. This could cause serious vascular or neurological damage.

The angle of the endplate should correspond to the anatomy and pathology of the respective patient. A wrong angle or size could cause the implant to subside into the bone, or prevent a successful correction.

Make sure that the central body is expanded with the corresponding spreader top.

Fig. 1

Mounting

When mounting the endplates, observe the following:

- The access to release the locking ring in the central body has to face the direction of the desired approach.
- Ensure that both endplates are in the same direction.
- When pressing on the endplates, ensure the endplate properly seats on the central body. This can be checked visually. If the endplate is not properly seated, there is a risk that it could detach from the central body (see Fig. 1).
- The modular XRL central body must never be implanted without cranial and caudal endplates properly secured with endplate screws.
- Please use the torque limiting handle according to the XRL technique guide to ensure appropriate screw locking.
- Only use the provided instruments as described in the XRL technique guide to reposition or to attach the endplates.

Implantation

Observe the following points when using the XRL implant:

- With the locking ring facing the instrument, slide the spreader top into the slots below the cranial endplate. Do not force spreader top onto implant. Slightly turn the T-driver clockwise until the detent on the fork of the spreader shaft engages the implant for a secure hold. Set the scale to zero.
- Completely insert the release tool through the XRL spreader and into the locking ring.
- The central body must not be excessively expanded. Otherwise the danger exists that the implant may be pressed into the neighboring vertebral body.
- When the implant is in its final position, verify the locking ring on the central body is closed. When the slot is approximately 1 mm (Fig. 2), the implant is locked and secured. If the slot is larger (Fig. 3), re-engage the implant with the spreader, and with the release tool disengaged, distract the implant slightly to close the locking ring.
- Check the implant for proper alignment after implantation.

Fig. 2


Fig. 3

For detailed instructions, please refer to the XRL technique guide.

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