Instructions for Use
Gebrauchsanweisung
Instrucciones de uso
Mode d'emploi
Istruzioni per l'uso
Instruções de utilização
Návod na použitie

09.804.500S–502S,
Vertebral Body Stent

09.804.600S–602S,
Vertebral Body Stent with Vertebral Body Balloon

PART OF THE JOHNSON & JOHNSON FAMILY OF COMPANIES

Description

The VBS System consists of the Vertebral Body Stent (VBS), the optional Vertebral Body Balloon (VBB), the access kit, and the Inflation System.

This leaflet contains information about the following products:

- 09.804.500S – 502S, VBS, containing: one stent (L605 cobalt-chromium alloy), one balloon-catheter and one stiffening wire (SST, Inox) and
- 09.804.600S – 602S, VBS with VBB, containing: one stent (L605 cobalt-chromium alloy), two balloon-catheters and two stiffening wires (SST, Inox)

The access kit (03.804.512S/03.804.612S) is used to prepare the operative access in the vertebral body. Subsequently, the Vertebral Body Stent is inserted into the vertebral body using a simultaneous bilateral approach. The Inflation System (03.804.413S) is then used to inflate the balloon, thereby expanding the stent. Once the vertebral body is restored to the desired height, the balloon is deflated and removed from the vertebral body. The stent remains in situ and stabilizes the cavity that has been created. The access kit (03.804.512S/03.804.612S) is then used to inject the PMMA based bone cement. As an option when using 09.804.600S-602S the enclosed VBB allows an in-situ preparation of the vertebral body prior to the use of VBS. Please read the surgical technique guide carefully prior to use of the Vertebral Body Stent System.

The Vertebral Body Stent System is supplied sterile in removable packaging. Do not use the product if the sterile packaging is damaged and inform the manufacturer.

The access kit and the Inflation System are designed specifically for use with the VBS System by Synthes. Alternate instrumentation must not be used with the VBS System.

Refer to the Instructions for Use for the access kit and Inflation System for additional details regarding these devices. Additionally, please adhere to the Instructions for Use for the specific PMMA based bone cement used during the procedure.

Intended use

The VBS System is intended for the reduction of painful vertebral compression fractures and/or creation of a void in cancellous bone in the spine for the treatment of levels ranging from Th5-L5. It is intended to be used in combination with a legally-marketed PMMA¹ based bone cement adequately indicated for use in vertebroplasty or kyphoplasty procedures.

Note: Refer to the manufacturer's directions accompanying the bone cement for specific information on its use, precautions and warnings.

Indications

- Painful osteoporotic vertebral compression fractures without posterior wall involvement. Classified after Genant, Grade 2 and Grade 3.
- Painful vertebral fractures classified after the AO classification:

A1.1 Endplate impaction

A1.2 Wedge impaction fracture

A1.3 Vertebral body collapse

A3.1 Incomplete burst fracture; matter of discretion (depending on the degree of posterior wall involvement, internal fixation must be used in addition)

In combination with internal fixation

A3.1 Incomplete burst fracture

A3.2 Burst-split fracture; matter of discretion (the extent of the gap width should not be too wide)

B1.2 Posterior disruption predominantly ligamentous associated with type A fracture of the vertebral body

B2.3 Posterior disruption predominantly osseous with type A fracture of the vertebral body.

- Palliative treatment of osteolytic lesions located within the vertebral body with intact cortical shell. Classified after Tomita Type 1.

¹ **Note:** Due to limited long-term efficacy data, the treating physician should weigh the benefits of the application of the PMMA based bone cement in younger patients against the potential risks.

Contraindications

- Lesions requiring open anterior column reconstruction
- Acute or chronic systemic or localized spinal infections

Cautions

- Before using the VBS System ensure that the size is suitable for the specific procedure.
- The surgeon is responsible for ensuring that the operation is carried out properly.
- The product may only be used by surgeons sufficiently trained in the application of the VBS System.
- General risks of surgery are not described in this document.
- The VBS System may only be used with an X-ray control with a device that offers a high image quality.
- The balloon pressure of the VBS and VBB may not exceed the maximum inflation pressure of 30 bar/atm. A manometer is used to monitor the pressure.
- The inflation balloon volumes of the VBS and VBB must not exceed the maximum volumes specified in the surgical technique guide.
- Only inflate the balloon with liquid, water-soluble, ionic or non-ionic contrast medium (VBS/VBB has been tested with a maximum iodine concentration of 320mg/ml). Contrast media may have different viscosity and precipitation levels that may influence inflation and deflation times, therefore a mixture ratio of contrast medium to saline solution of 1:2 is recommended.

It is essential to observe the manufacturer's instructions on the indications, use and safety measures for the contrast agent.

- The patient has to be checked for allergy to the contrast medium and stent material, i.e. any of the metal components of the CoCrW_{Ni} alloy.
- Do not use air or other gases to inflate the balloon-catheters of the VBS or the VBB.
- Never expose the balloon-catheter to organic solvents (e.g. alcohol).

- The efficacy of the balloon-catheter may be adversely affected if it comes into contact with bone splinters, bone cement and/or surgical instruments.
- The VBS System is designed for single use only and may not be resterilized and/or reused. It is expressly prohibited to recondition, reprocess, repair, modify or resterilize any component for further use. Even though the components may appear undamaged, they may have small defects and internal stress patterns that may cause material fatigue.
- The manufacturer is not liable for complications resulting from improper indication, diagnosis, incorrectly selected implants, incorrect combination of implant components and surgical techniques or limitations of the treatment methods or failure to ensure aseptic conditions.

Warnings

- Do not use the VBS System after the expiration date or if the packaging is opened or damaged.
- Do not leave the balloon implanted; the balloon material is not implant grade material.
- The timing for release of the bone cement is dependent on the selection of bone cement.

Possible side effects

- Deformation, loosening, wear or breakage of the implant components
- Transport of fat, thrombus or other materials which can in turn lead to a symptomatic pulmonary embolism or other clinical consequences
- Rupture and collapse of the inflated balloon-catheter and retention of a fragment in the vertebral body
- Rupture of the balloon-catheter with contrast agent exposure and possible subsequent allergic reaction
- Neurological damage
- Damage to vessels, nerve and organs
- Early and late infections
- Allergic reactions to the implant materials
- Tissue reactions to the implant materials

- Venous thromboses, pulmonary embolism and cardiac arrest
- Haematomas and impaired wound healing
- Rebounding fragments of the vertebral body that cause damage to the bone marrow or the nerve roots and can thus result in radiculopathy, paresis or paralysis
- Death

Disinfection, cleaning, care, storage and sterilization

- The implant components are EO sterilized.
- The implant components should be stored in their original packaging and should not be removed from the protective packaging until immediately before use. Store in a cool, dry place.
- The expiration date should be checked and it is also important to ensure that the sterile packaging is not damaged.
- Never use damaged or surgically removed implant components.
- The implants are not designed to be disinfected or cleaned by the user.

MRI information

Testing demonstrated that the VBS stent is MR conditional. A patient with this device can be scanned safely immediately after placement under the following conditions:

- Static magnetic field of 3-Tesla or less
 - Spatial gradient magnetic field of 720-Gauss/cm or less
 - Maximum whole body averaged specific absorption rate (SAR) of 3.0 W/kg for 15 minutes of MR imaging
- The VBS stent was determined to be MR conditional according to the terminology specified in the American Society for Testing and Materials (ASTM) International, Designation: F2503-05. Standard Practice for Marking Medical Devices and Other Items for Safety in the Magnetic Resonance Environment.

MR Imaging quality may be compromised if the area of interest is in the exact same area or relatively close to the position of the stent.

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