**Instructions for Use  
for chronOS Bone Void Filler Implants  
(chronOS Granules Bone Void Filler, chronOS  
Preforms Bone Void Filler, chronOS Vivify  
Bone Void Filler and chronOS Inserts Bone  
Void Filler)**

This instruction for use is not intended  
for distribution in the USA.

## English

This instruction for use is not intended for distribution in the USA. Please read these instructions for use, the Synthes "Important Information" and the corresponding surgical technique carefully before use. Ensure that you are familiar with the appropriate surgical technique. Sufficient experience in the field of the use of biomaterials is highly recommended before using chronOS Bone Void Filler implants.

### Material(s):

$\beta$ -tricalcium phosphate [ $\text{Ca}_3(\text{PO}_4)_2$ ]

### Standard(s):

ASTM-1088

### Intended Use

chronOS Bone Void Filler implants are a bone replacement material of synthetic, porous and resorbable  $\beta$ -tricalcium phosphate ( $\beta$ -TCP). In its function as temporary bone substitute, it serves to fill and bridge bone defects in children, adolescents and adults.

### Indications

chronOS Bone Void Filler implants should be used as bone void fillers or augmentation material in zones requiring cancellous rather than cortical bone. This includes the filling of bone defects after trauma, reconstruction or correction in non-load bearing indications only.

The porous structure of chronOS Bone Void Filler implants acts as a matrix for the ingrowth of bone. chronOS Bone Void Filler implants must always be applied by enosseal or subperiosteal implantation, i.e. by direct contact with the vital bone.

As a general rule, it is recommended that small granules are used in small defects and larger granules in larger bone defects:

- 0.5–0.7 mm granule size for filling defects up to 1.0 cc
- 0.7–1.4 mm granule size for filling defects up to 2.5 cc
- 1.4–2.8 mm and 2.8–5.6 mm granule size for filling defects from 2.5 cc

The implant must completely fill the bone defect (press-fit). However, it is essential to avoid overfilling to ensure tension-free wound closure. It is recommended that the filled defects are covered with a periost flap or a resorbable sheet to avoid soft tissue infiltration.

A mixture with autologous and/or allogeneous bone is highly recommended for larger defects (more than 20 cc). The size of the chronOS Bone Void Filler implant used depends on the size of the bone defect.

As a rule, chronOS Bone Void Filler implants are resorbed within 6 to 18 months and converted into bone; depending on patient conditions.

## **Contraindications**

chronOS Bone Void Filler implants should not be used in the following circumstances:

- Acute and chronic infections in the operation area, e.g. inflammation, bacterial bone diseases (posttraumatic or chronic osteomyelitis) and soft-tissue infections
- Malignant myeloma, Burkitt's lymphoma, and other lymphomas
- Defects and fractures in the region of an open epiphysis
- Osteoporosis
- Severe instability or deformation at the extraction point (harvesting site)
- Load bearing and unstable applications

## **Side Effects**

- Nonunion or delayed union, which may lead to a failure of the implant
- Pain, discomfort, abnormal sensation, or palpability due to the presence of the device
- Increased fibrous tissue response around fracture site and/or the implant

Apart from these possible adverse effects, there is also the risk of complications associated with any surgical procedure featuring bony defects such as, but not limited to, necrosis of bone, infection, nerve damage and pain which may not be related to the implant.

In general, good tissue response of  $\beta$ -tricalcium phosphate implants in bone is supported by experimental and clinical data.

Nevertheless, the following complications are possible:

- Fragment displacement as a result of use in inappropriate indications
- Neurovascular injuries caused by surgical trauma
- Foreign body reactions
- Allergic reactions
- Inflammatory reactions
- Infections that can lead to failure of the procedure
- General complications caused by invasive surgery

### **Sterile Device**

<b>STERILE R</b>	Sterilized using gamma-irradiation
------------------	------------------------------------

- Store implants in their original protective packaging, and do not remove them from the packaging until immediately before use.
- Prior to use, check the product expiration date and verify the integrity of the sterile packaging.
- Do not use if the package is damaged.

Do not re-sterilize

- This product is intended for single use and must not be re-used. Re-use or re-processing (e.g. cleaning and/or re-sterilization) may compromise the sterility or integrity of the device which may lead to infection, device failure and/or contamination.

### **Single Use Device**

Do not re-use

- This product is intended for single use and must not be re-used. Re-use or re-processing (e.g. cleaning and/or re-sterilization) may compromise the sterility or integrity of the device which may lead to infection, device failure and/or contamination.
- The contents of an opened or damaged outer protective packaging may no longer be used and must be destroyed.
- Once the packaging has been opened, the content must be used immediately. Any unused parts must be destroyed.
- Do not use after the expiry date.

## **Conditional Restrictions**

The use of chronOS Bone Void Filler implants has a more restricted indication in:

- Severe, endocrine-induced bone diseases (e.g. hyperparathyroidism)
- Current therapy with steroids and with drugs, which intervene in calcium metabolism (e.g. calcitonin)
- Severe, poorly controlled diabetes (diabetes mellitus) with bad wound healing tendencies
- Immunosuppressive therapy
- Poor bone quality

The use of chronOS Bone Void Filler implants is not known in:

- Cardiovascular diseases
- Pregnancy and lactation period
- Filling of bone defects in congenital deformations
- Immunosuppressive or radiation therapy
- Compromised immune system, compromised wound healing

## Warnings

- Do not use chronOS Bone Void Filler implants as a stand-alone product unless the cortical bone can bear the full load. Maximum mechanical stability in filling the defect with chronOS Bone Void Filler implants is a prerequisite for a good bony incorporation. If there is instability, use appropriate internal fixation to first stabilize the operating area.
- Do not use 2.8–5.6 mm granule size for filling of sinus lifts.
- Do not attempt to re-sterilize the unused contents of an opened pack, but dispose of such remnants. Re-sterilizing of chronOS Bone Void Filler implants can result in product not being sterile and/or meet performance specifications.
- Do not fix chronOS Bone Void Filler implants with screws or tacks directly in the bone.
- The perfusion pouch is intended for the impregnation of chronOS Vivify Preforms Bone Void Filler with bone marrow or blood.
- chronOS Inserts Bone Void Filler for Cervios must only be used in combination with Synthes Cervios Interbody Fusion Devices.
- Do not exceed the recommended upper volume limit as this may over pressurize the perfusion pouch.

## **Combination of the Medical Device with a Medicinal Product**

The effect of mixing chronOS Bone Void Filler implants with medicinal products is not known.

## **Interactions**

No negative interactions with bone marrow and autologous blood have been reported as yet.

## **Magnetic Resonance (MR) Environment**

chronOS Bone Void Filler implants are of non-metallic inorganic origin. This material is inherently diamagnetic and cannot be heated up or act as an antenna either by bringing the patients into the MRI magnet or during the MRI examinations. However, chronOS Bone Void Filler implants may be used in combination with metal hardware. Unless stated otherwise, devices have not been evaluated for safety and compatibility within the MR environment.

Please note that there are potential hazards which include but are not limited to:

- Heating or migration of the device
- Artifacts on MR image

## **Special Operating Instructions**

In removing the implant from its packaging, strictly observe the instructions concerning asepsis.

chronOS Bone Void Filler implants must always be applied by enosseal or subperiosteal implantation, i.e. by direct contact with the vital bone. The surgical procedure depends on the localization, nature and extent of the bone defect. To prepare the implant site, remove all inflamed necrotic tissue and bone particles.

The appropriate implant size is determined by the size of the bone defect, the vascularity and the size of the cancellous bone chips, if used. Keep the regions of enosseal vessel and nerve cords clear to avoid pressure sores.

### **chronOS Granules Bone Void Filler**

It is useful to mix chronOS Granules Bone Void Filler with autologous bone marrow or blood for improved handling. Introduce the resulting mixture into the defect. During application, remove individually displaced granules from the soft tissue. Based on experimental and clinical data, it is considered as safe to mix chronOS Granules Bone Void Filler with autologous or allogenic bone.

### **chronOS Inserts Bone Void Filler**

chronOS Inserts Bone Void Filler must be inserted into the Cervios Interbody Fusion Device in alignment with the orientation of the implant graft window.

## **chronOS Vivify Preforms Bone Void Filler**

chronOS Vivify Preforms Bone Void Filler are available in a perfusion pouch. The use of the perfusion pouch facilitates the enrichment of implants with the patient's own bone marrow or blood.

### **Perfusion Pouch**

chronOS Vivify Preforms Bone Void Filler are provided in a perfusion pouch with a Luer fitting to allow attachment of a syringe.

Steps for using the perfusion pouch:

1. Aspirate bone marrow using a standard bone marrow aspiration system or fill a syringe with the patient's blood.
2. Unscrew the protective cap on port of the perfusion pouch.
3. Attach the syringe containing autologous bone marrow or blood to the chronOS perfusion pouch using the port with the Luer fitting.
4. Fill the perfusion pouch with medium of choice following the recommended volume ranges to ensure adequate wetting of the chronOS Vivify Preforms Bone Void Filler surface.

### **chronOS Vivify Cylinders Bone Void Filler**

Diameter (mm)	Perfusion Pouch	Perfusion Volume Range (mL)
8.5	Small	5–8
9.5	Small	5–8
10.5	Small	5–8
12.5	Large	8–10
14.0	Large	8–10

### **chronOS Vivify Blocks Bone Void Filler**

Size (mm)	Perfusion Pouch	Perfusion Volume Range (mL)
5×5×10	Small	2–5
12.5×12.5×10	Small	5–8
20×20×10	Small	8–10

### **chronOS Vivify Wedges Bone Void Filler**

Angle	Perfusion Pouch	Perfusion Volume Range (mL)
10°	Small	5–8
14°	Small	8–10
18°	Small	8–10
22°	Large	8–10
26°	Large	8–10

### **chronOS Vivify Wedges Bone Void Filler, semi-circular**

Angle	Perfusion Pouch	Perfusion Volume Range (mL)
7°	Large	8–10
10°	Large	8–10
13°	Large	10–12

Do not exceed the recommended upper volume limit as this may over pressurize the perfusion pouch.

5. Perfuse chronOS Vivify Preforms Bone Void Filler by gently pumping the syringe plunger 10–12 times or until the implant is visibly wetted on all surfaces.

Allow a flow back of bone marrow or blood into the syringe before pumping it again into the pouch.

Do not pull the syringe plunger back to its maximum stroke as this may cause the plunger to disengage from the syringe barrel, resulting in the loss of syringe contents.

6. Remove the syringe to relieve any pressure from the perfusion pouch.
7. When ready to use the implant, remove it by grasping the implant through the perfusion pouch (to secure it) while tearing the designated notch.

8. Implant chronOS Vivify Preforms Bone Void Filler immediately or place in a sterile bowl for later use during the same surgical procedure.

## Device-related Storage and Handling Information

Keep dry

Keep away from sunlight

## Additional Device-specific Information

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© Synthes GmbH 2016. All rights reserved.
Synthes GmbH