

# Variable Angle LCP Dorsal Distal Radius Plate 2.4.

For fragment-specific fracture fixation with variable angle locking technology.

## Surgical Technique



This publication is not intended for distribution in the USA.

Instruments and implants approved by the AO Foundation.

---

 Image intensifier control

This description alone does not provide sufficient background for direct use of DePuy Synthes products. Instruction by a surgeon experienced in handling these products is highly recommended.

**Processing, Reprocessing, Care and Maintenance**

For general guidelines, function control and dismantling of multi-part instruments, as well as processing guidelines for implants, please contact your local sales representative or refer to:

<http://emea.depuysynthes.com/hcp/reprocessing-care-maintenance>

For general information about reprocessing, care and maintenance of Synthes reusable devices, instrument trays and cases, as well as processing of Synthes non-sterile implants, please consult the Important Information leaflet (SE\_023827) or refer to:

<http://emea.depuysynthes.com/hcp/reprocessing-care-maintenance>

# Table of Contents

---

Introduction	<b>Variable Angle LCP Dorsal Distal Radius Plate 2.4</b>	2
	<b>AO Principles</b>	4
	<b>Intended Use and Indications</b>	5
	<b>Clinical Cases</b>	6
	<b>Three-Column Theory</b>	7
<hr/>		
Surgical Technique	<b>Recommendations on Screw and Plate Insertion</b>	8
	Screw Insertion Techniques	8
	Screw Type Determination	9
	Plate Insertion Technique	10
	<b>Preparation</b>	12
	<b>Approach</b>	13
	<b>Plate Insertion</b>	14
	<b>Screw Insertion</b>	16
	<b>Postoperative Treatment / Implant Removal</b>	25
<hr/>		
Product Information	<b>Plates</b>	26
	<b>Screws</b>	28
	<b>Instruments</b>	30
<hr/>		
Bibliography		32
<hr/>		
MRI Information		33

## Variable Angle LCP Dorsal Distal Radius Plate 2.4. For fragment-specific fracture fixation with variable angle locking technology.

The low profile dorsal distal radius plates are intended for the double-plate technique. All implants are available in stainless steel and titanium.


### Variable angle locking

Holes allow up to 15° off-axis screw angulation in all directions in order to address the individual fracture patterns.




### Undercuts and bending notches


Undercuts and bending notches allow easy contouring of the plates while preserving the VA locking holes.



### Low profile construct


Rounded edges, polished surface and countersunk screws help reduce the risk of soft tissue irritation.





**Anatomical fit**

Anatomical pre-contoured plates for the radial and intermediate column minimize the need for bending.



**Oblong VA combi-hole**

Allows accurate plate positioning on the bone.



**Kirschner wire holes**

Enable preliminary plate fixation.

# AO Principles


In 1958, the AO formulated four basic principles, which have become the guidelines for internal fixation<sup>1,2</sup>.

## Anatomic reduction

Fracture reduction and fixation to restore anatomical relationships.

## Early, active mobilization

Early and safe mobilization and rehabilitation of the injured part and the patient as a whole.



## Stable fixation

Fracture fixation providing absolute or relative stability, as required by the patient, the injury, and the personality of the fracture.

## Preservation of blood supply

Preservation of the blood supply to soft tissues and bone by gentle reduction techniques and careful handling.

<sup>1</sup> Müller ME, Allgöwer M, Schneider R, Willenegger H. Manual of Internal Fixation. 3<sup>rd</sup> ed. Berlin, Heidelberg, New York: Springer. 1991.

<sup>2</sup> Rüedi TP, Buckley RE, Moran CG. AO Principles of Fracture Management. 2<sup>nd</sup> ed. Stuttgart, New York: Thieme. 2007.

# Intended Use and Indications

---

## **Intended Use**

The plate and screw implants included in the Radius Plate product family are intended for temporary fixation, correction or stabilization in the radius anatomical region.

## **Indications**

2.4 mm Variable Angle LCP Dorsal Distal Radius Plates are indicated for:

- Dorsally displaced fractures
- Extra-articular and intra-articular fractures with metaphyseal defect
- Open joint reconstruction
- Combination of distal radius with carpal and metacarpal fractures
- Corrective osteotomies after distal radius malunion

# Clinical Cases

Case 1  
69-year-old female  
with AO 23A2  
fracture



Preoperative, AP view



Preoperative, lateral view




Postoperative, AP view



Postoperative, lateral view

Case 2  
59-year-old male  
with AO 23C2  
fracture



Preoperative, AP view



Preoperative, lateral view



Postoperative, AP view



Postoperative, lateral view

# Three-Column Theory

The treatment of distal radius fractures requires a meticulous reconstruction of the joint surface, as well as stable internal fixation and early functional postoperative treatment.

Extra-articular fractures require both the restoration of the volar tilt and radial length to reduce the possibility of displacement. Any malalignment may result in limitations of movement, changes of load distribution, mid-carpal instability as well as the increased risk of osteoarthritis in the radio-carpal joint.


Intra-articular fractures with articular displacement of more than 2 mm in the radiocarpal joint inevitably result in osteoarthritis and functional impairment.

The distal radius and distal ulna form a three-column biomechanical construction<sup>3</sup>:

- The ulnar column is the distal ulna, the triangular fibrocartilage and the distal radio-ulnar joint.
- The intermediate column is the medial part of the distal radius, with the lunate fossa and the sigmoid notch.
- The radial column is the lateral part of the radius with the scaphoid fossa and the styloid process.

A dorsally displaced fracture of the distal radius indicates not only dorsiflexion in the sagittal plane, but also radial deviation in the frontal plane and supination in the transverse plane.

Following reduction, stabilization requires optimal fixation of the intermediate column as well as the radial column. In the case of a fractured distal ulna that compromises the distal radio-ulnar joint, the ulnar column must be stabilized as well.



- Radial column
- Intermediate column
- Ulnar column


<sup>3</sup> Rikli DA, Regazzoni P (1996) Fractures of the distal end of the radius treated by internal fixation and early function. A preliminary report of 20 cases. J Bone Joint Surg [Br] 78 (4): 588–592

# Recommendations on Screw and Plate Insertion

## Screw Insertion Techniques

Variable angle locking screws can be inserted using two different techniques:

- Variable angle technique
- Pre-defined nominal angle technique





### a) Variable angle technique

To drill variable angle holes up to 15° deviation from the nominal trajectory of the locking hole, insert the tip of the VA-LCP drill sleeve and key into the cloverleaf design of the VA locking hole. (1)

Use the funnel-shaped end of the VA-LCP drill sleeve to drill variable angle holes at the desired angle. (2)

Alternatively, use the freehand VA-LCP drill sleeve and insert it fully into the VA locking hole. (3) Drill variable angle holes at the desired angle. (4)



**Precaution:** It is important not to angulate more than 15° from the central axis of the screw hole. Overangulation could result in inappropriate screw-locking. Moreover, the screw head may not be fully countersunk.

### b) Pre-defined nominal angle technique

The fixed-angle end of the VA-LCP drill sleeve only allows the drill bit to follow the nominal trajectory of the VA locking hole.



## Screw Type Determination


Determine whether standard cortex screws or variable angle locking screws will be used for fixation.

The final screw placement and the use of VA locking and cortex screws are determined by the fracture pattern.

If a VA locking screw is inserted first, ensure that the plate is held securely against the bone to prevent the plate from spinning as the screw locks into the plate.

When using the pre-defined nominal angle technique standard locking screws can also be used instead of VA locking screws.

**Precaution:** The screw head is not completely countersunk if a cortex screw is inserted in a variable angle locking hole.



**VA locking hole:**

2.4 mm VA locking screw,  
1.8 mm VA locking buttress pin,  
2.4 mm locking screw  
(only nominal angle)  
or 2.4 mm cortex screw applicable

**Oblong VA combi-hole:**

2.4 mm cortex screw applicable in the  
compression portion (1), 2.4 mm VA  
locking screw, 1.8 mm VA locking  
buttress pin, 2.4 mm locking screw (only  
nominal angle) or 2.4 mm cortex screw  
applicable in the threaded portion (2)

## **Plate Insertion Technique**

### **1. Apply dorso-ulnar plate**

It is recommended to apply the dorso-ulnar plate first and fix it by inserting a 2.4 mm standard cortex screw in the oblong VA combi-hole in the proximal shaft. (See pages 16–17 for insertion of cortex screws.)

The plate supports the intermediate column and fixes the dorso-ulnar fragment.



### **2. Apply dorso-radial plate**

Apply the dorso-radial plate after provisional positioning of the dorso-ulnar plate by inserting a 2.4 mm standard cortex screw in the oblong VA combi-hole in the proximal shaft. It should form an angle of approximately 70° to the dorso-ulnar plate. (See pages 16 – 17 for insertion of cortex screws.)

Use the small notch (horse-shoe tip) in the distal end of the plate to position the dorso-radial plate properly.




### 3. Insert screws in dorso-ulnar plate

Insert a VA locking screw in the most proximal hole in the shaft of the dorso-ulnar plate (a). Complete internal fixation by inserting VA locking screws in the distal arm of the plate (b, c). (See pages 18 – 24 for insertion of variable angle locking screws.)



### 4. Insert screws in dorso-radial plate

Insert a VA locking screw in the most proximal hole in the shaft of the dorso-radial plate (d). Complete internal fixation by inserting VA locking screws in the distal arm of the plate (e, f). (See pages 18 – 24 for insertion of variable angle locking.)



# Preparation

---

## Select implants

Select the plates according to the fracture pattern and anatomy of the bone.

---

**Note:** This surgical technique describes the application of the VA-LCP Dorsal Distal Radius Plates 2.4 using the plates 04.115.151 and 04.115.540.

---


# Approach

Make a straight longitudinal incision over the dorsal distal radius extending 5 to 10 cm between the second and third dorsal extensor compartments. Open the extensor retinaculum by performing a longitudinal incision between the first and second extensor compartments.

Take care to elevate and mobilize the third compartment (extensor pollicis longus) proximally and distally, and translocate it radially for better access to the fracture site.

Elevate the second and fourth dorsal compartments subperiosteally to preserve their integrity.

For additional information on technique alternatives see Rikli (2005)<sup>4</sup>.



<sup>4</sup> Rikli DA, Businger A, Babst R (2005) Dorsal double-plate fixation of the distal radius. *Oper Orthop Traumatol* 17(6): 624 – 640

# Plate Insertion

## 1

### Reduce fracture

- Reduce the fracture under image intensifier control and, if necessary, fix with Kirschner wires or reduction forceps. The reduction method will be fracture-specific.

## 2

### Contour plate

#### Instrument

347.901	Pliers, flat-nosed, pointed for Plates 1.0 to 2.4
---------	---

If necessary, twist and bend the plate to suit anatomical conditions as indicated. Avoid repetitive bending. (1)

The anatomical pre-contoured plates (0X.115. 530 – 0X.115.641) do not usually require any contouring. (2)

Recommendation: Use non-serrated bending pliers for preservation of the plate's smooth finish.

#### Precautions:

- The design of the plate holes allows a certain degree of deformation. Undercuts help protect the threaded holes from distortion. Significant deformation of the VA locking holes reduces the locking effectiveness.
- Reverse bending or use of the incorrect instrumentation for bending may weaken the plate and lead to premature plate failure (e.g. breakage). Do not bend the plate beyond what is required to match the anatomy.



### 3

#### Position plate

---

#### Optional instruments

---

292.120	Kirschner Wire Ø 1.25 mm with trocar tip, length 150 mm, Stainless Steel
02.111.500.01(S)	Plate Reduction Wire Ø 1.25 mm, with thread, with Small Stop, length 150 mm, Stainless Steel
02.111.501.01(S)	Plate Reduction Wire Ø 1.25 mm, with thread, with Large Stop, length 150 mm, Stainless Steel
399.970	Reduction Forceps with Points, ratchet lock

---

Position the plate over the reduced fracture and, if necessary, fix provisionally with 1.25 mm Kirschner wires or reduction forceps. (1)

#### Option: Plate reduction wires



1.25 mm plate reduction wires can be used for preliminary plate fixation. (2)

These must be removed when no longer needed for temporary fixation.

---

**Precaution:** The plate reduction wires and Kirschner wires are single-use items, do not re-use.

---



# Screw Insertion

## Cortex screws

### 1

#### Drill screw hole for cortex screw

---

#### Instruments

---

310.509	Drill Bit $\varnothing$ 1.8 mm with marking, length 110/85 mm, 2-flute, for Quick Coupling
---------	--

---

323.202	Universal Drill Guide 2.4
---------	---------------------------

---

Drill the screw hole through the shaft of the plate using the 1.8 mm drill bit and the universal drill guide.



### 2

#### Determine screw length

---


#### Instrument

---

03.111.005	Depth Gauge for Screws $\varnothing$ 2.0 to 2.7 mm, measuring range up to 40 mm
------------	---

---

Determine the screw length with the depth gauge.



### 3

#### Insert cortex screw

---

#### Instruments

---

314.467	Screwdriver Shaft, Stardrive, T8, self-holding
---------	--

---

311.430	Handle with Quick Coupling, length 110 mm
---------	---

---


#### Optional instrument

---

314.453	Screwdriver Shaft, Stardrive, 2.4, short, self-holding, for Quick Coupling
---------	--

---

Insert the self-tapping cortex screw using the self-holding T8 Stardrive screwdriver shaft and quick coupling handle.



# Screw Insertion

## Variable angle locking screws

### 1a

#### Drill screw hole for VA locking screw using variable angle technique

---

#### Instruments


310.509	Drill Bit Ø 1.8 mm with marking, length 110/85 mm, 2-flute, for Quick Coupling
03.110.000	VA-LCP Drill Sleeve 2.4, for Drill Bits Ø 1.8 mm

---

#### Optional instruments

03.110.023	VA-LCP Drill Sleeve 2.4, conical, for Drill Bits Ø 1.8mm
03.111.004	VA-LCP Drill Sleeve 2.4, for Drill Bits Ø 1.8 mm, freehand useable

---



#### Drill using VA-LCP drill sleeve with funnel

Insert and lock the VA-LCP drill sleeve tip into the cloverleaf design of the VA locking hole.

Use the 1.8 mm drill bit to drill to the desired depth at the desired angle.

The funnel of the drill sleeve allows the drill bit to be angled up to 15° around the central axis of the locking hole.

---

#### Note:

- The drill guide inserts co-axially into the ole. Ensure that the tip of the drill guide remains fully seated in the hole while drilling.
  - When using the cone-end of the variable angle drill guide, measurement cannot be taken with the 1.8 mm drill bit with depth mark. The depth gauge must be used.
-

---

**Drill using VA-LCP drill sleeve for freehand use**

Alternatively, use the freehand VA-LCP drill sleeve. Fully extend it into the VA locking hole. Drill variable angle holes at the desired angle.

To ensure that the screw is locked correctly, do not angle it in excess of  $\pm 15^\circ$  from the nominal trajectory of the hole.

- To achieve the desired angle, verify the drill bit angle under image intensifier control. If necessary, drill at a different angle and verify again under image intensifier control.

---

**Note:** The previously inserted Kirschner wire can be used as a reference for the screw angulation by using the image intensifier.

---



## 1b

### Drill screw hole for VA locking screw using nominal angle technique

---

#### Instruments

---

310.509	Drill Bit $\varnothing$ 1.8 mm with marking, length 110/85 mm, 2-flute, for Quick Coupling
---------	--

---

03.110.000	VA-LCP Drill Sleeve 2.4, for Drill Bits $\varnothing$ 1.8 mm
------------	--

---

#### Optional instrument

---


03.110.024	VA-LCP Drill Sleeve 2.4, coaxial, for Drill Bits $\varnothing$ 1.8 mm
------------	---

---

The fixed-angle end of the drill sleeve only allows the drill bit to follow the nominal trajectory of the VA locking hole.

Use the 1.8 mm drill bit to drill to the desired depth.

Read the screw length directly from the laser mark on the drill bit. Alternatively, use the depth gauge to determine the screw length.



## 2

### Determine screw length

---


#### Instrument

---

03.111.005      Depth Gauge for Screws  $\varnothing$  2.0 to 2.7 mm,  
measuring range up to 40 mm

---

Determine the screw length with the depth gauge.



### 3

#### Insert VA locking screws

---

##### Instruments

---

314.467	Screwdriver Shaft, Stardrive T8, self-holding
---------	---

---

311.430	Handle with Quick Coupling, length 110 mm
---------	---

---

##### Optional instrument

---

314.453	Screwdriver Shaft, Stardrive, 2.4, short, self-holding, for Quick Coupling
---------	--

---

Insert the VA locking screws manually with the self-holding T8 Stardrive screwdriver shaft and quick coupling handle and tighten just enough for the screw head to be fully seated in the locking hole.

When using the pre-defined nominal angle technique, standard locking screws can also be used instead of VA locking screws.

---

**Note:** Do not over-tighten the screws. This allows the screws to be easily removed if they are not in the desired position.

---




## 4

### Ensure proper joint reconstruction

- After insertion of screws, ensure proper joint reconstruction, screw placement and screw length using the image intensifier. Verify that the distal screws are not in the joint by using additional views.

In an AP view, the dorso-ulnar plate should be projected almost antero-posteriorly, the dorso-radial plate almost laterally, and vice versa for the lateral view. If the plates appear to be parallel, the dorso-radial plate is positioned too far on the ulnar side.



## 5 Final fixation of VA locking screws

### Instruments

03.110.005	Handle for Torque Limiters 0.4/0.8/1.2 Nm
511.776	Torque Limiter, 0.8 Nm, with AO/ASIF Quick Coupling
314.467	Screwdriver Shaft, Stardrive, T8, self-holding


### Optional instrument

314.453	Screwdriver Shaft, Stardrive, 2.4, short, self-holding, for Quick Coupling
---------	---

**Precaution:** Use of the 0.8 Nm torque limiter (TLA) is mandatory when inserting locking screws into variable angle locking holes to ensure the adequate torque is applied (1). Final locking must be done manually using the TLA.

The torque limiter prevents over-tightening and ensures that the VA locking screws are securely locked into the plate. (2)

**Note:** For dense bone, visually inspect if the screw is countersunk after tightening with the torque limiter. If required, carefully tighten without the torque limiter until the screw head is flush with the plate surface.



# Postoperative Treatment / Implant Removal

---

## Postoperative treatment

Postoperative treatment with VA locking compression plates does not differ from conventional internal fixation procedures.

## Implant removal

---

### Instruments

311.430	Handle with Quick Coupling, length 110 mm
314.467	Screwdriver Shaft, Stardrive, T8, self-holding

---


### Optional instrument

314.453	Screwdriver Shaft, Stardrive, 2.4 short, self-holding, for Quick Coupling
---------	--

---

To remove locking screws, first unlock all screws from the plate; then remove the screws completely from the bone.

The last screw removed should be a non-locking screw on the shaft. This prevents the plate from spinning when locking screws are removed.



# Plates

## VA-LCP Dorsal Distal Radius Plate 2.4, Radial Column

Part number	Head holes	Length (mm)
OX.115.530	5	46
OX.115.540	6	57



## VA-LCP Dorsal Distal Radius Plate 2.4, Intermediate Column

**Note:** The plates for the right radius (OX.115.630 and OX.115.640) are left angled and the plates for the left radius (OX.115.631 and OX.115.641) are right angled.

Part number	Head holes	Shaft holes	Length (mm)	For radius
OX.115.630	2	3	41	Right
OX.115.631	2	3	41	Left
OX.115.640	2	4	49	Right
OX.115.641	2	4	49	Left



## VA-LCP Dorsal Distal Radius L-Plate 2.4

Part number	Head holes	Shaft holes	Length (mm)	Angle
OX.115.130	2	3	37	Right
OX.115.131	2	3	37	Left
OX.115.150	2	5	51	Right
OX.115.151	2	5	51	Left



---

**VA-LCP Dorsal Distal Radius L-Plate 2.4**

Part number	Head holes	Shaft holes	Length (mm)	Angle
OX.115.230	3	3	37	Right
OX.115.231	3	3	37	Left
OX.115.250	3	5	51	Right
OX.115.251	3	5	51	Left

**VA-LCP Dorsal Distal Radius L-Plate 2.4, oblique**

Part number	Head holes	Shaft holes	Length (mm)	Angle
OX.115.430	3	3	41	Right
OX.115.431	3	3	41	Left
OX.115.450	3	5	55	Right
OX.115.451	3	5	55	Left

**VA-LCP Dorsal Distal Radius T-Plate 2.4**

Part number	Head holes	Shaft holes	Length (mm)
OX.115.330	3	3	37
OX.115.350	3	5	51



All plates are also available sterile packed. Add suffix "S" to article number.

X = 2: Stainless steel (SS)

X = 4: Titanium (TiCp)

# Screws

---

## Variable Angle Locking Screws 2.4 mm

OX.210.108 – VA Locking Screw Stardrive  $\varnothing$  2.4 mm,  
OX.210.130 self-tapping, lengths 8 mm to 30 mm

For use in VA locking holes.



---

**Precaution:** For final locking, the 0.8 Nm torque limiter is required.

---

## Cortex Screws 2.4 mm

X01.756 – Cortex Screw Stardrive  $\varnothing$  2.4 mm,  
X01.780 self-tapping, lengths 6 mm to 30 mm

For use in VA locking holes or combi-holes.



All screws are also available sterile packed. Add suffix "S" to article number.

X = 2: Stainless steel (SS)

X = 4: Titanium Alloy (TAN)

---

## Optional

### Variable Angle Locking Buttress Pins 1.8 mm

---

OX.210.078 – VA-LCP Buttress Pin Stardrive  $\varnothing$  1.8 mm,  
OX.210.100 lengths 8 mm to 30 mm

For use in VA locking holes.

---

**Precaution:** For final locking, the 0.8 Nm torque limiter is required.

---



### Locking Screws 2.4 mm

---

X12.806 – Locking Screw Stardrive  $\varnothing$  2.4 mm,  
X12.830 self-tapping, lengths 6 mm to 30 mm

For use in VA locking holes but only in pre-defined angle using nominal angle technique.

---

**Precaution:** For final locking, the 0.8 Nm torque limiter is required.

---



All screws are also available sterile packed. Add suffix "S" to article number.

X = 2: Stainless steel  
X = 4: TAN

# Instruments

03.110.000	VA-LCP Drill Sleeve 2.4, for Drill Bits $\varnothing$ 1.8 mm	
323.202	Universal Drill Guide 2.4	
310.509	Drill Bit $\varnothing$ 1.8 mm with marking, length 110/85 mm, 2-flute, for Quick Coupling	
314.453	Screwdriver Shaft, Stardrive 2.4, short, self-holding, for Quick Coupling	
314.467	Screwdriver Shaft, Stardrive, T8, self-holding	
03.111.005	Depth Gauge for Screws $\varnothing$ 2.0 to 2.7 mm, measuring range up to 40 mm	
311.430	Handle with Quick Coupling, length 110 mm	
03.110.005	Handle for Torque Limiters 0.4/0.8/1.2 Nm	
511.776	Torque Limiter 0.8 Nm, with AO/ASIF Quick Coupling	
292.120(S)	Kirschner Wire $\varnothing$ 1.25 mm with trocar tip, length 150 mm, Stainless Steel	

---

**Optional Instruments**

---

03.111.038 Handle with Quick Coupling



---

03.110.023 VA-LCP Drill Sleeve 2.4, conical,  
for Drill Bits  $\varnothing$  1.8 mm



---

03.110.024 VA-LCP Drill Sleeve 2.4, coaxial,  
for Drill Bits  $\varnothing$  1.8 mm



---

03.111.004 VA-LCP Drill Sleeve 2.4,  
for Drill Bits  $\varnothing$  1.8 mm, freehand useable



---

02.111.500.01(S) Plate Reduction Wire  $\varnothing$  1.25 mm,  
with thread, with Small Stop,  
length 150 mm, Stainless Steel



---

02.111.501.01(S) Plate Reduction Wire  $\varnothing$  1.25 mm,  
with thread, with Large Stop,  
length 150 mm, Stainless Steel



## Bibliography

---

Blythe M et al. (2006) Volar Versus Dorsal Locking Plates With and Without Radial Styloid Locking Plates for the Fixation of Dorsally Comminuted Distal Radius Fractures: A Biomechanical Study in Cadavers. *J Hand Surg* 31A: 1587–1593

Letsch R et al. (2003) Surgical treatment of fractures of the distal radius with plates: a comparison of palmar and dorsal plate position. *Arch Orthop Trauma Surg* 123: 333–339

Lutsky K et al. (2009) Dorsal Fixation of Intra-articular Distal Radius Fractures Using 2.4mm Locking Plates. *Tech Hand Surg* 13: 187–196

Peine R et al. (2000) Comparison of Three Different Plating Techniques for the Dorsum of the Distal Radius: A Biomechanical Study. *J Hand Surg* 25A: 29–33

Rikli DA (2009) Dorsal Double Plating and Combined Palmar and Dorsal Plating for Distal Radius Fractures. In: Slutsky DJ, Osterman AL (2009) *Fractures and Injuries of the Distal Radius and Carpus*. Saunders Elsevier, 125–133

Rikli DA, Businger A, Babst R (2005) Dorsal double-plate fixation of the distal radius. *Oper Orthop Traumatol* 17(6): 624–640

Rikli DA, Regazzoni P (2000) The double plating technique for distal radius fractures. *Hand and upper extremity surgery* 4:101–114

Rikli DA, Regazzoni P (1996) Fractures of the distal end of the radius treated by internal fixation and early function. A preliminary report of 20 cases. *J Bone Joint Surg [Br]* 78-B (4): 588–592

Tavakolian JD, Jupiter JB (2005) Dorsal Plating for Distal Radius Fractures. *Hand Clin* 21: 341–346

---

## **Torque, Displacement and Image Artifacts according to ASTM F 2213-06, ASTM F 2052-06e1 and ASTM F 2119-07**

Non-clinical testing of worst case scenario in a 3 T MRI system did not reveal any relevant torque or displacement of the construct for an experimentally measured local spatial gradient of the magnetic field of 3.69 T/m. The largest image artifact extended approximately 169 mm from the construct when scanned using the Gradient Echo (GE). Testing was conducted on a 3 T MRI system.

## **Radio-Frequency-(RF-)induced heating according to ASTM F 2182-11a**

Non-clinical electromagnetic and thermal testing of worst case scenario lead to peak temperature rise of 9.5 °C with an average temperature rise of 6.6 °C (1.5 T) and a peak temperature rise of 5.9 °C (3 T) under MRI Conditions using RF Coils (whole body averaged specific absorption rate [SAR] of 2 W/kg for 6 minutes [1.5 T] and for 15 minutes [3 T]).

---

**Precautions:** The above mentioned test relies on non-clinical testing. The actual temperature rise in the patient will depend on a variety of factors beyond the SAR and time of RF application. Thus, it is recommended to pay particular attention to the following points:

- It is recommended to thoroughly monitor patients undergoing MR scanning for perceived temperature and/or pain sensations.
  - Patients with impaired thermoregulation or temperature sensation should be excluded from MR scanning procedures.
  - Generally, it is recommended to use a MR system with low field strength in the presence of conductive implants. The employed specific absorption rate (SAR) should be reduced as far as possible.
  - Using the ventilation system may further contribute to reduce temperature increase in the body.
-

